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LITERATURA

- Stamp D and Jenkins G. An overview of bile- acid synthesis, chemistry and function. In: Jenkins GJ, Hardie L, editors. Bile Acids: Toxicology and Bioactivity. Royal Society of Chemistry; 2008. p. 1-13.

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2. Hofmann AF. The continuing importance of bile acids in liver and intestinal disease. Arch. Intern. Med. 1999;159:2647– 2658.
3. Dufour DR. Liver Disease. In: Burtis CA, Ashwood ER, and Bruns DE eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 4th ed. Philadelphia: Elsevier Inc; 2006. p. 1782–1787.
4. Shima, T, Tada H, Morimoto M, et al. Serum total bile acid level as a sensitive indicator of hepatic histological improvement in chronic hepatitis C patients responding to interferon treatment. J. gastroenterol. Hepatol. 2000;15: 294-299.
5. Barnes S, Gallo GA, Trash DB, et al. Diagnostic value of serum bile acid estimations in liver disease. J Clin Pathol. 1975;28:506-9.
6. Skrede S, Solberg HE, Blomhoff JP, et al. Bile acids measured in serum during fasting as a test for liver disease. Clinical Chemistry.1978;24:1095-1099.
7. Pataia V, Dixon PH, Williamson C. Pregnancy and bile acid disorders. Am. J. Physiol. Gastrointest. Liver Physiol. 2017;313: G1–G6.
8. Neale G, Lewis B, Weaver V, et al. Serum bile acids in liver disease. Gut. 1971;12:145-152.
9. Charach G, Argov O, Geiger K, et al. Diminished bile acids excretion is a risk factor for coronary artery disease: 20-year follow up and long-term outcome. Therap. Adv. Gastroenterol. 2017;11:1-11.
10. Egan N et al. Reference standard for serum bile acids in pregnancy. BJOG. 2012;119:493–498.
11. Glantz A, Marschall HU, Mattson LA. Intrahepatic cholestasis of pregnancy: relationships between bile acid levels and fetal complication rates. Hepatology. 2004;40:467–74.
12. Brites D. Intrahepatic cholestasis of pregnancy: changes in maternal-fetal bile acid balance and improvement by ursodeoxycholic acid. Ann. Hepatol. 2002;1:20–28.
13. Ozkan S, Ceylan Y, Ozkan OV, et al. Review of a challenging clinical issue: Intrahepatic cholestasis of pregnancy. World J Gastroenterol. 2015;21:7134–41.
14. Schlesinger DP, Rubin SI. Serum bile acids and the assessment of hepatic function in dogs and cats. 1993;34:6.
15. Zhang GH, Cong AR, Xu GB, et al. An enzymatic cycling method for the determination of serum total bile acids with Total bile acids 21 FS – Page 3 844 2238 10 02 00 September 2022/3 recombinant 3alpha-hydroxysteroid dehydrogenase. Biochem Biophys Res Commun. 2005;326:87-92.
16. Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med. 2007;45:1240-3.
17. Brown RS. Use of Obeticholic Acid in Patients With Primary Biliary Cholangitis. Gastroenterol Hepatol. 2018;14:654-657.
18. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th ed. Volume 1 and 2. Washington DC: The American Association for Clinical Chemistry Press; 2000.
19. Young DS. Effects on Clinical Laboratory Tests - Drugs Disease, Herbs & Natural Products [Internet]. AACC Press and John Wiley and Sons, Inc; 2020 [cited 2020 Sept]. Available from: <https://clinfo.wiley.com/aaccweb/aacc/>.
20. Luo L, Aubrecht J, Li D, Warner RL, et al. Assessment of serum bile acid profiles as biomarkers of liver injury and liver disease in humans. PLoS One. 2018;13.
21. Jahnel J, Zöhrer E, Schamagel H, et al. Reference ranges of serum bile acids in children and adolescents. Clin Chem Lab Med. 2015;53:1807-13.

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Fabricante: DiaSys Diagnostic Systems GmbH