

Instructions for Use

For *in vitro* diagnostic use



SOLUÇÃO HEMOLISANTE

Hemolyzing solution

Anvisa 80115310147

ORDER INFORMATION

Cat. No.	Kit size
4190500K	1 x 500 mL

INTENDED USE

Hemolyzing solution for quantitative *in vitro* determination of hemoglobin A1c (HbA1c) in whole blood on photometric systems.

SUMMARY

Hemolyzing Solution is used for sample preparation of whole blood, calibrators, and controls for determination of HbA1c with Kovalent HbA1c Kit.

COMPONENTS AND CONCENTRATIONS

Water
Stabilizer.

STORAGE INSTRUCTIONS AND STABILITY

The solution is stable up to the end of the indicated month of expiry, if stored at 2 – 25°C and contamination is avoided. Do not freeze the solution.

WARNINGS AND PRECAUTIONS

Take necessary precautions when handling laboratory reagents.

PREPARATION

The reagent is ready to use.

SAMPLE PREPARATION PROCEDURE

Hemolyzing solution: 1000µL
Sample/Calibrator/Control: 20µL
Mix and let stand for 5 minutes or until complete hemolysis is evident.

WASTE MANAGEMENT

Follow the requirements of the current guidelines about technical regulation for the management of healthcare service waste, as well as other equivalent biosafety practices.

WARRANTY

These instructions for use should be read carefully before using the product and the information contained therein should be strictly adhered to. The reliability of the test results cannot be guaranteed if the instructions are not followed.

CONSUMER INFORMATION

Symbols used:

	Manufacturer
	Temperature limit
	In vitro diagnostic device
	Caution
	Operating instructions
	Recycling material
	Do not discard directly into the environment
	Batch code
	Date of manufacture
	Use by date
	Biological hazards
	Highly toxic
	Corrosive
	Harmful

Manufacturer:

Kovalent do Brasil Ltda.

Rua Cristóvão Sardinha, 110 – Jd. Bom Retiro
São Gonçalo – RJ – CEP 24722-714
www.kovalent.com.br
CNPJ: 04.842.199/0001-56

Customer service: sac@kovalent.com.br - (21) 3907-2534 / 0800 015 1414

Expiration date and Lot no.: See label