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GLICOSE GOD-PAP

Anvisa 80115310038

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LITERATURA

- Hallbach J. Klinische Chemie und Hämatologie – Biomedizinische Analytik für MTLA und Studium. 3rd ed. Stuttgart: Georg Thieme Verlag KG; 2011. p. 170-171.
- Thomas L. Clinical Laboratory Diagnostics [Internet]. Prof. Lothar Thomas; 2023 [cited 2023 Nov 21]. Available from: <https://www.clinical-laboratory-diagnostics.com/>
- Barham D, Trinder P. An improved color reagent for the determination of blood glucose by the oxidase system. Analyst 1972; 97: 142-5.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. ClinChemLabMed 2007;45(9):1240-1243.
- Guder WG et al. Die Qualität diagnostischer Proben – Empfehlung der Arbeitsgruppe Präanalytik der Deutschen Vereinten Gesellschaft für Klinische Chemie und Laboratoriumsmedizin. 7th ed. Heidelberg: BD Diagnostics Preanalytical Systems; 2012. p. 46-47, p. 68-69.
- Sacks DB. Carbohydrates. In: Burtis CA, Ashwood ER, Bruns DE, editors. Tietz Textbook of Clinical Chemistry and Molecular Diagnosis. 4th ed. St. Louis, Missouri: Elsevier Saunders Company; 2006. p. 837-901.
- Sacks DB, Bruns DE, Goldstein DE, MacLaren NK, McDonald JM, Parrott M. Guidelines and Recommendations for Laboratory Analysis in the Diagnosis and Management of Diabetes Mellitus. ClinChem 2002; 48: 436-472.
- Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press 2000.
- Young DS. Effects on Clinical Laboratory Tests - Drugs Disease, Herbs & Natural Products, <https://clinf.wiley.com/aaccweb/aacc/>, accessed in February 2024. Published by AACC Press and John Wiley and Sons, Inc.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem. 2001 Jul;38:376-85.

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