

# Instructions for Use

For *in vitro* diagnostic use



## CREATININA WS

### Creatinine WS

Anvisa 80115310205

#### ORDER INFORMATION

| Cat. No.    | Kit size                        |
|-------------|---------------------------------|
| 1030500KWS  | R1 2 x 200 mL + R2 1 x 100 mL   |
| 1030150MWS  | R1 4 x 30 mL + R2 2 x 15 mL     |
| 1030200RWS  | R1 4 x 38,6 mL + R2 4 x 11,4 mL |
| 1030100MKWS | R1 2 x 40 mL + R2 2 x 10 mL     |
| 1030250KWS  | R1 1 x 200 mL + R2 1 x 50 mL    |

#### INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of creatinine in serum, plasma or urine on photometric systems.

#### SUMMARY

Creatinine is a waste product excreted by the kidneys mainly by glomerular filtration. The concentration of creatinine in plasma of a healthy individual is fairly constant, independent from water intake, exercise and rate of urine production. Therefore, increased plasma creatinine values always indicate decreased excretion, i.e. impaired kidney function. Creatinine clearance is a good indicator for the glomerular filtration rate (GFR) which allows better detection of kidney diseases and monitoring of renal function. For this purpose, creatinine is measured simultaneously in serum and urine collected over a defined time period. [1,2]

#### METHOD

Kinetic test without deproteinization according to the Jaffé method

#### PRINCIPLE

Creatinine forms a colored orange-red complex in an alkaline picrate solution. The difference in absorbance at fixed times during conversion is proportional to the concentration of creatinine in the sample.

Creatinine + Picric acid  $\longrightarrow$  Creatinine picrate complex

#### REAGENTS

##### Components and Concentrations

|    |                  |           |
|----|------------------|-----------|
| R1 | Sodium hydroxide | 0.2 mol/L |
| R2 | Picric acid      | 20 mmol/L |

#### WARNINGS AND PRECAUTIONS

1. Reagent 1: Warning! May be corrosive to metals. It causes skin irritation. Causes serious eye irritation, if irritation persists, seek medical advice. Keep only in original packaging. Wash hands and face thoroughly after handling. Wear protective gloves/protective clothing/eye protection. IF ON SKIN: Wash with plenty of water/soap. If skin irritation occurs: Get medical advice/attention. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do, and continue rinsing. Absorb spillage to prevent material damage
2. Reagent 2: Warning! May be corrosive to metals. Keep only in original packaging. Wear protective gloves/protective clothing/eye protection. Absorb spillage to prevent material damage.
3. High homogentisic acid concentrations in urine samples lead to false results.
4. In very rare cases, samples of patients with gammopathy might give falsified results [11].
5. Eltrombopag medication (active ingredient of the drug Revolade®) leads to falsely low or high results in patient samples.
6. Please refer to the safety data sheets (SDS) and take the necessary precautions for the use of laboratory reagents. For diagnostic purposes, the results should always be assessed with the patient's medical history, clinical examinations and other findings.
7. For professional use only.

#### STORAGE AND STABILITY

Reagents are stable up to the date of expiry indicated on the kit, if stored at 2 - 8°C, protected from light and if contamination is avoided. Do not freeze the reagents!

#### WASTE MANAGEMENT

Follow the requirements of the current guidelines about technical regulation for the management of healthcare service waste, as well as other equivalent biosafety practices.

#### REAGENT PREPARATION

##### Starting with Substrate

The reagents are ready to use

##### Starting with Sample

Mix 4 parts of R1 with 1 part of R2  
(e.g. 20 mL R1 + 5 mL R2) = monoreagent

Stability: 5 hours at 15-25°C

#### MATERIALS REQUIRED BUT NOT PROVIDED

1. NaCl solution 9 g/L.
2. General laboratory equipment.

#### SPECIMEN

Serum, heparin plasma or urine

Stability [5]:

|               |                   |    |            |
|---------------|-------------------|----|------------|
| Serum/plasma: | 7 days            | at | 4 - 25 °C  |
|               | at least 3 months | at | -20 °C     |
| Urine:        | 2 days            | at | 20 - 25 °C |
|               | 6 days            | at | 4 - 8 °C   |
|               | 6 months          | at | -20 °C     |

Dilute urine 1 + 49 with distilled water; multiply the result by 50.

Discard contaminated specimens! Only freeze once!

#### ASSAY PROCEDURE

Applications for automatic systems are available upon request or on our website: [www.kovalent.com.br](http://www.kovalent.com.br)

|              |                         |
|--------------|-------------------------|
| Wavelength   | Hg 492nm (490 - 510 nm) |
| Optical path | 1 cm                    |
| Temperature  | 20 - 25 °C / 37 °C      |
| Measurement  | Against reagent blank   |

##### Starting with Substrate

|  | Blank   | Sample or calibrator |
|--|---------|----------------------|
| Sample or calibrator   | -       | 50 µL                |
| Distilled water  | 50 µL   | -                    |
| Reagent 1  | 1000 µL | 1000 µL              |
| Mix, incubate for 0 - 5 min, then add:   |         |                      |
| Reagent 2  | 250 µL  | 250 µL               |
| Mix and read the absorbance A1 after 60 seconds, then read the absorbance A2 after the next 120 seconds. |         |                      |

$$\Delta A = (A_2 - A_1)_{\text{sample or calibrator}}$$

##### Starting with Sample

|  | Blank   | Sample or calibrator |
|--|---------|----------------------|
| Sample or calibrator   | -       | 50 µL                |
| Distilled water  | 50 µL   | -                    |
| Monoreagent  | 1000 µL | 1000 µL              |
| Mix and read the absorbance A1 after 60 seconds, then read the absorbance A2 after the next 120 seconds. |         |                      |

$$\Delta A = (A_2 - A_1)_{\text{sample or calibrator}}$$

#### CALCULATION

##### With calibrator

##### Serum/plasma

$$\text{Creatinine [mg/dL]} = \frac{\Delta A_{\text{Sample}}}{\Delta A_{\text{Cal}}} \times \text{Conc Cal. [mg/dL]}$$

##### Urine

$$\text{Creatinine [mg/dL]} = \frac{\Delta A_{\text{Sample}}}{\Delta A_{\text{Cal}}} \times \text{Conc Cal. [mg/dL]} \times 50$$

##### Creatinine Clearance [mL/min/1.73 m<sup>2</sup>] [7]

$$= \frac{\text{mg Creatinine} / 100 \text{ mL Urine} \times \text{mL Urine}}{\text{mg Creatinine} / 100 \text{ mL Serum} \times \text{min (urine collection time)}}$$

# Instructions for Use

For *in vitro* diagnostic use

The calculated creatinine clearance refers to the average body surface of an adult (1.73 m<sup>2</sup>).

## Conversion Factor

Creatinine [mg/dL] x 88.4 = Creatinine [μmol/L]

## COMPENSATED METHOD [3,4]

Picric acid which forms the colored complex reacts unspecifically with interfering serum components, so-called pseudo-creatinines. This leads to falsely elevated creatinine values in serum and plasma samples especially in the low measuring range. To compensate these interferences, the calibrator value for the compensated method indicated in the value sheet of Topkal U has to be used for calculation. Additionally, 0.3 mg/dL has to be subtracted from the calculated creatinine value. For use of the compensated method, calibration with the calibrator Topkal U is strictly recommended. The method is applicable only for serum and plasma samples. The compensated method is traceable to GC-IDMS.

## CALIBRATORS AND CONTROLS

For calibration in automated photometric systems, Kovalent Topkal U calibrator is recommended. Use Kovalent Topkon N and P for internal quality control. Each laboratory should establish corrective action in case of deviations in control recovery.

## WARRANTY

These instructions for use should be read carefully before using the product and the information contained therein should be strictly adhered to. The reliability of the test results cannot be guaranteed if the instructions are not followed.

## PERFORMANCE CHARACTERISTICS

### Measuring range

The test is designed to determine creatinine concentrations within a measuring range of 0.2 - 15 mg/dL (18 - 1330 μmol/L). When values exceed this range, samples should be diluted 1 + 1 with NaCl solution (9 g/L) and the result multiplied by 2.

### Specificity / Interferences

No interference was observed by ascorbic acid up to 30 mg/dL, hemoglobin up to 500 mg/dL and lipemia up to 2000 mg/dL of triglycerides. Bilirubins interfere from a concentration of 4 mg/dL. For more information on interfering substances, see Young DS [10].

### Sensitivity / Limit of Detection

The lowest detection limit is 0.2 mg/dL (17.7 μmol/L).

### Precision

| Within run<br>n = 10 | Mean<br>[mg/dL] | SD<br>[mg/dL] | CV<br>[%] |
|----------------------|-----------------|---------------|-----------|
| Normal control       | 1.44            | 0.04          | 3.13      |
| Pathological control | 4.089           | 0.06          | 1.50      |
| Between day<br>n = 9 | Mean<br>[mg/dL] | SD<br>[mg/dL] | CV<br>[%] |
| Normal control       | 1.51            | 0.02          | 1.22      |
| Pathological control | 4.26            | 0.08          | 1.93      |

### Method comparison

Method comparison between Kovalent Creatinina (y) and a commercially available test with Jaffé method (x) using 30 samples of human serum demonstrated the following results:  
y = 0.9939x + 0.0067; R<sup>2</sup> = 0.9839

## REFERENCE RANGE

### Serum or plasma – Jaffé method, not compensated

|                       | mg/dL     | μmol/L   |
|-----------------------|-----------|----------|
| <b>Adults [1]</b>     |           |          |
| Women                 | 0.6 – 1.1 | 53 – 97  |
| Men                   | 0.7 – 1.3 | 62 – 115 |
| <b>Children [2,8]</b> |           |          |
| Neonate               | 0.5 – 1.2 | 44 – 106 |
| Infant                | 0.4 – 0.7 | 35 – 62  |
| Child                 | 0.5 – 1.2 | 44 – 106 |

### Serum or plasma – Jaffé method, compensated

|                   | mg/dL     | μmol/L   |
|-------------------|-----------|----------|
| <b>Adults [3]</b> |           |          |
| Women             | 0.5 – 0.9 | 44 – 80  |
| Men               | 0.7 – 1.2 | 62 – 106 |

### Children [9]

|         |             |         |
|---------|-------------|---------|
| Neonate | 0.24 – 1.04 | 21 – 92 |
| Infant  | 0.17 – 0.42 | 15 – 37 |
| Child   | 0.24 – 0.87 | 21 – 77 |

### 1st urine of the morning [3] – Jaffé method, compensated

|                      | mg/dL      | μmol/L       |
|----------------------|------------|--------------|
| Women                | 28 – 217   | 2470 – 19200 |
| Men                  | 39 – 259   | 3460 – 22900 |
|                      | mg/kg/24 h | μmol/kg/24 h |
| <b>Urine 24h [1]</b> |            |              |
| Women                | 11 – 20    | 97 – 177     |
| Men                  | 14 – 26    | 124 – 230    |

### Albumin/creatinine ratio (early morning urine) [12]:

< 30 mg/g Creatinine

### Creatinine Clearance [2]

|       | mL/min/1.73 m <sup>2</sup> |
|-------|----------------------------|
| Women | 95 – 160                   |
| Men   | 98 – 156                   |

Each laboratory should check if the reference ranges are transferable to its own patient population and determine its own reference ranges if necessary.

## LITERATURE















- Newman DJ, Price CP. Renal function and nitrogen metabolites. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1204-.
- Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 366-74.
- Mazzachi BC, Peake MJ, Ehrhardt V. Reference Range and Method Comparison Studies for Enzymatic and Jaffé Creatinina Assays in Plasma and Serum and Early Morning Urine. Clin. Lab. 2000; 46: 53-55
- Swanson AF, Swartzentruber M, Nolen PA et al. Multicenter Evaluation of the Boehringer Mannheim Compensated, Rate-Blanked Creatinine/Jaffe Application on BM/Hitachi Systems. Advances in Clinical Diagnostics. 1993. Boehringer Mannheim Corporation.
- Guder WG, Zawta B. Recommendations of the Working group on Preanalytical Quality of the German Society for Clinical Chemistry and the German Society for Laboratory Medicine: The Quality of Diagnostic Samples. 1<sup>o</sup> ed Darmstadt: GIT Verlag 2001; p. 24-5,50-1.
- Levey AS, Coresh J, Greene T, Marsh J et al: Expressing the Modification of Diet in Renal Disease Study Equation for Estimating Glomerular Filtration Rate with Standardized Serum Creatinine Values. Clin Chem 2007; 53 (4): 766-72.
- Junge W, Wilke B, Halabi A, Klein G. Determination of reference intervals for serum creatinine, creatinine excretion and creatinine clearance with an enzymatic and a modified Jaffé method. Clin Chim Acta 2004; 344: 137-148.
- Soldin SJ, Brugnara C, Wong EC, editores. Pediatric Reference Intervals. 6<sup>o</sup> ed. AACC Press, 2007: p. 77-78.
- Schlebusch H, Liappis N, Klein G. Ultrasensitive CRP and Creatinine: Reference intervals from infancy to childhood. Clin Chem Lab Med. 2001; 39 Special supplement pp S1-S448; Maio 2001. PO-T042.
- Young DS. Effects of Drugs on Clinical Laboratory Tests. 5<sup>th</sup> ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press 2000.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. ClinChemLabMed 2007;45(9):1240-1243.
- Dati F, Metzmann E. Proteins-Laboratory testing and clinical use. 1st ed. Holzheim: DiaSys Diagnostic Systems; 2005: p. 93.

# Instructions for Use

For *in vitro* diagnostic use

## CONSUMER INFORMATION

Symbols used:

|  |  |
|--|--|
|   | Manufacturer                                 |
|   | Temperature limit                            |
|   | <i>In vitro</i> diagnostic device            |
|   | Caution                                      |
|   | Operating instructions                       |
|   | Recycling material                           |
|   | Do not discard directly into the environment |
|   | Batch code                                   |
|   | Date of manufacture                          |
|   | Use by date                                  |
|   | Biological hazards                           |
|   | Highly toxic                                 |
|   | Corrosive                                    |
|  | Harmful                                      |

### Manufacturer:

#### Kovalent do Brasil Ltda.

Rua Cristóvão Sardinha, 110 – Jd. Bom Retiro

São Gonçalo – RJ – CEP 24722-414 – Brasil

[www.kovalent.com.br](http://www.kovalent.com.br)

CNPJ: 04.842.199/0001-56

Kit sizes variations on demand:

| Anvisa No.  | Kit size              |
|-------------|-----------------------|
| 80115310205 | R1 4x40mL + R2 4x10mL |
| 80115310205 | R1 1x40mL + R2 1x10mL |
| 80115310205 | R1 3x40mL + R2 3x10mL |

Customer service: [sac@kovalent.com.br](mailto:sac@kovalent.com.br) - (21) 3907-2534 / 0800 015 1414

Expiration date and Lot no.: See label