(EN)Hemoglobin A1c Calibrator**PRODUCT**

Hemoglobin A1c Calibrator

PACKAGE SIZE

REF	Package size
105-007857-00	CAL-1: 2mLx1, CAL-2: 2mLx1
105-007858-00	CAL-1: 2mLx3, CAL-2: 2mLx3

INTENDED USE

The product is used to calibrate HbA1c in whole blood to guarantee that the test result is accurate.

PRINCIPLE

The measurement system is calibrated against calibrator reference values following the calibration procedure. A calibration function is also established. Based on the calibration function, the measurement system calculates the corresponding sample result against measured reaction degree, the units being the same of those for reference values. The precision of the results is determined by the traceability of the calibrators and the performance of the measurement system.

ACTIVE INGREDIENTS

Freeze-dried powder of whole blood matrix of people.

* The content concentrations of products of different batches differ slightly.

The calibrator of HbA1c can be traced to the reference methods of IFCC and DCCT according to the following method.

In the standard test system that uses secondary reference materials for calibration and the IFCC reference method to determine values, determine the reference value of the calibrator according to the standard valuation procedure.

For the IFCC and NGSP reference values of the calibrator, see the reference value table.

STORAGE AND EXPIRATION DATE

Products that are not decapped should be stored in the temperature 2°C to 8°C for 1 year. After the products are re-dissolved and decapped, the products should be stored in the temperature 2°C to 8°C for 7 days.

APPLICABLE INSTRUMENT

This product applies to H50/H50P/H-120/H-120P Automated Glycohemoglobin Analyzer produced by Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

SAMPLE REQUIREMENT

Not applicable.

MATERIALS REQUIRED BUT NOT PROVIDED

The following materials are required but not provided with the product: Mindray-manufactured measurement instruments and matched reagents, controls, as well as general laboratory devices.

For the required reagents and control, refer to the table below.

For H50 and H50P models:

Reagents and control	REF
Hemolysis Solution	105-007856-00
Eluent A	105-007852-00
Eluent B	105-007853-00 105-007854-00 105-007855-00
Hemoglobin A1c Control	105-007859-00 105-007862-00

For H-120 and H-120P models:

Reagents and control	REF
Hemoglobin A1c Hemolysis Solution	105-048455-00
Hemoglobin A1c Eluent A	105-048453-00
Hemoglobin A1c Eluent B	105-048454-00 105-048464-00
Hemoglobin A1c Control	105-007859-00 105-007862-00

INSTRUCTION FOR USE

- Take the product from the environment of 2°C to 8°C and wait until the product reaches the using temperature(10°C to 32°C);
- Gently knock the bottle cap to guarantee that the freeze-dried sample completely falls onto the bottom. Open the bottle cap and rubber stopper with caution to avoid loss of the content;
- Accurately measure 2mL Hemolysis Solution/pure water*, slowly fill it into the bottle, and put the rubber stopper on the bottle. Gently turn over the small bottle for several seconds and place the product for 30 minutes in a dark place;
- Gently turn over the bottom to guarantee that the content is completely dissolved and avoid foam;
- After the content is re-dissolved, select an appropriate centrifuge tube or sample cup to hold the content and mark it for later use;
- Refer to the calibration procedure described in the Operator's Manual. Use the specific values listed in the reference value sheet to set calibration parameters and execute calibration.

* For H50 and H50P models, use Hemolysis Solution or pure water; for H-120 and H-120P models, use only pure water.

CUT-OFF VALUE/REFERENCE INTERVAL

Not applicable.

RESULT ELABORATION

When run a Mindray control on a calibrated Mindray measurement system, the mean of results should fall into the reference range. If it does not, check the measurement system. For instance, check

the expiration dates of reagents, controls and calibrators; and the analyzer's performance. The calibrator values are lot-specific with the matched analyzer models.

LIMITATIONS

Not applicable.


PRODUCT SPECIFICATIONS

The calibrator is pale red brown or pale brown solid before it is re-dissolved. After the calibrator is re-dissolved and placed for a while, the calibrator is red liquid without sediment, suspended matters, or floccules.

PRECAUTIONS

- For professional use in vitro diagnosis only.
- The calibrator values are lot-specific. Please confirm before use.
- Perform calibration after report of outliers, reagent lot change or execution of specific maintenance procedure. Determine the calibration frequency based on laboratory QC rules.
- Perform scheduled maintenance and standard operation including calibration and analysis to assure the analyzer's performance.
- The calibrator was tested and was shown to be negative for antibody to human immunodeficiency virus (HIV), hepatitis B surface antigen (HBsAg), antibody to hepatitis C virus (HCV) and antibody to Treponema pallidum (TP). However, as no testing method can rule out the potential risk of infection with absolute certainty, this material should be handled as a patient sample to avoid biological risk.
- Take the necessary precautions for the use of reagents. Do not swallow. Avoid contact with skin and mucous membranes. If you accidentally take the reagent into your mouth, or the reagents accidentally spill on your skin or into your eyes, wash them off with plenty of water and go see a doctor if necessary.
- Disposal of waste liquid and materials should be in accordance with local guidelines.
- The Material Safety Data Sheet (SDS) is available upon request.
- All identified risks have been reduced as far as possible by generally acknowledged state of art, and the overall residual risk is acceptable.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

GRAPHICAL SYMBOLS

		
Batch code	Use-by date	Manufacturer
		
Unique Device Identifier	Catalogue number	Temperature limit
		
Biological risks	Consult instructions for use	In vitro diagnostic medical device
		
Authorized representative in the European Community	European Conformity	

REFERENCES

Biosafety in Microbiological and Biomedical Laboratories; U.S. Department of Health and Human Services; US Government Printing Office; Washington: 2007.

COMPANY CONTACT

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen 518057, P. R. China

E-mail Address: service@mindray.com

Tel: +86 755 81888998

Fax: +86 755 26582680

AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY

EC-Representative: Shanghai International Holding Corp. GmbH(Europe)

Address: Eiffestraße 80, 20537 Hamburg, Germany

Tel: 0049-40-2513175

Fax: 0049-40-255726

© 2021-2025 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved.

APPROVAL DATE OF THE INSTRUCTION FOR USE

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Fabricante Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Representante da CE: Shanghai International Holding Corp. GmbH(Europe)
© 2012-2025 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Todos os direitos reservados.