REF

H-3 105-048461-00

H-6 105-048462-00

H-10 105-048463-00

Instructions for Use

Hemoglobin A1c Analytical Column

(EN)Hemoglobin A1c Analytical Column	1
(ES)Columna analítica para Hemoglobina A1c	7
(FR)Colonne analytique pour hémoglobine A1c	14
(IT)Colonna analitica per emoglobina A1c	21
(TR)Hemoglobin A1c Analitik Kolonu	28
(RU)Аналитическая колонка для гемоглобина A1c	34
(BRPT)Coluna analítica para hemoglobina A1c	41

Find IFUs in more languages on our website:

<https://www.mindray.com/en/resources-center/ivd-edocs>

**COMPANY CONTACT / KONTAKT MIT DEM
UNTERNEHMEN/CONTACTO CON LA EMPRESA /
COORDONNEES DE LA SOCIETE/INFORMAZIONI
RELATIVE ALLA SOCIETÀ / ŞİRKET İLETİŞİM / CONTATO
DA EMPRESA / КОНТАКТНАЯ ИНФОРМАЦИЯ**

	Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
	Mindray Building, Keji 12 th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China
	www.mindray.com
	service@mindray.com
Tel /Tél/Тел:	+86 755 81888998
Fax/Télécopie/ Faks/Факс:	+86 755 26582680

<table border="1"> <tr> <td>EC</td> <td>REP</td> </tr> </table>	EC	REP	Shanghai International Holding Corp. GmbH(Europe)
	EC	REP	
Eiffestraße 80, Hamburg 20537, Germany			
Tel /Tél/Тел:	0049-40-2513175		
Fax/ Télécopie/ Faks/Факс:	0049-40-255726		

Revision bar indicates update to previous version / Das Revisionsbarsymbol weist auf eine Aktualisierung der vorigen Version / La barra de revisión indica una actualización a la versión anterior / La barre de révision indique la mise à jour vers la version précédente / La barra di revisione indica l'aggiornamento alla versione precedente / Değişiklik çubuğu, önceki sürüme yapılan güncellemeyi gösterir / A barra de revisão indica a atualização para a versão anterior / Строка обновления указывает на обновление предыдущей версии

(EN) Hemoglobin A1c Analytical Column

Product Configuration

Table 1. Configurations of different models

Model	Accessories
H-3	At least one filter
H-6	At least two filters
H-10	At least four filters

Packing Specifications

1 piece / box

Intended Use

Based on ion-exchange high-performance liquid chromatography (IEX HPLC) method, the Hemoglobin A1c Analytical Column is clinically intended to be used with matching reagents and analyzer for quantitative determination of hemoglobin A1c (HbA1c) in human whole blood samples.

Strongly correlated to the chronic complication of diabetes, the HbA1c level is commonly recognized as the "golden index" of glycemic control, and also a good tool assessing the efficacy of the therapy for diabetes. Besides, HbA1c is now acknowledged by the World Health Organization (WHO) and the diabetes associations of many countries as an independent diagnostic criterion for diabetes.

Test Principle

Used together with the Automated Glycohemoglobin Analyzer, Hemoglobin A1c Analytical Column (hereafter, "the column") works based on the ion-exchange high-performance liquid chromatography (IEX HPLC) method to build up a chromatogram and calculate the related parameters.

The blood sample containing various types of hemoglobin is loaded onto the column. In the presence of the low ionic strength eluent (Hemoglobin A1c Eluent A), HbA1c hardly carries any positive charge, whereas HbA0 retains its positive charges. Therefore, HbA1c is first eluted from the Hemoglobin A1c Analytical Column. Then under the action of the high ionic strength eluent (Hemoglobin A1c Eluent B), HbA0 is eluted. Through this process, a chromatogram is developed for the purpose of calculating the ratio between the area of the HbA1c peak in relation to the total Hb area.

Components and Structure

The column is composed of the ion-exchange stationary phase, the

column tubular that is made of stainless steel, the sealing caps and accessories. The accessories include filters.

Storage Conditions and Stability

- The product may be used up to the expiry date indicated on the label if stored under sealed condition at 2°C to 8°C.
- When the column is stored without being sealed, the material will dry out and become unusable.
- When not used, store the column sealed with the sealing cap under refrigerated condition at 2°C to 8°C. Avoid freezing.
- Use the product at the relative humidity of 20%RH to 85%RH. The working temperature range of the product is consistent with that of its applicable instruments.
- The production date is indicated on the product package or label.
- Refer to the table below for the maximum number of tests for each model.

Table 2. Maximum number of tests for each model

Model	Number of tests
H-3	3,000 tests
H-6	6,000 tests
H-10	10,000 tests

Applicable Instruments

This product applies to H-120 and H-120P Automated Glycohemoglobin Analyzers manufactured by Shenzhen Mindray Bio-Medical Electronics Co., Ltd..

Material Required but not Provided

The following materials are required but not provided with the product: Mindray manufactured measurement instruments and matched reagents, as well as general laboratory devices.

Sample Collection and Preparation

- Fresh human venous whole blood samples anticoagulated with EDTA K₂, EDTA K₃, lithium heparin, sodium citrate, potassium oxalate or sodium fluoride.
- Immediately analyze the whole blood samples once they have been collected. If this is not feasible, the samples may be stored at 18°C to 25°C for one day, or at 2°C to 8°C for seven days before analysis.
- For more information regarding the requirements for sample collection and preparation, refer to the Operator's Manual of the

applicable instruments.

Test Methods

- Before every day test, follow the instructions in the Operator's Manual of the applicable instrument to perform the quality control program.
- For the test method, refer to the Operator's Manual of the applicable instruments.

Installation and Replacement

- Replace the column when it reaches its shelf life or is broken.
- When replacing the column, wear proper protective clothing (gloves, glasses, masks, etc.).
- Confirm the integrity of the package before use. Do not use the product if the package is damaged. The results cannot be assured when the reagents are stored at inappropriate conditions.
- The product is supplied with a snap-fit pre-installed. Use the column with the snap-fit together. It is unnecessary to detach the product from the snap-fit.
- Install the column in the correct direction; otherwise, the column cannot be installed into the column compartment of the analyzer.
- If any error occurs during the installation process, follow the software instruction to remove the error and try again. If the error remains, contact our customer service department.

Follow the instructions below to replace / install the column:

1. Open the package of the column and take out the column. Wait for the column to return to its service temperature.
2. Make sure the analyzer is in standby state, and open the front cover of the analyzer.
3. Flip open the column compartment door, and remove the old column (if there is one) with its snap-fit together.
4. Remove the two sealing caps from both ends of the new column.
5. Installed the column together with its snap-fit into the column compartment, and then close the compartment door.
6. Close the front cover the analyzer.
7. Follow the instructions in the analyzer's Operator's Manual and the software guidance to register the new column. Then the analyzer will automatically prime the tubing of the column.
8. After column installation / replacement, calibrate the analyzer before testing.

For detailed information, refer to the Operator's Manual of the applicable instrument.

Cut-off Value or Reference Interval

- DCCT/UKPDS recommends 4.0% - 6.0% (20 mmol/mol-42mmol/

mol) as the reference interval for HbA1c in venous whole blood samples from non-diabetic healthy adults.

- It is recommended that each laboratory to establish its own reference interval based on its own patient population.

Result Elaboration

- Any variations in the red blood cell lifespan may affect the patient's HbA1c concentration level. Take care when interpreting the test results of such patients.
- In case of insufficient sample volume or very low hemoglobin concentration, the analyzer may report the message "Chromatogram area too small". Take care when reporting such results.
- Test results should be interpreted based on the individual's medical history, diagnosis, and the results of other tests.

Limitation

- A low concentration of hemoglobin in a sample can result in a smaller total chromatogram area.
- When tested at or below the specified concentrations, the following substances do not present significant interference to the results (interference is lower than 7%).

Table 3. Possible interfering substances

Substances	Concentration
Free bilirubin or conjugated bilirubin	20 mg/dL
Triglyceride	3300 mg/dL
Glucose	1000 mg/dL
Aspirin (medicine)	65 mg/dL
Sodium cyanate (among patients with kidney diseases)	25 mg/dL
Ascorbic acid	200 mg/dL
Acetaldehyde (produced due to alcohol)	25 mg/dL

Performance Specifications

Appearance

The column body is intact and undamaged, with no signs of leakage.

Column pressure

The pressure of the column should not be higher than 15 MPa.

Accuracy

Test certified reference material or manufacturer-provided reference material on an applicable analyzer, the relative deviation between the test result and the reference value should be within the range of $\pm 5.0\%$.

Repeatability

Test samples with HbA1c concentration at 4.00%- 6.50% (NGSP) or 20.20 mmol/mol – 47.50 mmol/mol in (IFCC) on an applicable instrument for multiple times, the coefficient variation (CV) between the results should not exceed 1.0%.

Linearity

When the sample's concentration is within the range of 3.00%- 20.10% (NGSP), the linear correlation coefficient (r) is no less than 0.9900.

Warnings and Precautions

- For *in-vitro* diagnostic use only. For laboratory professional use.
- Read the instruction for use carefully before use, and use the product according to the instructions.
- This column should only be used together with the Auto Glycohemoglobin Analyzers and applicable reagents manufactured by Shenzhen Mindray Bio-Medical Electronics Co., Ltd.. For more information of the applicability, refer to the Operator's Manual of the analyzers.
- Disposal of waste liquid and materials should be in accordance with local guidelines.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the country in which the user and/or the patient is established.
- This instruction for use is intended to be read by clinical laboratory professionals to:
 - a. learn about the product;
 - b. perform daily operating task;
 - c. perform maintenance and troubleshooting.

References

American Diabetes Association. Standards of Medical Care for Patients with Diabetes Mellitus. Diabetes Care 2002; 25 (Suppl. 1), S33-49.

Symbols

Batch code

Use-by date

Temperature limit

Consult instructions
for use

Manufacturer

In vitro diagnostic
medical device

Serial number

Unique device
identifier

European
Conformity

Authorized
representative
in the European
Community

Approval Date of the Instructions for Use

2025-09

The summary of safety and performance is available from the eudamed database:

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