(EN) Hemoglobin A1c Control**[PRODUCT]**

Hemoglobin A1c Control

[PACKAGE SIZE]

REF	Package size
105-007859-00	CRL-1: 2mLx1, CRL-2: 2mLx1
105-007862-00	CRL-1: 2mLx3, CRL-2: 2mLx3

[INTENDED USE]

The product is used for quality control over the test on HbA1c in whole blood to monitor and assess the precision of the test result.

[PRINCIPLE]

The component of this product is similar to the test samples. Perform statistical analysis on the test result or compare the test result with the reference value (reference range) to determine whether it meets requirements.

[ACTIVE INGREDIENTS]

Freeze-dried powder of whole blood matrix of people.

* The content concentrations of products of different batches differ slightly.

* For products of different batches, the range of reference values is different. Therefore, a separate reference value table is provided for each batch of products.

[STORAGE AND EXPIRATION DATE]

Products that are not decapped should be stored in the temperature 2°C to 8°C for 1 year. After the products are re-dissolved and decapped, the products should be stored in the temperature 2°C to 8°C for 14 days.

[APPLICABLE INSTRUMENT]

This product applies to H50/H50P/H-120/H-120P Automated Glycohemoglobin Analyzer produced by Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

[SAMPLE REQUIREMENT]

Not applicable.

[MATERIALS REQUIRED BUT NOT PROVIDED]

The following materials are required but not provided with the product: Mindray-manufactured measurement instruments and matched reagents, calibrators, as well as general laboratory devices.

For the required reagents and calibrator, refer to the table below.

For H50 and H50P models:

Reagents and calibrator	REF
Hemolysis Solution	105-007856-00
Eluent A	105-007852-00
Eluent B	105-007853-00 105-007854-00 105-007855-00
Hemoglobin A1c Calibrator	105-007857-00 105-007858-00

For H-120 and H-120P models:

Reagents and calibrator	REF
Hemoglobin A1c Hemolysis Solution	105-048455-00
Hemoglobin A1c Eluent A	105-048453-00
Hemoglobin A1c Eluent B	105-048454-00 105-048464-00
Hemoglobin A1c Calibrator	105-007857-00 105-007858-00

[INSTRUCTION FOR USE]

- Take the product from the environment of 2°C to 8°C and wait until the product reaches the using temperature(10°C to 32°C);
 - Gently knock the bottle cap to guarantee that the freeze-dried sample completely falls onto the bottom. Open the bottle cap and rubber stopper with caution to avoid loss of the content.
 - Accurately measure 2mL Hemolysis Solution/pure water*, slowly fill it into the bottle, and put the rubber stopper on the bottle. Gently turn over the small bottle for several seconds and place the product for 30 minutes in a dark place;
 - Gently turn over the bottom to guarantee that the content is completely dissolved and avoid foam.
 - After the content is re-dissolved, select an appropriate centrifuge tube or sample cup to hold the content and mark it for later use.
- * For H50 and H50P models, use Hemolysis Solution or pure water; for H-120 and H-120P models, use only pure water.

[CUT-OFF VALUE/REFERENCE INTERVAL]

Not applicable.

[RESULT ELABORATION]

The mean of the control results measured by the measurement system should fall into the reference range. If it does not, check the measurement system for the control's expiration date, storage conditions; calibration method, and the analyzer's performance and working status. Laboratories may also establish their own reference values or ranges for the measurement systems in accordance with their QC rules.

[LIMITATIONS]

Not applicable.

[PRODUCT SPECIFICATIONS]


The reference value and reference range of the quality control materials are determined according to the standard valuation procedure. The reference value is the average value obtained through repeated test on

the quality control materials. For the reference value, reference range and test method, see the reference value table.

[PRECAUTIONS]

- For professional use in vitro diagnosis only.
- The control values are lot-specific. Please confirm before use.
- Perform scheduled maintenance and standard operation including calibration and analysis to assure the analyzer's performance.
- The control was tested with a method confirmed by a national management authority and is shown to be negative for antibody to human immunodeficiency virus (HIV), hepatitis B surface antigen (HBsAg), antibody to hepatitis C virus (HCV) and antibody to Treponema pallidum (TP). However, as no testing method can rule out the potential risk of infection with absolute certainty, this material should be handled as a patient sample to avoid biological risk.
- Take the necessary precautions for the use of reagents. Do not swallow. Avoid contact with skin and mucous membranes. If you accidentally take the reagent into your mouth, or the reagents accidentally spill on your skin or into your eyes, wash them off with plenty of water and go see a doctor if necessary.
- Besides the reference values and ranges listed in the reference value sheet for Mindray's measurement system, laboratories may also establish their own QC rules in accordance with their specific conditions, and determine the reference values or ranges accordingly.
- Disposal of waste liquid and materials should be in accordance with local guidelines.
- Do not inhale and/or ingest. In case of inhalation or ingestion, seek medical advice immediately.
- Avoid skin contact. In case of skin contact, flush the affected area with plenty of water immediately.
- Avoid eye contact. In case of eye contact, flush the affected area with plenty of water and seek medical advice immediately.
- The Material Safety Data Sheet (SDS) is available upon request
- All identified risks have been reduced as far as possible by generally acknowledged state of art, and the overall residual risk is acceptable.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

[GRAPHICAL SYMBOLS]

		
Batch code	Use-by date	Temperature limit
		
<i>In vitro</i> diagnostic medical device	Manufacturer	Consult Instruction for Use
		
Authorized representative in the European Community	Catalogue number	Biological risks
		
European Conformity	Unique device identifier	Control

[REFERENCES]

Biosafety in Microbiological and Biomedical Laboratories; U.S. Department of Health and Human Services; US Government Printing Office; Washington: 2007

[COMPANY CONTACT]

Manufacturer Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Address Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

Website www.mindray.com
E-mail address service@mindray.com

Tel +86 755 81888998

Fax +86 755 26582680

Authorized Representative in the European Community

EC-Representative Address Shanghai International Holding Corp. GmbH(Europe)
Tel Eiffestraße 80, 20537 Hamburg, Germany 0049-40-2513175
Fax 0049-40-255726

© 2021-2025 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved.

[APPROVAL DATE OF THE OPERATOR'S MANUAL]
2025-07

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REF	Paket boyutu
105-007859-00	CRL-1: 2mLx1, CRL-2: 2mLx1
105-007862-00	CRL-1: 2mLx3, CRL-2: 2mLx3

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Fabricante Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Representante da CE: Shanghai International Holding Corp. GmbH(Europe)
© 2012-2025 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Todos os direitos reservados.