1.5 L × 1: 105-048453-00

Instructions for Use

Hemoglobin A1c Eluent A

Find IFUs in more languages on our website:
<https://www.mindray.com/en/resources-center/ivd-edocs>

COMPANY CONTACT / KONTAKT MIT DEM UNTERNEHMEN/CONTACTO CON LA EMPRESA / COORDONNEES DE LA SOCIETE/INFORMAZIONI RELATIVE ALLA SOCIETÀ / ŞİRKET İLETİŞİM / CONTATO DA EMPRESA / КОНТАКТНАЯ ИНФОРМАЦИЯ

	Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
	Mindray Building, Keji 12 th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China
	www.mindray.com
	service@mindray.com
Tel /Tél/Тел:	+86 755 81888998
Fax/Télécopie/ Faks/Факс:	+86 755 26582680

	Shanghai International Holding Corp. GmbH(Europe)
	Eiffestraße 80, Hamburg 20537, Germany
Tel /Tél/Тел:	0049-40-2513175
Fax/Télécopie/ Faks/Факс:	0049-40-255726

Revision bar indicates update to previous version / Das Revisionsbarsymbol weist auf eine Aktualisierung der vorigen Version / La barra de revisión indica una actualización a la versión anterior / La barre de révision indique la mise à jour vers la version précédente / La barra di revisione indica l'aggiornamento alla versione precedente / Değişiklik çubuğu, önceki sürüme yapılan güncellemeyi gösterir / A barra de revisão indica a atualização para a versão anterior / Строка обновления указывает на обновление предыдущей версии

Hemoglobin A1c Eluent A

Packing Specifications

Table 1. Packing specification

Model	Packing specification
H-12A	1.5 L × 1

Intended Use

This product applies to Mindray H-120/H-120P Automated Glycohemoglobin Analyzer. The product participates in the measurement of HbA1c in blood samples.

Test Principle

Hemoglobin A1c Eluent A (hereafter Eluent A) and Hemoglobin A1c Eluent B (hereafter Eluent B) work together in HbA1c analysis.

In the presence of the low ionic strength eluent (Eluent A), HbA1c hardly carries any positive charge, whereas HbA0 retains its positive charges. Therefore, HbA1c is first eluted from the Hemoglobin A1c Analytical Column. Then under the action of the high ionic strength eluent (Eluent B), HbA0 is eluted. In this way, the analyzer builds up a complete chromatogram.

Major Components

Citric acid salt..... ≤ 5.9 g/L
 Isothiazolinones antimicrobial agent..... ≤ 0.099 g/L

Storage Conditions and Stability

- The product may be used up to the expiry date indicated on the label if stored unopened at 2°C to 30°C, in a well-ventilated area, away from corrosive gas and direct sunlight.
- The working temperature range of the product is consistent with that of its applicable instruments. The open-vial validity is 100 days.
- For the production date and expiry date, see the package or the label.

Applicable Instruments

This product applies to H-120 and H-120P Automated Glycohemoglobin Analyzers manufactured by Shenzhen Mindray Bio-Medical Electronics Co., Ltd..

Material Required but not Provided

The following materials are required but not provided with the product: Mindray manufactured measurement instruments and matched reagents, as well as general laboratory devices.

Sample Collection and Preparation

- Fresh human venous whole blood samples anticoagulated with EDTA K₂, EDTA K₃, lithium heparin, sodium citrate, potassium oxalate or sodium fluoride.
- Immediately analyze the whole blood samples once they have been collected. If this is not feasible, the samples may be stored at 18°C to 25°C for one day, or at 2°C

to 8°C for seven days before analysis.

- For more information regarding the requirements for sample collection and preparation, refer to the Operator's Manual of the applicable instruments.

Reagent Preparation

The product is a ready-for-use reagent.

Test Procedure

1. Open the external packing of the product, and restore the product to usage temperature.
2. When the instrument is not in use, open the front cover, remove the old reagent container (if there is), and install the new container to the correct reagent compartment.
3. Close the instrument front cover.
4. Follow the instructions in the Operator's Manual of the applicable instrument to replace the reagent in the instrument, and dispose the old reagent container properly.

For detailed information, refer to the Operator's Manual of the applicable instruments.

Cut-off Value or Reference Intervals

- DCCT/UKPDS recommends 4.0% - 6.0% (20 mmol/mol - 42mmol/mol) as the reference interval for HbA1c in venous whole blood samples from non-diabetic healthy adults.
- It is recommended that each laboratory to establish its own reference interval based on its own patient population.

Result Elaboration

- Any variations in the red blood cell lifespan may affect the patient's HbA1c concentration level. Take care when interpreting the test results of such patients.
- In case of insufficient sample volume or very low hemoglobin concentration, the analyzer may report the message "Chromatogram area too small". Take care when reporting such results.
- Test results should be interpreted based on the individual's medical history, diagnosis, and the results of other tests.

Limitation

- A low concentration of hemoglobin in a sample can result in a smaller total chromatogram area.
- When tested at or below the specified concentrations, the following substances do not present significant interference to the results (interference is lower than 7%).

Table 2. Possible interfering substances

Substances	Concentration
Free bilirubin or conjugated bilirubin	20 mg/dL
Triglyceride	3,300 mg/dL

Substances	Concentration
Glucose	1,000 mg/dL
Aspirin (medicine)	65 mg/dL
Sodium cyanate (among patients with kidney diseases)	25 mg/dL
Ascorbic acid	200 mg/dL
Acetaldehyde (produced due to alcohol)	25 mg/dL

Performance Characteristics

Appearance

The product is transparent liquid without deposits, suspended grains or flocks.

Accuracy

When testing certified reference material or manufacturer-provided reference material on an applicable analyzer, the relative deviation between the test result and the reference value should be within the range of $\pm 5.0\%$.

Repeatability

When testing samples or HbA1c controls with HbA1c concentration at 4.00%- 6.50% (NGSP) or 20.20 mmol/mol - 47.50 mmol/mol in (IFCC); and at 9.00% - 12.00% (NGSP) or 74.87 mmol/mol - 107.66 mmol/mol (IFCC) on an applicable instrument for multiple times, the coefficient variation (CV) between the results should not exceed 1.0%.

Linearity

When the sample's concentration is within the range of 3.00%- 20.10% (NGSP), the linear correlation coefficient (r) is no less than 0.9900.

Warnings and Precautions

- For in-vitro diagnostic use only. For laboratory professional use.
- Read the instruction for use carefully before use, and use the product according to the instructions.
- The product is transparent liquid without deposits, suspended grains or flocks. If it gets polluted or affected by other factors and becomes abnormal, stop using it and replace it with a normal one.
- Do not mix reagents of different batches.
- If the reagent is frozen, thaw the frozen reagent at room temperature and ensure it is thoroughly mixed before using it.
- The sample analysis may be impacted by various factors, including the use of expired or ineffective reagents, mixing reagents with third-party products, combining residue from an old container with newly-opened reagents, and using reagents outside of their specified conditions.
- Do not take the product into mouth. If you accidentally ingest the product, seek medical treatment immediately.

- Avoid exposure to skin. If you accidentally spill the product on your skin, wash it off with plenty of water immediately.
- Avoid exposure to eyes. If you accidentally spill the product into your eyes, wash it off with plenty of water immediately, and seek medical treatment if necessary.
- Confirm the integrity of the package before use. Do not use the reagents with damaged packages. The results cannot be assured when the reagents are stored at inappropriate conditions.
- If the reagents are opened unintentionally before use, they shall be used as soon as possible.
- Disposal of waste liquid and materials should be in accordance with local guidelines.
- The Material Safety Data Sheet (SDS) is available upon request.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the country in which the user and/or the patient is established.
- This product contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:


	
Warning	
Hazard statement(s)	
H317	May cause an allergic skin reaction.
Precautionary statement(s)	
EUH210	Safety data sheet available on request.
Precautionary statement(s) Prevention	
P280	Wear protective gloves and protective clothing.
P261	Avoid breathing mist/vapours/spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
Precautionary statement(s) Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice/attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Precautionary statement(s) Disposal	

P501	Dispose of contents/container to authorised hazardous or special waste collection point in accordance with any local regulation.
------	--

References

American Diabetes Association. Standards of Medical Care for Patients with Diabetes Mellitus. Diabetes Care 2002; 25 (Suppl. 1), S33-49.

Symbols

		
Batch code	Use-by date	Temperature limit
		
Authorized representative in the European Community	Manufacturer	<i>In vitro</i> diagnostic medical device
		
Consult instructions for use	Keep away from sunlight (Indicates a medical device that needs protection from light source)	Catalogue number
		
European Conformity	Unique device identifier	

Approval Date of the Instructions for Use

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