MR Urine Quality Control

Order Information

Cat. No.	Package size
105-001691-00	High: 1 × 100 mL, Low: 1 × 100 mL

Intended Purpose

MR Urine Quality Control is an electrolyte reference urine control used to monitor the urine mode performance of ion-selective electrode on ISE Modules in Mindray BS series Chemistry Analyzers. Daily use of this urine control provides quality control data for confirming the precision and accuracy of urine sample tests.

Summary

MR Urine Quality Control applies to Mindray BS series chemistry analyzers. Please use the MR Urine Quality Control as instructed by the operator's manual of the analyzer.

Components

MR Urine Quality Control is mainly composed of sodium chloride (<1.0%), potassium chloride (<1.0%) and formaldehyde (<0.1%).

Storage and stability

MR Urine Quality Control is stable up to expiration date indicated on the label when stored in unopened vial at 2-10°C. The control vial should be kept away from light. After each use, promptly replace the cap and ensure the vial is capped tightly.

Once opened, this product is stable for 8 weeks when stored tightly capped at 2-10°C.

Usage

Use this product on the chemistry analyzer as instructed in the instrument manual.

Materials required but not provided

1. Mindray Sodium electrode, potassium electrode, chloride electrode and reference electrode for ISE module of chemistry analyzers; Mindray MR Urine Standard. The catalogue number of electrodes and Urine Standard can be obtained from the corresponding Mindray electrode and Urine Standard instruction for use.
2. Mindray BS series chemistry analyzers, and General laboratory equipment.

Performance Characteristics

■ Assay Values

The control assay values (the target value and range) determined by Mindray standard transfer procedure and routine method are listed in the value sheet. The target value was obtained from Mindray measurement system and the

range was calculated as the target values ± 3 standard deviations. The control values are specific for different lot numbers, so please check the lot number before use.

Mindray systems (BS series chemistry analyzers/Mindray electrodes) are calibrated by Mindray Urine Standard, then test the control on Mindray systems, the assayed values should be within the defined range.

■ Repeatability

Items	Repeatability
Na ⁺	CV \leq 3.5%
K ⁺	CV \leq 3.5%
Cl ⁻	CV \leq 3.5%

Result interpretation

If the control result falls outside the range, the measurement system should be checked. For instance, position of sample placed on the analyzer; expiration date or storage condition of the Urine Standard and control; the parameter setting in software; or performance of the analyzer; or the calibration process. Each laboratory should establish its own internal quality control scheme and procedure for corrective action.

Limitations

The product only applies to the ISE module of Mindray BS series chemistry analyzers. In addition, the warranted results only come from processing tests strictly as instructed in the instrument manual and storing the product correctly.

Warnings and precaution

1. For in vitro diagnostic use only. For laboratory professional use.
2. Please take the necessary precautions for handling all laboratory reagents.
3. Please confirm the integrity of the package before use. Do not use the controls with damaged packages.
4. If unintentionally opened before used, store the controls tightly capped at 2-10°C and protected from light, and the stability is equally to in-use stability.
5. Do not mix controls with different lots and bottles.
Do not use the control beyond the expiration date and the in-use date.
Do not mix new opened control with in-use control.
Avoid the formation of foam.
6. If any appearance change occurs, like a cloudy or turbid appearance shows up, do not use this vial of the product any more.
7. Do not swallow. Avoid contact with skin and mucous membranes.
8. When the control accidentally enters the eyes and mouth, or contact with

the skin, immediately wash with plenty of water. If necessary, visit the doctor for further medical treatment.

9. Safety data sheet is available for professional user on request.
10. Disposal of all waste material should be in accordance with local guidelines.
11. All human material should be considered potentially infectious¹.
12. All identified risks have been reduced as far as possible without adversely affecting the benefit-risk ratio, and the overall residual risk is acceptable.
13. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

References

1. Biosafety in Microbiological and Biomedical Laboratories; U.S. Department of Health and Human Services; US Government Printing Office; Washington:2007.

Graphical symbols

In Vitro Diagnostic
medical device

Unique device
identifier

European
Conformity

Consult Instructions
For use

Use-by
date

Authorized representative in
the European Community

Batch Code

Temperature
limit

Manufacturer

Catalogue
number

Control

Keep away from sunlight
Indicates a medical device that needs
protection from light sources

© 2022-2026 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

E-mail Address: service@mindray.com

Website: www.mindray.com

Tel: +86-755-81888998; **Fax:** +86-755-26582680

EC-Representative: Shanghai International Holding Corp. GmbH(Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

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© 2022-2026 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Tous droits réservés
Fabricant : Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Représentant en Europe : Shanghai International Holding Corp. GmbH (Europe)