# T3

Total Triiodothyronine (CLIA)

INSTRUCTIONS FOR USE

**mindray**

## T3

### Total Triiodothyronine (CLIA)

#### Order Information

Catalog No.	Package Size
105-004210-00	2×50 tests
105-004237-00	2×100 tests

#### Intended Use

The CL-series T3 assay is a Chemiluminescent Immunoassay (CLIA) for the quantitative determination of total triiodothyronine (T3) in human serum or plasma.

Measurements obtained by this device are used as an aid in the diagnosis and treatment of thyroid disease.

#### Summary

3,5,3' Triiodothyronine (T3) is a thyroid hormone with a molecular weight of 651 daltons and a half-life in serum of 1.5 days.<sup>1</sup> It is an iodide of tyrosine residues, synthesized and secreted by thyroid gland. The degree of their biological effect, are controlled by the hypothalamo-hypophyseal-thyroid axis. Thyrotropin releasing hormone (TRH) is released by the hypothalamus in response to circulating concentration of free thyroid hormones. TRH produce a marked effect on the adenohypophysis and initiates an intracellular cascade of events that result in the production and release of thyroid-stimulating hormone (hTSH). The target organ for hTSH is the thyroid gland, which synthesis and secret T4 and T3 as the response of hTSH stimulation. Only 0.2-0.4% of the total T3 and T4 is present in solution as unbound or free, as almost all is bound by serum proteins, which is not available to incite biological activity. Free and bound thyroid hormone can transform and maintain dynamic equilibrium.<sup>3</sup> It is only the free thyroid hormone that is readily available to bind its receptor, and stimulate a response from the target organ or tissues.<sup>5,6</sup> Only a small proportion of total T3 in serum comes from direct secretion from the thyroid gland. The remaining fraction of total T3 is derived from enzymatic monodeiodination of T4 to T3 by the peripheral tissues. They are transferred into cell nucleus through energy and carrier and combined with intranuclear acceptor.<sup>4</sup>

Serum T3 determination can be a valuable component of a thyroid screening panel in the diagnosis of thyroid disorders. In some serious and chronicity thyroid illness, the concentration of free T4 increase or decrease with high TSH and low T3.<sup>2</sup> Normal concentration of free T4 and elevated concentration of T3 will give rise to hyperthyroidism.<sup>7</sup> Depression concentration of T3 is caused

by primary hypothyroidism or secondary cases hypothyroidism caused by loss of hypothalamo or hypophyseal functions, such as Hashimoto's thyroiditis.<sup>5B</sup>

### Assay Principle

The CL-series T3 assay is a competitive binding assay to determine the level of total T3.

In the first step, sample, stripping agent and monoclonal anti-T3 antibody-alkaline phosphatase conjugate are added into a reaction cuvette. After incubation, the stripping agent dissociates T3 from the binding proteins.

In the second step, the streptavidin coated microparticle and biotinylated T3 are added into the reaction cuvette. T3 in the sample competes with the biotinylated T3 for anti-T3 alkaline phosphatase conjugate. The biotinylated T3-antibody complexes are bound to the streptavidin coated microparticle. Microparticle is magnetically captured while other unbound substances are removed by washing.

In the third step, the substrate solution is added to the reaction cuvette. It is catalyzed by anti-T3 antibody-alkaline phosphatase conjugate in the immunocomplex retained on the microparticle. The resulting chemiluminescent reaction is measured as relative light units (RLUs) by a photomultiplier built into the system. The amount of total T3 present in the sample is inversely proportional to RLUs generated during the reaction. The total T3 concentration can be determined via a calibration curve.

### Reagent Components

Ra	Paramagnetic microparticles coated with streptavidin. Minimum concentration: 0.1 g/L. HEPES <sup>a)</sup> buffer: 50 mmol/L. Preservatives: 0.05% ProClin 300 and 0.09% sodium azide.
Rb	Monoclonal anti-T3 antibody-alkaline phosphatase conjugate. Minimum concentration: 45 µg/L. MES <sup>b)</sup> buffer: 50 mmol/L. Preservatives: 0.05% ProClin 300 and 0.09% sodium azide.
Rc	Biotinylated T3. Minimum concentration: 1.0 µg/L. HEPES buffer: 50 mmol/L. Preservatives: 0.05% ProClin 300 and 0.09% sodium azide.

Rd	<p>8-Anilino-1-Napthalenesulfonic Acid (ANS).          Minimum concentration: 1 g/L.          HEPES buffer: 50 mmol/L.          Preservatives: 0.05% ProClin 300 and 0.09% sodium azide.</p>
----	--

- a) HEPES=2-[4-(2-Hydroxyethyl)-1-piperazinyl]ethanesulfonic acid
- b) MES=2-(N-morpholino)ethanesulfonic acid

### Storage and Stability

The unopened T3 (CLIA) reagent kit is stable up to the stated expiration date when stored at 2-8°C.

The T3 (CLIA) reagent kit can be stored onboard and used for a maximum of 56 days after opening at 2-8°C.

### Reagent Preparation

The reagents in the kit have been assembled into an intact unit that cannot be separated.

### Materials Required but not Provided

Mindray CL-series Chemiluminescence Immunoassay Analyzer

Cat.No.105-004279-00: Total T3 Calibrators, 1x2.0 mL for each of calibrator C0, C1 and C2.

Cat.No.105-007383-00: Thyroid Function Multi Control (L), 1x5.0 mL

Cat.No.105-007384-00: Thyroid Function Multi Control (H), 1x5.0 mL

Cat.No.105-007379-00: Thyroid Function Multi Control (L), 3x5.0 mL

Cat.No.105-007380-00: Thyroid Function Multi Control (H), 3x5.0 mL

Cat.No.105-007371-00: Thyroid Function Multi Control (L), 6x5.0 mL

Cat.No.105-007372-00: Thyroid Function Multi Control (H), 6x5.0 mL

Cat.No.105-007375-00: Thyroid Function Multi Control (L), 12x5.0 mL

Cat.No.105-007376-00: Thyroid Function Multi Control (H), 12x5.0 mL

Cat.No.105-004552-00: Wash Buffer, 1x10 L

Cat.No.105-009044-00: Substrate Solution, 4x75 mL

Cat.No.105-004274-00: Substrate Solution, 4x115mL

Reaction Cuvette.

**Applicable Instrument**

Mindray CL-series Chemiluminescence Immunoassay Analyzer

**Specimen Collection and Preparation****Specimen Types**

- Human serum or plasma collected in sodium heparin and lithium heparin are recommended for this assay.
- Blood collection tubes from various manufacturers may contain additives which could affect the test results in some cases. Not all available tubes on the market were tested by Mindray. Each laboratory should determine the acceptability of different blood collection tubes and serum/plasma separation products. Variations in these products may exist among manufacturers and from lot-to-lot.
- When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

**Specimen Conditions**

- Do not use:
  - heat-inactivated specimens
  - grossly hemolyzed specimens
  - specimens with obvious microbial contamination
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.

**Preparation for Analysis**

- Follow blood collection tube manufacturer's recommendations for centrifugation. Centrifuge the specimens after clot formation is complete. Ensure residual fibrin and cellular matter has been removed prior to analysis.
- For optimal results, inspect all samples for bubbles. Remove bubbles with a pipette tip prior to analysis. Specimens must be mixed thoroughly after thawing. Thawed samples should be centrifuged prior to use.
- If the sample was covered with lipid layer after centrifugation, the sample should be transferred to a clean tube and centrifuged before testing. Do not transfer the lipid layer. Handle carefully to prevent cross contamination.

**Specimen Storage**

- Specimens should be tested timely after sample collection. If the assay cannot be completed within 8 hours, refrigerate the samples at 2 to 8°C. If testing will be delayed for more than 72

hours, specimens should be frozen at  $-20^{\circ}\text{C}$  or colder. The specimens can be stored at  $-20^{\circ}\text{C}$  for up to 30 days.

- The number of repeated freezing and thawing shall not exceed 5 times.

### Assay Procedure

For optimal performance of this assay, operators should read the related system operation manual carefully, to get sufficient information such as operation instructions, sample preservation and management, safety precaution, and maintenance. Prepare all required materials for the assay as well.

Before loading the T3 (CLIA) reagent kit on the machine for the first time, unopened reagent bottle should be inverted gently for at least 30 times to resuspend the microparticles that have settled during shipment or storage. Visually inspect the bottle to ensure the microparticles have been resuspended. If the microparticles remain adhered to the bottle, continue inverting until the microparticles have been completely resuspended. If the microparticles cannot be resuspended, it is recommended not to use this bottle of reagent. Contact Mindray Customer Service for help. Do not invert opened reagent bottle.

This assay requires 20  $\mu\text{L}$  of sample for a single test. This volume does not include the dead volume of the sample container. Additional volume is required when performing additional tests from the same sample. Operators should refer to the system operation manual and specific requirement of the assay to determine the minimum sample volume.

### Calibration

CL-series T3 (CLIA) has been standardized against the IDMS (isotope dilution mass spectrometry) reference measurement procedure.

The specific information of master calibration curve of T3 (CLIA) reagent kit is stored in the two-dimensional barcode attached in the reagent pack. It's used together with calibrators for the calibration of the specific reagent lot. When performing calibration, scan the information of master calibration curve from the barcode into the system first, and then use the calibrators at three levels. Valid calibration curve is required before any T3 test. Recalibration is recommended every 4 weeks, or when a new reagent lot is used, or the quality controls are out of specified range. For detailed instruction of calibration, refer to the system operation manual.

### Quality Control

It is recommended to run quality controls once every 24 hours if the tests are in use, or after every calibration. The quality control frequency should be adapted to each laboratory's individual

requirements. The recommended two levels of quality controls for this assay are Mindray Thyroid Function Multi Control (L) and Thyroid Function Multi Control (H). In addition, other suitable control material can be used.

Quality control results should be within the acceptable ranges. If a control is out of its specified range, the associated test results are invalid and the samples must be retested. Recalibration may be required. Examine the assay system referring to the system operation manual. If the quality control results are still out of the specified range, please contact Mindray Customer Service for help.

### Calculation

The analyzer automatically calculates the analyte concentration of each sample on the master calibration curve read from the barcode, and the RLUs generated from the three level calibrators of defined concentration values. The results are shown in the unit of ng/mL or nmol/L (selectable).

Conversion factors:  $\text{ng/mL} \times 1.536 = \text{nmol/L}$

$\text{nmol/L} \times 0.651 = \text{ng/mL}$

### Dilution

Samples with T3 concentrations above the upper limit can be diluted with Mindray Sample Diluent. The recommended dilution is 1:2 (either automatically by the analyzer or manually). The concentration of the diluted sample must be  $> 2 \text{ ng/mL}$ . After manual dilution, multiply the result by the dilution factor. If use automated dilution on the analyzer, the system automatically multiplies the result by the dilution factor when calculating the sample concentration.

### Expected values

A normal range of 0.58 ng/mL to 1.62 ng/mL (central 95% interval) was obtained by testing serum specimens from 530 non-pregnant individuals determined as normal by MINDRAY TSH and Free T4 assays.

Representative T3 ranges during normal pregnancy are summarized below.

Pregnancy Gestation	N	Central 95% Range
First Trimester	156	0.71 -1.60 ng/mL
Second Trimester	160	0.85 -1.95 ng/mL
Third Trimester	154	0.72 -1.85 ng/mL

It is highly recommended that each laboratory establishes its own normal range, which may be unique to the population it serves depending upon geography, race, sex, age, disease, dietary, or environmental factors.

**Limitation**

The upper limit of this assay is 8.0 ng/mL. A specimen with a total T3 concentration lower than the upper limit can be quantitatively determined, while specimen with a concentration higher than the upper limit will be reported as >8.0 ng/mL or retested by diluting the samples with Mindray Sample Diluent.

The concentration of total T3 in a given specimen, determined with assays from different manufacturers, can vary due to differences in assay methods, calibration, and reagent specificity. The assay results should be used in conjunction with other data, such as symptoms, results of other tests, clinical history, etc.

Specimens of hyperlipidemia, icterus, hemolysis or high total protein may result in incorrect results.

Heterophilic antibodies in samples can react with immunoglobulin in kit and interfere with the immunoassay results.<sup>9</sup> Patients who have been regularly exposed to animals or have been treated with immunoglobulins or immunoglobulin fragments may produce such antibodies, resulting in abnormal results. For example, specimens from individuals who have been exposed to mouse monoclonal antibodies may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values with assay kits employing mouse monoclonal antibodies.<sup>10,11</sup> There are anti-interference components in Mindray T3 (CLIA) reagent kit, which can effectively reduce the influence caused by HAMA in samples, but there may still be interference problems in a few samples. Clinical examinations, medical history, or other relevant information are needed to determine the diagnosis of the patient.

In rare cases, interference due to extremely high titers of antibodies to analyte-specific antibodies or streptavidin can occur.

**Performance Characteristics****Lower limits of measurement**

Limit of Blank, Limit of Detection and Limit of Quantitation

Limit of Blank = 0.1 ng/mL

Limit of Detection = 0.2 ng/mL

Limit of Quantitation = 0.36 ng/mL

The Limit of Blank, Limit of Detection and Limit of Quantitation were determined in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP17-A2 requirements.

The Limit of Blank is the 95th percentile value from  $n \geq 60$  measurements of analyte-free samples over several independent series. The Limit of Blank corresponds to the concentration below which analyte-free samples are found with a probability of 95 %.

The Limit of Detection is determined based on the Limit of Blank and the standard deviation of low concentration samples. The Limit of Detection corresponds to the lowest analyte concentration which can be detected (value above the Limit of Blank with a probability of 95 %).

The Limit of Quantitation is the lowest analyte concentration that can be reproducibly measured with an intermediate precision CV of  $\leq 20$  %.

### Measuring Range

0.2-8.0 ng/mL (defined by the Limit of Detection and the maximum of the master calibration curve). Values below the Limit of Detection are reported as  $< 0.2$  ng/mL. Values above the measuring range are reported as  $> 8$  ng/mL (or up to 16.0 ng/mL for 2-fold diluted samples).

### Analytical Specificity

Hemoglobin up to 500 mg/dL, bilirubin up to 10 mg/dL, triglycerides up to 1800 mg/dL, total protein up to 10 g/dL, rheumatoid factor up to 200 IU/mL or antinuclear antibody will not interfere with the CL-series T3 assay. Criterion: recovery within  $\pm 10$  % of initial value.

In vitro tests were performed on 24 commonly used pharmaceuticals. No interference with the assay was found.

No interference was observed from the following substances at the levels indicated below. Criterion: recovery within  $\pm 10$  % of initial value.

Samples should not be taken from patients receiving high dose biotin supplements. If medically possible, these patients should discontinue biotin use before sample drawing to minimize the risk of interference.

Substance	Concentration
Phenylbutazone	1000 ng/mL
Methimazole	1.00 mg/dL
Biotin	20 ng/mL
Sodium salicylate	1000 ng/mL
6-N-propyl-2-thiouracil	1000 ng/mL

No obvious cross reactivity was observed when Mindray Total T3 Calibrator C0 was spiked with the T3 analogue at specific levels indicated in the table below. Criterion: Cross reactivity  $\leq 3$ % (Tetraiodothyroacetic acid: Cross reactivity  $\leq 5$ %). The results are summarized in the table below\*.

Substance	Cross-reactant Concentration	Cross Reactivity
3,5-L-diiodothyronine	100 ng/mL	1.59%
Reverse T3	500 ng/mL	0.01%
Tetraiodothyroacetic acid	25 ng/mL	0.33%
L-thyroxine	100 ng/mL	-0.01%
Monoiodotyrosine	1000 ng/mL	0.01%
Diiodo-L-Tyrosine	1000 ng/mL	0.01%

\*Representative data; results in individual laboratories may vary.

### Accuracy

Two trueness controls with traceable and defined values were used to verify the accuracy of this assay. The results showed that the relative deviation was within  $\pm 10\%$ . The results are listed in the following table\*.

Sample	Measured T3 Value (ng/mL)	Defined T3 Value (ng/mL)	Relative Deviation
Level 1	1.65	1.64	0.66%
Level 2	3.73	3.59	3.86%

\*Representative data; results in individual laboratories may vary.

### Precision

Precision was determined by following Clinical and Laboratory Standards Institute (CLSI) Protocol EP5-A2<sup>12</sup>. Two levels of quality controls were tested in duplicate in two separate runs per day, for a total of 20 days, using a single lot of reagents and a single calibration curve. The precision data are summarized in the table below\*.

Sample	Mean T3 (ng/mL)	Within-run CV	Between-run CV	Within-Device CV
1	2.05	1.26%	1.63%	2.16%
2	4.18	1.66%	1.20%	2.64%

\*Representative data; results in individual laboratories may vary.

### Linearity

A high concentration T3 sample (approximately 8.0 ng/mL) was mixed with a low concentration sample ( $< 0.2$  ng/mL) at different ratios, generating a series of dilutions. The total T3 of each dilution was determined using the Mindray CL-Series T3 Assay. Linearity was demonstrated in the range of 0.2 ng/mL to 8.0

ng/mL, the correlation coefficient  $r$  is  $\geq 0.990$ . The linearity data are summarized in the table below\*.

Concentration (ng/mL)	1	2	3	4	5	6
Expected T3	0.00	1.79	3.59	5.38	7.17	8.96
Measured T3	0.00	1.83	3.46	4.93	6.66	8.96

\*Representative data; results in individual laboratories may vary.

### Method Comparison

The Mindray CL-Series T3 Assay was compared to a commercially available diagnostic kit in a correlation study with about 304 specimens. The statistical data obtained by linear regression are shown in the table below.

Concentration Range(ng/mL)	Slope	Intercept	Correlation Coefficient
0.2-8	1.0066	-0.002	0.9892

### Warnings and Precautions

- For in vitro diagnostic use only. For laboratory professional use.
- Follow all the rules in handling laboratory reagents and take necessary safety precautions.
- The concentration of total T3 in a given specimen determined with different manufacturers can vary due to differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the T3 assay used. Values obtained with different assay methods cannot be used interchangeably.
- Do not use reagent kits beyond the expiration date.
- Do not use reagents mixed from different reagent lots.
- Always keep the reagent pack in the upright position to ensure no microparticle has been lost prior to use.
- Reagent pack opened for more than 56 days is not recommended for use.
- Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
- All the specimen and reaction wastes should be considered potentially biohazard. The handling of specimens and reaction wastes should be in accordance with the local regulations and guidelines.
- The Material Safety Data Sheet (MSDS) is available upon

request.

11. Please confirm the integrity of the package before use. Do not use the reagents with damaged packages.
12. If the reagents are opened unintentionally before use, they shall be used as soon as possible.
13. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the local competent authority.
14. Instability or deterioration should be suspected if there are visible signs of leakage, turbidity, precipitates or microbial growth.
15. Do not freeze. The results can't be assured when the reagents are stored at inappropriate condition.
16. This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

#### **Warning**

H317 May cause an allergic skin reaction.

H412 Harmful to aquatic life with long lasting effects.

#### **Prevention:**

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P272 Contaminated work clothing should not be allowed out of the workplace.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P273 Avoid release to the environment.

#### **Response:**

P302 + P352 IF ON SKIN: Wash with plenty of water.

P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.

P362 + P364 Take off contaminated clothing and wash it before reuse.

#### **Disposal:**

P501 Dispose of contents/ container in accordance with local regulation.

## Graphical Symbols

In vitro  
diagnostic  
medical  
device

Authorized  
representative  
in the European  
Community

Consult  
instructions  
for use

European  
Conformity

Catalogue  
number

Temperature  
limit

Manufacturer

Use-by date

Caution

Batch code

This way up

Unique device  
identifier

## References

1. Larsen PR. Triiodothyronine: Review of Recent Studies of Its Physiology and Pathophysiology in Man. *Metabolism* 1972; 21: 1073-1092.
2. Klee GG. Clinical Usage Recommendations and Analytic Performance Goals for Total and Free Triiodothyronine Measurements. *Clinical Chemistry* 1996; 42: 155-159.
3. Fisher DA. Physiological Variations in Thyroid Hormones; Physiological and Pathophysiological Considerations. *Clinical Chemistry* 1996; 42: 135-139.
4. Ekins, R. Measurement of Free Hormones in Blood. *Endocrinol. rev.*, 1986; 11: 5-46.
5. Gornall AG, Luxton AW, Bhavnani BR. Endocrine Disorders. In *Applied Biochemistry of Clinical Disorders*. 1986; 305-318. Edited by Gornall, AG Philadelphia, PA: J. B. Lippincott Co.
6. White GH. Recent Advances in Routine Thyroid Function Testing. *CRC - Critical Reviews in Clinical Laboratory Sciences*. 1986; 24: 315-362.
7. Kaplan MM, Larsen PR, Crantz FR, Deau VJ, Rossing TH. Prevalence of Abnormal Thyroid Function Test Results in Patients with Acute Medical Illnesses. *Am J Med* 1982; 72: 9-16.

8. Wahner HW, Gorman CA. Interpretation of Serum Tri-Iodothyronine Levels Measured by the Sterling Technic. N Engl J Med 1971; 2824: 225-230.
9. Boscatto LM, Stuart MC. Heterophilic Antibodies: a Problem for All Immunoassays. Clin Chem 1988; 34: 27-33.
10. Kricka L. Interferences in Immunoassays - Still a Threat. Clin Chem 2000; 46: 1037-1038.
11. Bjerner J, et al. Immunometric Assay Interference: Incidence and Prevention. Clin Chem 2002; 48: 613-621.
12. CLSI. EP5-A2: Vol. 24, No. 25, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline – Second Edition.

© 2013-2023 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

All rights Reserved

**Manufacturer:** Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

**Address:** Mindray Building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R.China.

**E-mail Address:** service@mindray.com

**Website:** www.mindray.com

**Tel:** +86-755-81888998

**Fax:** +86-755-26582680

**EC-Representative:** Shanghai International Holding Corp. GmbH (Europe)

**Address:** Eiffestraße 80, Hamburg 20537, Germany

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1. Larsen PR. Triiodothyronine: Review of Recent Studies of Its Physiology and Pathophysiology in Man. *Metabolism* 1972; 21: 1073-1092.
2. Klee GG. Clinical Usage Recommendations and Analytic Performance Goals for Total and Free Triiodothyronine Measurements. *Clinical Chemistry* 1996; 42: 155-159.
3. Fisher DA. Physiological Variations in Thyroid Hormones; Physiological and Pathophysiological Considerations. *Clinical Chemistry* 1996; 42: 135-139.
4. Ekins, R. Measurement of Free Hormones in Blood. *Endocrinol. rev.*, 1986; 11: 5-46.
5. Gornall AG, Luxton AW, Bhavnani BR. Endocrine Disorders. In *Applied Biochemistry of Clinical Disorders*. 1986; 305-318. Edited by Gornall, AG Philadelphia, PA: J. B. Lippincott Co.
6. White GH. Recent Advances in Routine Thyroid Function Testing. *CRC - Critical Reviews in Clinical Laboratory*

Sciences. 1986; 24: 315-362.

7. Kaplan MM, Larsen PR, Crantz FR, Deau VJ, Rossing TH. Prevalence of Abnormal Thyroid Function Test Results in Patients with Acute Medical Illnesses. Am J Med 1982; 72: 9-16.
8. Wahner HW, Gorman CA. Interpretation of Serum Tri-Iodothyronine Levels Measured by the Sterling Technic. N Engl J Med 1971; 2824: 225-230.
9. Boscato Lm, Stuart MC. Heterophilic Antibodies: a Problem for All Immunoassays. Clin Chem 1988; 34: 27-33.
10. Kricka L. Interferences in Immunoassays - Still a Threat. Clin Chem 2000; 46: 1037-1038.
11. Bjerner J, et al. Immunometric Assay Interference: Incidence and Prevention. Clin Chem 2002; 48: 613-621.
12. CLSI. EP5-A2: Vol. 24, No. 25, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline – Second Edition.

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### Referências

1. Larsen PR. Triiodothyronine: Review of Recent Studies of Its Physiology and Pathophysiology in Man. *Metabolism* 1972; 21: 1073-1092.
2. Klee GG. Clinical Usage Recommendations and Analytic Performance Goals for Total and Free Triiodothyronine Measurements. *Clinical Chemistry* 1996; 42: 155-159.
3. Fisher DA. Physiological Variations in Thyroid Hormones; Physiological and Pathophysiological Considerations. *Clinical Chemistry* 1996; 42: 135-139.
4. Ekins, R. Measurement of Free Hormones in Blood. *Endocrinol. rev.*, 1986; 11: 5-46.
5. Gornall AG, Luxton AW, Bhavnani BR. Endocrine Disorders. In *Applied Biochemistry of Clinical Disorders*. 1986; 305-318. Edited by Gornall, AG Philadelphia, PA: J. B. Lippincott Co.
6. White GH. Recent Advances in Routine Thyroid Function Testing. *CRC - Critical Reviews in Clinical Laboratory Sciences*. 1986; 24: 315-362.
7. Kaplan MM, Larsen PR, Crantz FR, Deau VJ, Rossing TH. Prevalence of Abnormal Thyroid Function Test Results in Patients with Acute Medical Illnesses. *Am J Med* 1982; 72: 9-16.
8. Wahner HW, Gorman CA. Interpretation of Serum Tri-Iodothyronine Levels Measured by the Sterling Technic. *N Engl J Med* 1971; 2824: 225-230.
9. Boscato Lm, Stuart MC. Heterophilic Antibodies: a Problem for All Immunoassays. *Clin Chem* 1988; 34: 27-33.
10. Kricka L. Interferences in Immunoassays - Still a Threat. *Clin Chem* 2000; 46: 1037-1038.
11. Bjerner J, et al. Immunometric Assay Interference: Incidence and Prevention. *Clin Chem* 2002; 48: 613-621.
12. CLSI. EP5-A2: Vol. 24, No. 25, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline - Second Edition.

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**Fabricante:** Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

**Endereço:** Mindray Building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R.China.

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### Referencias

1. Larsen PR. Triiodothyronine: Review of Recent Studies of Its Physiology and Pathophysiology in Man. *Metabolism* 1972; 21: 1073-1092.
2. Klee GG. Clinical Usage Recommendations and Analytic Performance Goals for Total and Free Triiodothyronine Measurements. *Clinical Chemistry* 1996; 42: 155-159.
3. Fisher DA. Physiological Variations in Thyroid Hormones; Physiological and Pathophysiological Considerations. *Clinical Chemistry* 1996; 42: 135-139.
4. Ekins, R. Measurement of Free Hormones in Blood. *Endocrinol. rev.*, 1986; 11: 5-46.
5. Gornall AG, Luxton AW, Bhavnani BR. Endocrine Disorders. In *Applied Biochemistry of Clinical Disorders*. 1986; 305-318. Edited by Gornall, AG Philadelphia, PA: J. B. Lippincott Co.
6. White GH. Recent Advances in Routine Thyroid Function Testing. *CRC - Critical Reviews in Clinical Laboratory Sciences*. 1986; 24: 315-362.
7. Kaplan MM, Larsen PR, Crantz FR, Deau VJ, Rossing TH. Prevalence of Abnormal Thyroid Function Test Results in Patients with Acute Medical Illnesses. *Am J Med* 1982; 72: 9-16.
8. Wahner HW, Gorman CA. Interpretation of Serum Tri-Iodothyronine Levels Measured by the Sterling Technic. *N Engl J Med* 1971; 2824: 225-230.
9. Boscatto Lm, Stuart MC. Heterophilic Antibodies: a Problem for All Immunoassays. *Clin Chem* 1988; 34: 27-33.
10. Kricka L. Interferences in Immunoassays - Still a Threat. *Clin Chem* 2000; 46: 1037-1038.

11. Bjerner J, et al. Immunometric Assay Interference: Incidence and Prevention. Clin Chem 2002; 48: 613-621.
12. CLSI. EP5-A2: Vol. 24, No. 25, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline – Second Edition.

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**Riferimenti**

1. Larsen PR. Triiodothyronine: Review of Recent Studies of Its Physiology and Pathophysiology in Man. *Metabolism* 1972; 21: 1073-1092.
2. Klee GG. Clinical Usage Recommendations and Analytic Performance Goals for Total and Free Triiodothyronine

- Measurements. Clinical Chemistry 1996; 42: 155-159.
3. Fisher DA. Physiological Variations in Thyroid Hormones; Physiological and Pathophysiological Considerations. Clinical Chemistry 1996; 42: 135-139.
  4. Ekins, R. Measurement of Free Hormones in Blood. Endocrinol. rev., 1986; 11: 5-46.
  5. Gornall AG, Luxton AW, Bhavnani BR. Endocrine Disorders. In Applied Biochemistry of Clinical Disorders. 1986; 305-318. Edited by Gornall, AG Philadelphia, PA: J. B. Lippincott Co.
  6. White GH. Recent Advances in Routine Thyroid Function Testing. CRC - Critical Reviews in Clinical Laboratory Sciences. 1986; 24: 315-362.
  7. Kaplan MM, Larsen PR, Crantz FR, Deau VJ, Rossing TH. Prevalence of Abnormal Thyroid Function Test Results in Patients with Acute Medical Illnesses. Am J Med 1982; 72: 9-16.
  8. Wahner HW, Gorman CA. Interpretation of Serum Tri-Iodothyronine Levels Measured by the Sterling Technic. N Engl J Med 1971; 2824: 225-230.
  9. Boscato Lm, Stuart MC. Heterophilic Antibodies: a Problem for All Immunoassays. Clin Chem 1988; 34: 27-33.
  10. Kricka L. Interferences in Immunoassays - Still a Threat. Clin Chem 2000; 46: 1037-1038.
  11. Bjerner J, et al. Immunometric Assay Interference: Incidence and Prevention. Clin Chem 2002; 48: 613-621.
  12. CLSI. EP5-A2: Vol. 24, No. 25, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline - Second Edition.

© 2013-2023 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.  
Tutti i diritti riservati

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**Indirizzo:** Mindray Building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R.China.

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Seri kodu

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#### Referanslar

1. Larsen PR. Triiodothyronine: Review of Recent Studies of Its Physiology and Pathophysiology in Man. *Metabolism* 1972; 21: 1073-1092.
2. Klee GG. Clinical Usage Recommendations and Analytic Performance Goals for Total and Free Triiodothyronine Measurements. *Clinical Chemistry* 1996; 42: 155-159.
3. Fisher DA. Physiological Variations in Thyroid Hormones; Physiological and Pathophysiological Considerations. *Clinical Chemistry* 1996; 42: 135-139.
4. Ekins, R. Measurement of Free Hormones in Blood.

- Endocrinol. rev., 1986; 11: 5-46.
5. Gornall AG, Luxton AW, Bhavnani BR. Endocrine Disorders. In Applied Biochemistry of Clinical Disorders. 1986; 305-318. Edited by Gornall, AG Philadelphia, PA: J. B. Lippincott Co.
  6. White GH. Recent Advances in Routine Thyroid Function Testing. CRC - Critical Reviews in Clinical Laboratory Sciences. 1986; 24: 315-362.
  7. Kaplan MM, Larsen PR, Crantz FR, Deau VJ, Rossing TH. Prevalence of Abnormal Thyroid Function Test Results in Patients with Acute Medical Illnesses. Am J Med 1982; 72: 9-16.
  8. Wahner HW, Gorman CA. Interpretation of Serum Tri-iodothyronine Levels Measured by the Sterling Technic. N Engl J Med 1971; 2824: 225-230.
  9. Boscato Lm, Stuart MC. Heterophilic Antibodies: a Problem for All Immunoassays. Clin Chem 1988; 34: 27-33.
  10. Kricka L. Interferences in Immunoassays - Still a Threat. Clin Chem 2000; 46: 1037-1038.
  11. Bjerner J, et al. Immunometric Assay Interference: Incidence and Prevention. Clin Chem 2002; 48: 613-621.
  12. CLSI. EP5-A2: Vol. 24, No. 25, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline - Second Edition.

© 2013-2023 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

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**Üretici:** Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

**Adres:** Mindray Building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R.China.

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Haut

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**Références**

1. Larsen PR. Triiodothyronine: Review of Recent Studies of Its Physiology and Pathophysiology in Man. *Metabolism* 1972; 21: 1073-1092.
2. Klee GG. Clinical Usage Recommendations and Analytic Performance Goals for Total and Free Triiodothyronine Measurements. *Clinical Chemistry* 1996; 42: 155-159.
3. Fisher DA. Physiological Variations in Thyroid Hormones; Physiological and Pathophysiological Considerations. *Clinical Chemistry* 1996; 42: 135-139.
4. Ekins, R. Measurement of Free Hormones in Blood. *Endocrinol. rev.*, 1986; 11: 5-46.
5. Gornall AG, Luxton AW, Bhavnani BR. *Endocrine Disorders*.

- In Applied Biochemistry of Clinical Disorders. 1986; 305-318. Edited by Gornall, AG Philadelphia, PA: J. B. Lippincott Co.
- White GH. Recent Advances in Routine Thyroid Function Testing. CRC - Critical Reviews in Clinical Laboratory Sciences. 1986; 24: 315-362.
  - Kaplan MM, Larsen PR, Crantz FR, Deau VJ, Rossing TH. Prevalence of Abnormal Thyroid Function Test Results in Patients with Acute Medical Illnesses. Am J Med 1982; 72: 9-16.
  - Wahner HW, Gorman CA. Interpretation of Serum Tri-Iodothyronine Levels Measured by the Sterling Technic. N Engl J Med 1971; 2824: 225-230.
  - Boscato Lm, Stuart MC. Heterophilic Antibodies: a Problem for All Immunoassays. Clin Chem 1988; 34: 27-33.
  - Kricka L. Interferences in Immunoassays - Still a Threat. Clin Chem 2000; 46: 1037-1038.
  - Bjerner J, et al. Immunometric Assay Interference: Incidence and Prevention. Clin Chem 2002; 48: 613-621.
  - CLSI. EP5-A2: Vol. 24, No. 25, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline - Second Edition.

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**Fabricant** : Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
**Représentant en Europe** : Shanghai International Holding Corp. GmbH (Europe)