Total T4 Calibrators

Order Information

Catalog No.	Package Size
105-004280-00	C0:1×2.0 mL/vial C1:1×2.0 mL/vial C2:1×2.0 mL/vial

Intended Use

Mindray Total T4 Calibrators (T4 CAL) are intended to calibrate the quantitative total thyroxine (T4) assay on Mindray CL-series Chemiluminescence Immunoassay Analyzer.

Summary

Mindray CL measurement system is composed of Mindray CL-series Chemiluminescence Immunoassay Analyzer, Mindray reagent kits, calibrators, and controls. Total T4 Calibrators contain thyroxine (T4) analyte. The mathematical relationship between the measured responses and the known analyte concentration establishes the calibration curve. It is used to convert Relative Light Units (RLUs) measured from samples to quantitative analyte concentration.

Components

C0, C1, C2	T4 (synthetic) at three concentration levels with 0.05% ProClin 300 and 0.09% sodium azide as preservatives in human serum
------------	--

Warnings and Precautions

- For in vitro diagnostic use only. For laboratory professional use.
- Laboratory professionals must be trained to take all necessary precautions while handling laboratory reagents. Do not use reagents beyond the expiration date.
- The calibrators' values are lot-specific with the matched set listed in the Calibration Card.
- Calibration is recommended after each reagent lot change, special maintenance, or troubleshooting procedure.
- Do not mix calibrators from different lots.
- Please perform maintenance and standard operation as scheduled, including calibration and analysis to assure the performance of the measurement system.
- Recalibrate the system when controls are out of specified range.
- The calibrators were tested with CE-approved kits for HIV-1/2 Antibody, HCV Antibody, Treponema pallidum Antibody and HBsAg. All results are negative. However, as no testing method can rule out the potential risk of infection with absolute certainty, this material should be handled as potential biohazards.¹
- If there is any evidence of microbial contamination or excessive turbidity in the product, discard the vial.

- Material safety data sheet is available upon request.
- Disposal of all waste material should be in accordance with local guidelines.
- Please confirm the integrity of the package before use. Do not use the reagents with damaged packages. The results can't be assured when the reagents are stored at inappropriate condition.
- If the reagents are opened unintentionally before use, they shall be tightly capped and stored at 2-8°C, and then used within in-use stability time.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the local competent authority.
- This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

Warning

H317 May cause an allergic skin reaction.
H412 Harmful to aquatic life with long lasting effects.
Prevention:
P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
P272 Contaminated work clothing should not be allowed out of the workplace.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P273 Avoid release to the environment.
Response:
P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.
P302 + P352 IF ON SKIN: Wash with plenty of water.
P362 + P364 Take off contaminated clothing and wash it before reuse.
Disposal:
P501 Dispose of contents/ container in accordance with local regulation.

Testing Method

- Refer to Mindray System Operations Manual and / or customer service to acquire information of calibration principle, use of calibrators, and validity of data of calibrators.
- Mix contents by gently inverting the vial no less than 10 times before use. Avoid bubble formation.
- Three levels of T4 calibrators are provided, which are 0, ~4, ~30 µg/dL. Test each calibrator twice for calibration.
- Refer to calibration card for accurate information about lot-specific T4 concentration.

Preparation

- Take the vial out of refrigerator and let it stand for no less than 30 minutes at room temperature.
- According to the calibration procedure described in the instrument manual, use the Calibrators' lot-specific values to set calibration parameters.

Storage and Stability

The calibrators are stable up to the expiration date as indicated on the label when stored in unopened vial at 2-8°C. Once opened, it is stable for 30 days at 2-8°C or 90 days at -20°C.

Materials Required but not Provided

- Mindray CL-series Chemiluminescence Immunoassay Analyzer, reagent kits and controls.
- General laboratory equipment.

Traceability of Calibrator

The traceability process is based on EN ISO 17511: 2021,² the measurand (analyte) in the calibrators is traceable to the IDMS (isotope dilution mass spectrometry) reference measurement procedure. The concentrations of calibrators are lot-specific on Mindray CL-series Chemiluminescence Immunoassay Analyzer, and cannot be used on the other system or method, as bias may exist among different systems and methods. Details about traceability are available upon request.

Quality Control

Use Mindray controls as the samples to monitor the performance of system. The control result should be within the defined range as shown in the control value sheet. If the controls are out of specified range, the measurement system should be checked. For instance, expiration date or storage condition of the calibrators, reagent and control; contamination of reagent; position of reagent or sample placed on the analyzer; the parameter setting in software; or performance of the analyzer. The assigned values of controls are fitted to Mindray CL-series Chemiluminescence Immunoassay Analyzer and are lot-specific.

Graphical Symbols

References

- HHS Publication, 5th ed., December 2009. Biosafety in Microbiological and Biomedical Laboratories.
- EN ISO 17511:2021 In vitro diagnostic medical devices – Requirements for establishing metrological traceability of values assigned to calibrators, trueness control materials and human samples.

© 2022-2023 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

All rights Reserved

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen, 518057 P.R.China

E-mail Address: service@mindray.com

Website: www.mindray.com

Tel: +86-755-81888998

Fax: +86-755-26582680

EC-Representative: Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

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Fabricant : Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Représentant en Europe : Shanghai International Holding Corp. GmbH (Europe)