LH

Luteinizing Hormone (CLIA)

Order information

Catalog No.	Package Size
LH111	2×50 tests
LH112	2×100 tests
LH113	2×30 tests

Intended Use

The CL-series LH assay is a Chemiluminescent Immunoassay (CLIA) for the quantitative determination of luteinizing hormone (LH) in human serum or plasma.

Summary

Human luteinizing hormone (hLH) is a glycoprotein consisting of two subunits (α - and β -chains), with 219 amino acids in total.^{1,2} Besides of pituitary, there are LH-like materials in rat hypothalamus, which exert a negative feedback on the release of LH. There are also materials similar with HCG-LH in 6~8d rabbit blastocyst and 18d cattle conceptus, which stimulate corpus luteum cells to synthesize progesterone.³ The α -subunit of LH is essentially identical to the α -subunits of follicle stimulating hormone (FSH, follitropin), thyroid stimulating hormone (TSH, thyrotropin), and human chorionic gonadotropin (HCG) in molecular structure. The β -subunit is considerably different from that of FSH and TSH. However, the β -subunits of LH and HCG are similar. The different β -subunits determine their immune and physiological specificities.⁴ LH, together with FSH, is secreted by the gonadotroph cells in the pituitary in response to the secretion of the gonadotropin releasing hormone (GnRH) from the medial basal hypothalamus.^{3,5}

LH and FSH regulate and stimulate the growth and function of the gonads (ovaries and testes) synergistically. LH regulate the menstrual cycle in females with FSH synergistically, causing the release of ovum and the formation of corpus luteum, which secretes progesterone and estradiol. The biosynthesis of LH was controlled by the hypothalamus gonadotropin releasing hormone (GnRH). Meanwhile, the ovarian steroids, exert an either a positive or negative feedback on the release of LH. The determination of LH in conjunction with FSH is utilized for clarifying the causes of amenorrhea (primary amenorrhea or secondary amenorrhea) in women and clarifying the causes of testicular failure in men. The surge of LH has a compact relation with the ovulation in menstrual cycle. The Ovulation is always occurred in 24~36 hours after the surge of LH. So the determination and monitoring of the concentration of serum LH in menstrual cycle can define the optimal period of fertility. Series detections of LH for 6-8 days in mid-cycle of women with regular menstrual cycles can make out the release surge. Determination of LH is useful in the diagnosis and treatment of

infertility and pituitary and gonadal disorders.⁶

Assay Principle

The CL-series LH assay is a two-site sandwich assay to determine the level of LH.

In the first step, sample and paramagnetic microparticles coated with mouse monoclonal anti-LH antibody are added into a reaction vessel. After incubation, LH present in the sample binds to the mouse monoclonal anti-LH antibody coated microparticles. Microparticles are magnetically captured while other unbound substances are removed by washing.


In the second step, mouse monoclonal anti-LH antibody alkaline phosphatase labeled conjugate was added to the reaction vessel. The mouse monoclonal anti-LH antibody alkaline phosphatase labeled conjugate binds to another reaction site of the LH that has bound to the paramagnetic microparticles coated with LH antibody to form a sandwich complex. Microparticles are magnetically captured while other unbound substances are removed by washing.

In the third step, the substrate solution is added to the reaction vessel. It is catalyzed by mouse monoclonal anti-LH antibody-alkaline phosphatase conjugate in the immunocomplex retained on the microparticle. The resulting chemiluminescent reaction is measured as relative light units (RLUs) by a photomultiplier built into the system. The amount of LH present in the sample is proportional to relative light units (RLUs) generated during the reaction. The LH concentration can be determined via a calibration curve.

Reagent Components

Ra	Paramagnetic microparticles coated with monoclonal anti-LH antibody (mouse) in TRIS buffer with preservative
Rb	Monoclonal anti-LH antibody (mouse)-alkaline phosphatase conjugate in MES buffer with preservative.
Rc	MES Buffer with preservative.

The position of each reagent component is shown in the figure below (front view on the left

and top view on the right):

Storage and Stability

The unopened LH (CLIA) reagent kit is stable up to the stated expiration date when stored at 2-8°C.

The LH (CLIA) reagent kit can be stored onboard and used for a maximum of 28 days after opening at 2-8°C.

Reagent Preparation

Ra: Ready to use

Rb: Ready to use

Rc: Ready to use

Materials Required But Not Provided

Mindray CL-series Chemiluminescence Immunoassay Analyzer

Cat.No.LH211: Mindray LH Calibrators, 1×2.0 mL for each of calibrator C0, C1 and C2.

Cat.No.LH212: Mindray LH Calibrators, 1×1.2 mL C0, 1×1.0 mL C1 and 1×1.0 mL C2.

Cat.No.RHL311: Mindray Reproductive Multi Control (L), 6×5.0 mL.

Cat.No.RHH311: Mindray Reproductive Multi Control (H), 6×5.0 mL.

Cat.No.RHL312: Mindray Reproductive Multi Control (L), 12×5.0 mL.

Cat.No.RHH312: Mindray Reproductive Multi Control (H), 12×5.0 mL

Cat.No.RHL313: Mindray Reproductive Multi Control (L), 1×5.0 mL.

Cat.No.RHH313: Mindray Reproductive Multi Control (H), 1×5.0 mL.

Cat.No.RHL314: Mindray Reproductive Multi Control (L), 3×5.0 mL.

Cat.No.RHH314: Mindray Reproductive Multi Control (H), 3×5.0 mL

Cat.No.WB411: Mindray Wash Buffer, 1×10 L.

Cat.No.CS511: Mindray Substrate Solution, 4×115 mL.

Cat.No.CS512: Mindray Substrate Solution, 4×75 mL.

Mindray Reaction Vessels.

Applicable Instrument

Mindray CL-series Chemiluminescence Immunoassay Analyzer

Specimen Collection and Preparation

Human serum or plasma collected in EDTA, sodium heparin and lithium heparin are recommended for this assay.

Centrifuge the specimens after clot formation is complete. Transfer the supernatants into tubes for storage or test within two hours after centrifugation.

Specimens should be tested as soon as possible after sample collection. If testing is not completed within 8 hours, specimens should be tightly capped and refrigerated at 2-8°C. If testing will be delayed for more than 72 hours, specimens should be frozen at -20°C or below.

Avoid repeated freeze and thaw cycles.

Assay Procedure

For optimal performance of this assay, operators should read the related system operation manual carefully, to get sufficient information such as operation instructions, sample preservation and management, safety precaution, and maintenance. Prepare all

required materials for the assay as well.

Before loading the LH (CLIA) reagent kit on the machine for the first time, unopened reagent bottle should be inverted gently for at least 30 times to resuspend the microparticles that have settled during shipment or storage. Visually inspect the bottle to ensure the microparticles have been resuspended. If the microparticles remain adhered to the bottle, continue inverting until the microparticles have been completely resuspended. If the microparticles cannot be resuspended, it is recommended not to use this bottle of reagent. Contact Mindray Customer Service for help. Do not invert opened reagent bottle.

This assay requires 15 μ L of sample for a single test. This volume does not include the dead volume of the sample container. Additional volume is required when performing additional tests from the same sample. Operators should refer to the system operation manual and specific requirement of the assay to determine the minimum sample volume.

Calibration

CL-series LH (CLIA) has been standardized against the WHO International Standard Luteinizing Hormone, Human, Pituitary (NIBSC code: 80/552).

The specific information of master calibration curve of LH (CLIA) reagent kit is stored in the two-dimensional barcode attached in the reagent pack. It's used together with calibrators for the calibration of the specific reagent lot. When performing the calibration, scan the information of master calibration curve from the barcode into the system first, and then use the calibrators at three levels. Valid calibration curve is required before any LH test. Recalibration is recommended every 4 weeks, or when a new reagent lot is used, or the quality controls are out of specified range. For detailed instruction of calibration, refer to the system operation manual.

Quality control

It is recommended that quality controls should be run once every 24 hours if the tests are in use, or after every calibration. The quality control frequency should be adapted to each laboratory's individual requirements. The recommended two levels of quality controls for this assay are Mindray Reproductive Multi Control (L) and Reproductive Multi Control (H).

Quality control results should be within the acceptable ranges. If a control is out of its specified range, the associated test results are invalid and the samples must be retested. Recalibration may be required. Examine the assay system referring to the system operation manual. If the quality control results are still out of the specified range, please contact Mindray Customer Service for help.

Calculation

The analyzer automatically calculates the analyte concentration of each sample on the

master calibration curve read from the barcode, and a 4-Parameter Logistic Curve Fitting (4PLC) with the relative light units (RLUs) generated from three level calibrators of defined concentration values. The results are shown in the unit of mIU/mL.

Conversion factors: mIU/mL x 1 = IU/L

Dilution

Samples with LH concentrations above the upper limit can be diluted with Mindray Sample Diluent. The recommended dilution is 1:4 (either automatically by the analyzer or manually). The concentration of the diluted sample must be > 2 mIU/mL. After manual dilution, multiply the result by the dilution factor. After automated dilution by the analyzers, the system automatically multiply the result by the dilution factor when calculating the sample concentration.

Expected values

An extended study on a cohort of 565 healthy individuals (135 males and 430 females) has determined the reference range of CL-series LH assay.

Category		N	5-95 th Percentile
Male		135	1.4-7.7 mIU/mL
Female	follicular phase	125	1.9-11.6 mIU/mL
	ovulation period	43	12.9-105.2 mIU/mL
	luteal phase	128	0.8-10.5 mIU/mL
	post-menopause	134	6.6-66.4 mIU/mL

Due to the variation in geography, race, sex, and age, it is highly recommended that each laboratory should establish its own reference range.

Limitation

The upper limit of this assay is 250 mIU/mL.

A specimen with a LH concentration lower than the upper limit can be quantitatively determined, while specimen with a concentration higher than the upper limit will be reported as >250 mIU/mL or diluting the samples with Mindray Sample Diluent.

The concentration of LH in a given specimen, determined with assays from different manufacturers, can vary due to differences in assay methods, calibration, and reagent specificity. The assay results should be used in conjunction with other data, such as symptoms, results of other tests, clinical history, etc.

Specimen from individuals who have been exposed to mouse monoclonal antibodies may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values with assay kits employing mouse monoclonal antibodies. However, no obvious interference of HAMA was observed in this assay.

Performance Characteristics

Analytical Sensitivity/Limit of Detection

The LH (CLIA) reagent kit has an analytical

sensitivity of ≤ 0.2 mIU/mL. Analytical sensitivity is defined as the lowest concentration of analyte that can be differentiated from a sample that contains no analyte. It is defined as the LH concentration at two standard deviations above the mean RLU from 20 measurements of an analyte-free sample.

Reportable Range

Reportable range is defined by the analytical sensitivity and the upper limit of the master calibration curve. The reportable range of LH (CLIA) reagent kit is 0.2-250 mIU/mL (or the upper limit is up to 1000 mIU/mL for 4-fold diluted samples).

Specificity

Hemoglobin up to 500 mg/dL, bilirubin up to 10 mg/dL, triglycerides up to 1800 mg/dL, and total protein up to 10.0 g/dL will not interfere with the CL-series LH assay. These substances show less than 10% interferences at indicated concentration.

No obvious interference was observed from rheumatoid factor up to 1500 IU/mL or antinuclear antibody up to 4000 U/L.

No significant cross-reactivity was observed when human chorionic gonadotropin (HCG), follicle stimulating hormone (FSH), and thyroid stimulating hormone (TSH) were added to the LH Calibrator C0, at specific levels indicated in the table below. The results showed no more than 2.0 mIU/mL in the CL-series LH assay at the levels indicated below.

Substance	Interference Substance Concentration
HCG	1000 mIU/mL
FSH	200 mIU/mL
TSH	100 μ IU/mL

Accuracy

A sample of WHO International Standard Luteinizing Hormone, Human, Pituitary (NIBSC code: 80/552) with traceable and defined value was used to verify the accuracy of this assay. The results showed that the relative deviation was less than $\pm 10\%$. The results are listed in the following table.

Sample	Measured LH Value (mIU/mL)	Defined LH Value (mIU/mL)	Relative Deviation
WHO LH	50.32	50.89	-1.13%

Precision

The CL-series LH assay is designed to have a precision of $\leq 10\%$ (within-device CV). Precision was determined by following National Committee for Clinical Laboratory Standards (NCCLS) Protocol EP5-A2. 7 Two levels of quality controls were tested in duplicate in two separate runs per day, for a total of 20 days, using a single lot of reagents and a single calibration curve. The precision data are summarized in the table below.

Sample	Mean LH (mIU/mL)	Within-run CV	Between-run CV	Within-Device CV
--------	------------------	---------------	----------------	------------------

1	10,21	3,29%	7,51%	8,20%
2	50,37	3,56%	5,20%	7,00%

Linearity

A high concentration LH sample (approximately 100 mIU/mL) was mixed with a low concentration sample (< 0.4 mIU/mL) at different ratios, generating a series of dilutions. The LH of each dilution was determined using the Mindray CL-Series LH Assay. Linearity was demonstrated in the range of 0.4 mIU/mL to 100 mIU/mL. The correlation coefficient r is ≥ 0.990 . The linearity data are summarized in the table below.

Concentration (mIU/mL)	1	2	3	4	5	6
Expected LH	0.31	27.38	54.45	81.52	108.59	135.66
Measured LH	0.31	27.64	54.91	82.82	107.68	135.66

Method Comparison


The Mindray CL-Series LH Assay was compared to a commercially available diagnostic kit in a correlation study with about 422 specimens. The statistical data obtained by Deming computing mode are shown in the table below.

Concentration Range (mIU/mL)	Slope	Intercept	Correlation Coefficient
0.15-210.85	1.1719	-0.9300	0.997

Warnings and Precautions

- For in vitro diagnostic use only.
- Follow all the rules in handling laboratory reagents and take necessary safety precautions.
- Due to the differences in methodology and antibody specificity, test results of the same sample may be different when using reagent kits from different manufacturers on Mindray system, or using Mindray reagent kits on other systems.
- Do not use reagent kits beyond the expiration date.
- Do not use reagents mixed from different reagent lots.
- Always keep the reagent pack in the upright position to ensure no microparticle has been lost prior to use.
- Reagent pack opened for more than 28 days is not recommended for use.
- Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
- All the specimen and reaction wastes should be considered potentially biohazard. The handling of specimens and reaction wastes should be in accordance with the local regulations and guidelines.
- The Material Safety Data Sheet (MSDS) is available upon request.

Graphical symbols

References

- Bishop WH, Nureddin A, Ryan RJ. Pituitary Luteinizing and Follicle Stimulating Hormones. In: Parsons JA, editor. Peptide Hormones. Baltimore: University Park Press, 1976; 273-98.
- Keutmann HT, Williams RM, Ryan RJ. Structure of Human Luteinizing Hormone Beta Subunit: Evidence for a Related Carboxyl-Terminal Sequence among Certain Peptide Hormones. Biochem Biophys Res Commun 1979; 90:842-8.
- Schroff RW, Foon KA, Beatty SM, Oldham RK, Morgan AC Jr. Human Anti-Murine Immunoglobulin Responses in Patients Receiving Monoclonal Antibody Therapy. Cancer Res 1985; 45:879-85.
- Vaitukaitis JL, Ross GT. Antigenic similarities among the human glycoprotein hormones and their subunits. In Gonadotropins. Edited by Saxena BB, Gandy HM, Beling CG. 1972; New York: John Wiley and Sons, 435-443.
- Beastall GH, Ferguson KM, O'Reilly DSJ, Seth J, Sheridan B. Assays for follicle stimulating hormone and luteinizing hormone: Guidelines for the provision of a clinical biochemistry service. Ann Clin Biochem. 1987; 24:246-262.
- Ross GT. Disorders of the Ovary and Female Reproductive Tract. In: Wilson JD and Foster DW, editors. Williams Textbook of Endocrinology. Philadelphia: Saunders, 1985; 206-58.
- NCCLS. EP5-A2 Evaluation of precision performance of quantitative measurement method; approved guideline-second edition. 2004.

© 2013 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

All rights Reserved

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, Shenzhen, 518057 P.R.China

E-mail Address: service@mindray.com

Website: www.mindray.com

Tel: +86-755-26582888

Fax: +86-755-26582680

EC-Representative: Shanghai International Holding Corp. GmbH(Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

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- Bishop WH, Nureddin A, Ryan RJ. Pituitary Luteinizing and Follicle Stimulating Hormones. In: Parons JA, editor. Peptide Hormones. Baltimore: University Park Press, 1976; 273-98.
- Keutmann HT, Williams RM, Ryan RJ. Structure of Human Luteinizing Hormone Beta Subunit: Evidence for a Related Carboxyl-Terminal Sequence among Certain Peptide Hormones. Biochem Biophys Res Commun 1979; 90:842-8.
- Schroff RW, Foon KA, Beatty SM, Oldham RK, Morgan AC Jr. Human Anti-Murine Immunoglobulin Responses in Patients Receiving Monoclonal Antibody Therapy. Cancer Res 1985; 45:879-85.
- Vaitukaitis JL, Ross GT. Antigenic similarities among the human glycoprotein hormones and their subunits. In: Gonadotropins. Edited by Saxena BB, Gandy HM, Beling CG. 1972; New York: John Wiley and Sons, 435-443.
- Beastall GH, Ferguson KM, O'Reilly DSJ, Seth J, Sheridan B. Assays for follicle stimulating hormone and luteinizing hormone: Guidelines for the provision of a clinical biochemistry service. Ann Clin Biochem. 1987; 24:246–262.
- Ross GT. Disorders of the Ovary and Female Reproductive Tract. In: Wilson JD and Foster DW, editors. Williams Textbook of Endocrinology. Philadelphia: Saunders, 1985; 206-58.
- NCCLS. EP5-A2 Evaluation of precision performance of quantitative measurement method; approved guideline—second edition. 2004.

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
Símbolos gráficos

Referências

1. Bishop WH, Nureddin A, Ryan RJ. Pituitary Luteinizing and Follicle Stimulating Hormones. Em: Parons JA, editor. Peptide Hormones. Baltimore: University Park Press, 1976; 273-98.
2. Keutmann HT, Williams RM, Ryan RJ. Structure of Human Luteinizing Hormone Beta Subunit: Evidence for a Related Carboxyl-Terminal Sequence among Certain Peptide Hormones. Biochem Biophys Res Commun 1979; 90:842-8.
3. Schroff RW, Foon KA, Beatty SM, Oldham RK, Morgan AC Jr. Human Anti-Murine Immunoglobulin Responses in Patients Receiving Monoclonal Antibody Therapy. Cancer Res 1985; 45:879-85.
4. Vaitukaitis JL, Ross GT. Antigenic similarities among the human glycoprotein hormones and their subunits. Em Gonadotropins. Editado por Saxena BB, Gandy HM, Belling CG. 1972; Nova York: John Wiley and Sons, 435-443.
5. Beastall GH, Ferguson KM, O'Reilly DSJ, Seth J, Sheridan B. Assays for follicle stimulating hormone and luteinizing hormone: Guidelines for the provision of a clinical biochemistry service. Ann Clin Biochem. 1987; 24:246-262.
6. Ross GT. Disorders of the Ovary and Female Reproductive Tract. Em: Wilson JD e Foster DW, editores. Williams Textbook of Endocrinology. Philadelphia: Saunders, 1985; 206-58.
7. NCCLS. EP5-A2 Evaluation of precision performance of quantitative measurement method; approved guideline-second edition. 2004.

© 2013 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

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Fabricante: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Endereço: Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, ShenZhen 518057, P.R. China.

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Símbolos gráficos

Referencias

- Bishop WH, Nureddin A, Ryan RJ. Pituitary Luteinizing and Follicle Stimulating Hormones. In: Parsons JA, editor. Peptide Hormones. Baltimore: University Park Press, 1976; 273-98.
- Keutmann HT, Williams RM, Ryan RJ. Structure of Human Luteinizing Hormone Beta Subunit: Evidence for a Related Carboxyl-Terminal Sequence among Certain Peptide Hormones. Biochem Biophys Res Commun 1979; 90:842-8.
- Schroff RW, Foon KA, Beatty SM, Oldham RK, Morgan AC Jr. Human Anti-Murine Immunoglobulin Responses in Patients Receiving Monoclonal Antibody Therapy. Cancer Res 1985; 45:879-85.
- Vaitukaitis JL, Ross GT. Antigenic similarities among the human glycoprotein hormones and their subunits. In Gonadotropins. Edited by Saxena BB, Gandhi HM, Belling CG. 1972; New York: John Wiley and Sons, 435-443.
- Beastall GH, Ferguson KM, O'Reilly DSJ, Seth J, Sheridan B. Assays for follicle stimulating hormone and luteinizing hormone: Guidelines for the provision of a clinical biochemistry service. Ann Clin Biochem. 1987; 24:246-262.
- Ross GT. Disorders of the Ovary and Female Reproductive Tract. En: Wilson JD and Foster DW, editors. Williams Textbook of Endocrinology. Philadelphia: Saunders, 1985; 206-58.
- NCCLS. EP5-A2 Evaluation of precision performance of quantitative measurement method; approved guideline-second edition. 2004.

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Manufacturer - Fabbricante

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- Bishop WH, Nureddin A, Ryan RJ. Pituitary Luteinizing and Follicle Stimulating Hormones. In: Parons JA, editor. Peptide Hormones. Baltimore: University Park Press, 1976; 273-98.
- Keutmann HT, Williams RM, Ryan RJ. Structure of Human Luteinizing Hormone Beta Subunit: Evidence for a Related Carboxyl-Terminal Sequence among Certain Peptide Hormones. Biochem Biophys Res Commun 1979; 90:842-8.
- Schroff RW, Foon KA, Beatty SM, Oldham RK, Morgan AC Jr. Human Anti-Murine Immunoglobulin Responses in Patients Receiving Monoclonal Antibody Therapy. Cancer Res 1985; 45:879-85.
- Vaitukaitis JL, Ross GT. Antigenic similarities among the human glycoprotein hormones and their subunits. In: Gonadotropins. Edited by Saxena BB, Gandy HM, Belling CG. 1972; New York: John Wiley and Sons, 435-443.
- Beastall GH, Ferguson KM, O'Reilly DSJ, Seth J, Sheridan B. Assays for follicle stimulating hormone and luteinizing hormone: Guidelines for the provision of a clinical biochemistry service. Ann Clin Biochem. 1987; 24:246-262.
- Ross GT. Disorders of the Ovary and Female Reproductive Tract. In: Wilson JD and Foster DW, editors. Williams Textbook of Endocrinology. Philadelphia: Saunders, 1985; 206-58.
- NCCLS. EP5-A2 Evaluation of precision performance of quantitative measurement method; approved guideline-second edition. 2004

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Indirizzo: Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, ShenZhen 518057, Cina.

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
Grafiksel Semboller

Referanslar

- Bishop WH, Nureddin A, Ryan RJ. Pituitary Luteinizing and Follicle Stimulating Hormones. In: Parsons JA, editor. Peptide Hormones. Baltimore: University Park Press, 1976; 273-98.
- Keutmann HT, Williams RM, Ryan RJ. Structure of Human Luteinizing Hormone Beta Subunit: Evidence for a Related Carboxyl-Terminal Sequence among Certain Peptide Hormones. Biochem Biophys Res Commun 1979; 90:842-8.
- Schroff RW, Foon KA, Beatty SM, Oldham RK, Morgan AC Jr. Human Anti-Murine Immunoglobulin Responses in Patients Receiving Monoclonal Antibody Therapy. Cancer Res 1985; 45:879-85.
- Vaitukaitis JL, Ross GT. Antigenic similarities among the human glycoprotein hormones and their subunits. In: Gonadotropins. Edited by Saxena BB, Gandy HM, Belling CG. 1972; New York: John Wiley and Sons, 435-443.
- Beastall GH, Ferguson KM, O'Reilly DSJ, Seth J, Sheridan B. Assays for follicle stimulating hormone and luteinizing hormone: Guidelines for the provision of a clinical biochemistry service. Ann Clin Biochem. 1987; 24:246-262.
- Ross GT. Disorders of the Ovary and Female Reproductive Tract. In: Wilson JD and Foster DW, editors. Williams Textbook of Endocrinology. Philadelphia: Saunders, 1985; 206-58.
- NCCLS. EP5-A2 Evaluation of precision performance of quantitative measurement method; approved guideline-second edition. 2004.

© 2013 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

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Üretici: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Adres: Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, Shenzhen, 518057 Çin

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
Symboles graphiques

Références

- Bishop WH, Nureddin A, Ryan RJ. Pituitary Luteinizing and Follicle Stimulating Hormones. In: Parons JA, editor. Peptide Hormones. Baltimore: University Park Press. 1976; 273-98.
- Keutmann HT, Williams RM, Ryan RJ. Structure of Human Luteinizing Hormone Beta Subunit: Evidence for a Related Carboxyl-Terminal Sequence among Certain Peptide Hormones. Biochem Biophys Res Commun 1979; 90:842-8.
- Schroff RW, Foon KA, Beatty SM, Oldham RK, Morgan AC Jr. Human Anti-Murine Immunoglobulin Responses in Patients Receiving Monoclonal Antibody Therapy. Cancer Res 1985; 45:879-85.
- Vaitukaitis JL, Ross GT. Antigenic similarities among the human glycoprotein hormones and their subunits. In: Gonadotropins. Edited by Saxena BB, Gandy HM, Beling CG. 1972; New York: John Wiley and Sons. 435-443.
- Beastall GH, Ferguson KM, O'Reilly DSJ, Seth J, Sheridan B. Assays for follicle stimulating hormone and luteinizing hormone: Guidelines for the provision of a clinical biochemistry service. Ann Clin Biochem. 1987; 24:246-262.
- Ross GT. Disorders of the Ovary and Female Reproductive Tract. In: Wilson JD and Foster DW, editors. Williams Textbook of Endocrinology. Philadelphia: Saunders. 1985; 206-58.
- NCCLS. EP5-A2 Evaluation of precision performance of quantitative measurement method; approved guideline-second edition. 2004.

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Fabricant : Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Représentant EC : Shanghai International Holding Corp. GmbH (Europe)