Anti-Tg

Antibody to thyroglobulin (CLIA)

Order Information

Catalog No.	Package Size
Anti-Tg111	2×50 tests
Anti-Tg112	2×100 tests
Anti-Tg113	2×30 tests

Intended Use

The CL-series Anti-Tg assay is a Chemiluminescent Immunoassay (CLIA) for the quantitative determination of Antibody to thyroglobulin (Anti-Tg) in human serum or plasma.

Summary

The level of thyroglobulin autoantibodies are often rising in patients with autoimmune thyroid diseases¹. High concentrations of anti-Tg together with anti-TPO are indicative of chronic lymphocytic-infiltrative thyroiditis (Hashimoto's disease). The positive rate of anti-Tg in subjects with autoimmune-thyroiditis is approx.70–80%, and approx.30% in individuals with Graves' disease. The anti-Tg assay is important for use in the course of Hashimoto's disease and for the differential diagnosis of suspected autoimmune-thyroiditis of unknown origin with negative anti-TPO results, Graves' disease without lymphocytic infiltration². The anti-Tg assay is also valuable to rule out interference by thyroglobulin autoantibodies in the Tg test.

Assay Principle

The CL-series Anti-Tg assay is a two-site sandwich assay to determine the level of Antibody to thyroglobulin.

In the first step, sample, sample treatment reagent, paramagnetic microparticles coated with monoclonal anti-biotin antibody (sheep) and biotinylated Tg are added into a reaction cuvette. After incubation, anti-Tg antibody present in the sample binds to biotinylated Tg and biotinylated Tg binds to monoclonal anti-biotin antibody coated microparticles. Microparticles are magnetically captured while other unbound substances are removed by washing.


In the second step, the Tg-alkaline phosphatase conjugate is added to the reaction cuvette. The Tg-alkaline phosphatase conjugate binds to the anti-Tg antibody complexed with Tg on microparticles. Microparticles are magnetically captured while other unbound substances are removed by washing.

Then the substrate solution is added to the reaction cuvette. It is catalyzed by Tg-alkaline phosphatase conjugate in the immunocomplex retained on the microparticles. The resulting chemiluminescent reaction is measured as relative light units (RLUs) by a photomultiplier built inside the system. The amount of anti-Tg antibody present in the sample is proportional to the relative light units (RLUs) generated during the reaction. The anti-Tg antibody concentration can be determined via a calibration curve.

Reagent Components

Ra	Paramagnetic microparticles coated with monoclonal anti-biotin antibody (sheep) in MES buffer with preservative.
Rb	Tg-alkaline phosphatase conjugate in PB buffer with preservative.
Rc	biotinylated Tg in PB buffer with preservative.
Rd	Sample treatment reagent with preservative.

The position of each reagent component is shown in the figure below (front view on the left and top view on the right):

Storage and Stability

The unopened Anti-Tg (CLIA) reagent kit is stable up to the stated expiration date when stored at 2-8°C.

The Anti-Tg (CLIA) reagent kit can be stored onboard and used for a maximum of 56 days after opening at 2-8°C.

Reagent Preparation

Ra: Ready to use

Rb: Ready to use

Rc: Ready to use

Rd: Ready to use

Materials Required but not Provided

Mindray CL-series Chemiluminescence Immunoassay Analyzer.

Cat.No.Anti-Tg 211: Anti-Tg Calibrators, 1×2.0 mL for each of calibrator C0, C1 and C2.

Cat.No.Anti-Tg 212: Anti-Tg Calibrators, 1×1.2 mL C0, 1×1.0 mL C1 and 1×1.0 mL C2.

Cat.No. ATABL311: Anti-thyroid Antibodies Control (L), 3×2.0 mL.

Cat.No. ATABH311: Anti-thyroid Antibodies Control (H), 3×2.0 mL.

Cat.No. ATABL311: Anti-thyroid Antibodies Control (L), 6×2.0 mL.

Cat.No. ATABH311: Anti-thyroid Antibodies Control (H), 6×2.0 mL.

Cat.No.WB411: Wash Buffer, 1×10 L.

Cat.No.CS511: Substrate Solution, 4×115 mL.

Reaction Cuvette

Applicable Instrument

Mindray CL-series Chemiluminescence Immunoassay Analyzer

Specimen Collection and Preparation

Human serum or plasma collected in EDTA, sodium

heparin and lithium heparin are recommended for this assay.

Collect all blood samples observing routine precautions for venipuncture. Follow blood collection tube manufacturer's recommendations for centrifugation. Centrifuge the specimens after clot formation is complete. Some specimens, especially those from patients receiving anticoagulant therapy may exhibit increased clotting time. Please ensure that residual fibrin and cellular matter have been removed prior to analysis.

For optimal results, inspect all samples for bubbles. Remove bubbles with a pipette tip prior to analysis. Specimens must be mixed thoroughly after thawing. Thawed samples should be centrifuged prior to use. If the sample was covered with lipid layer after centrifugation, the sample should be transferred to a clean tube before testing. Do not transfer the lipid layer. Handle carefully to prevent cross contamination. Do not use grossly hemolyzed specimens. Do not use heat-inactivated specimens.

Specimens should be tested as soon as possible after sample collection. If testing is not completed within 8 hours, specimens should be stored at 2-8°C or colder. Specimens would be stable for 48 hours at 2-8°C, 3 months at -20°C. Avoid more than five freeze cycles.

Assay Procedure

For optimal performance of this assay, operators should read the related system operation manual carefully, to get sufficient information such as operation instructions, sample preservation and management, safety precaution, and maintenance. Prepare all required materials for the assay as well.

Before loading the Anti-Tg (CLIA) reagent kit on the machine for the first time, unopened reagent bottle should be inverted gently for at least 30 times to resuspend the microparticles that have settled during shipment or storage. Visually inspect the bottle to ensure the microparticles have been resuspended. If the microparticles remain adhered to the bottle, continue inverting until the microparticles have been completely resuspended. If the microparticles cannot be resuspended, it is recommended not to use this bottle of reagent. Contact Mindray Customer Service for help. Do not invert opened reagent bottle.

This assay requires 10 µL of sample for a single test. This volume does not include the dead volume of the sample container. Additional volume is required when performing additional tests from the same sample. Operators should refer to the system operation manual and specific requirement of the assay to determine the minimum sample volume.

Calibration

CL-series Anti-Tg (CLIA) has been standardized against the WHO International Standard Antibody to thyroglobulin 65/093.

The specific information of master calibration curve of Anti-Tg (CLIA) reagent kit is stored in the two-dimensional barcode attached in the reagent pack. It's used in combination with Anti-Tg calibrators for the calibration of the specific reagent lot. When

performing the calibration, scan the information of master calibration curve from the barcode into the system first, and then use the Anti-Tg calibrators at three levels. Valid calibration curve is required before any Anti-Tg test. Recalibration is recommended every 4 weeks, or when a new reagent lot is used, or the quality controls are out of specified range. For detailed instruction of calibration, refer to the system operation manual.

Quality Control

It is recommended that quality controls should be run once every 24 hours if the tests are in use, or after every calibration. The quality control frequency should be adapted to each laboratory's individual requirements. The recommended two levels of quality controls for this assay are Anti-thyroid Antibodies Control (L) and Anti-thyroid Antibodies Control (H).

Quality control results should be within the acceptable ranges. If a control is out of its specified range, the associated test results are invalid and the samples must be retested. Recalibration may be required. Examine the assay system referring to the system operation manual. If the quality control results are still out of the specified range, please contact Mindray Customer Service for help.

Calculation

The analyzer automatically calculates the analyte concentration of each sample on the master calibration curve read from the barcode, and a 4-Parameter Logistic Curve Fitting (4PLC) with the relative light units (RLUs) generated from Anti-Tg calibrators of defined concentration values. The results are shown in the unit of IU/mL.

Expected values

An study on a cohort of 480 apparently healthy individuals(The screening criteria included TSH levels between 0.5-2.0 IU/mL,no thyroid disease and absence of non-thyroid autoimmune disease) confirmed the reference value of 4 IU/ml; this value corresponds to the 95th percentile.

Due to the variation in geography, race, sex, and age, it is highly recommended that each laboratory should establish its own reference range.

Limitation

The measuring range of this assay is 0.9 IU/mL ~ 2500 IU/mL(defined by the lower detection limit and the upper limit).A specimen with an Anti-Tg concentration lower than the upper limit can be quantitatively determined, while the specimen with a concentration higher than the upper limit will be reported as >2500 IU/mL.

The concentration of Anti-Tg in a given specimen, determined with assays from different manufacturers, can vary due to differences in assay methods, calibration, and reagent specificity. The assay results should be used in conjunction with other data, such as symptoms, results of other tests, clinical history, etc.

Specimen from individuals who have been exposed to mouse monoclonal antibodies may contain human anti-mouse antibodies (HAMA)³. Such specimens may show either falsely elevated or falsely depressed

values with assay kits employing mouse monoclonal antibodies^{4,5}. However, no obvious interference of HAMA has been observed in this assay.

Performance Characteristics

Analytical Sensitivity/Limit of Detection

The Anti-Tg (CLIA) reagent kit has an analytical sensitivity of ≤ 0.9 IU/mL. Analytical sensitivity is defined as the lowest concentration of analyte that can be differentiated from a sample that contains no analyte. It is defined as the Anti-Tg concentration at two standard deviations above the mean RLU from 20 measurements of an analyte-free sample.

Reportable Range

Reportable range is defined by the analytical sensitivity and the upper limit of the master calibration curve. The reportable range of Anti-Tg (CLIA) reagent kit is 0.9-2500 IU/mL.

Specificity

Hemoglobin up to 500 mg/dL, bilirubin up to 40 mg/dL, triglycerides up to 3000 mg/dL, and total protein up to 10 g/dL will not interfere with the CL-series Anti-Tg assay. These substances show less than 10% interferences at indicated concentration.

No obvious interference was observed from rheumatoid factor (up to 1500 IU/mL) or antinuclear antibody.

Anti-Tg Calibrator C0 was spiked with anti-TPO. No obvious cross reactivity was observed. The results are summarized in the table below.

Substance	Cross-reactant Concentration	Reported Anti-Tg (IU/mL)	Acceptance Criteria
Anti-TPO	1400 IU/mL	1.63	Reported Anti-Tg \leq 5.0 IU/mL

High Dose Hook

For the CL-series Anti-Tg assay, no high dose hook effect was observed with samples containing up to approximately 350000 IU/mL of Anti-Tg were assayed.

Accuracy

Two samples traced to the European Community Bureau of Reference (BCR) CRM 457 thyroglobulin standard were used to verify the accuracy of this assay. The results showed that the relative deviation was less than $\pm 10\%$.

Sample	Measured Anti-Tg Value (IU/mL)	Defined Anti-Tg Value (IU/mL)	Relative Deviation
L	11.27	11.05	1.99%
H	248.27	250.36	-0.83%

Precision

The CL-series Anti-Tg assay is designed to have a

duplicate in two separate runs per day, for a total of 20 days, using a single lot of reagents and a single calibration curve. The precision data are summarized in the table below.

Sample	Mean Anti-Tg (IU/mL)	Within-run CV	Between-run CV	Within-Device CV
1	8,29	4,65%	6,40%	8,99%
2	261,32	4,70%	5,23%	8,94%

Method Comparison

The Mindray CL-Series Anti-Tg Assay was compared to a commercially available diagnostic kit in a correlation study with about 336 specimens. The statistical data obtained by Deming computing mode are summarized in the table below.

Concentration Range (IU/mL)	Slope	Intercept	Correlation Coefficient
0.01~2500	1.0327	-1.0641	0.9950

Warnings and Precautions

- For in vitro diagnostic use only.
- Follow all the rules in handling laboratory reagents and take necessary safety precautions.
- Due to the differences in methodology and antibody specificity, test results of the same sample may be different when using reagent kits from different manufacturers on Mindray system, or using Mindray reagent kits on other systems.
- Do not use reagent kits beyond the expiration date.
- Do not use reagents mixed from different reagent lots.
- Always keep the reagent pack in the upright position to ensure no microparticle has been lost prior to use.
- Reagent pack opened for more than 56 days is not recommended for use.
- Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
- All the specimen and reaction wastes should be considered potentially biohazard. The handling of specimens and reaction wastes should be in accordance with the local regulations and guidelines.
- The Material Safety Data Sheet (MSDS) is available upon request.

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practice guidelines – Laboratory support for the diagnosis and monitoring of thyroid disease. The National Academy of Clinical Biochemistry. Доступно на сайте <http://www.cdc.gov/niosh>.

- Boscatto LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27–33.
- Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
- Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.

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1. F
e
l
d
t

-Rasmussen, U. Analytical and clinical performance goals for testing autoantibodies to thyroperoxidase, thyroglobulin and thyrotropin receptor. Clinical Chemistry: 42:1 pp. 160-3. (1996).

2. Demers LM, Spencer CA. Laboratory medicine

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
Símbolos gráficos

Referências

- Feldt-Rasmussen, U. Analytical and clinical performance goals for testing autoantibodies to thyroperoxidase, thyroglobulin and thyrotropin receptor. Clinical Chemistry: 42:1 pp. 160-3. (1996).
- Demers LM, Spencer CA. Laboratory medicine practice guidelines – Laboratory support for the diagnosis and monitoring of thyroid disease. The National Academy of Clinical Biochemistry. Available <http://www.cdc.gov/niosh>.
- Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
- Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
- Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.

© 2015 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Todos os direitos reservados
Fabricante: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Endereço: Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, Shenzhen, 518057 P.R.China

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
- Demers LM, Spencer CA. Laboratory medicine practice guidelines – Laboratory support for the diagnosis and monitoring of thyroid disease. The National Academy of Clinical Biochemistry. Available <http://www.cdc.gov/niosh>.
- Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
- Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
- Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.

© 2015 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Todos los derechos reservados

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1. Feldt-Rasmussen, U. Analytical and clinical performance goals for testing autoantibodies to

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Tel: 0049-40-2513175

Faks: 0049-40-255726

- nd clinical performance goals for testing autoantibodies to thyroperoxidase, thyroglobulin and thyrotropin receptor. Clinical Chemistry: 42:1 pp. 160-3. (1996).
2. Demers LM, Spencer CA. Laboratory medicine practice guidelines – Laboratory support for the diagnosis and monitoring of thyroid disease. The National Academy of Clinical Biochemistry. [Http://www.cdc.gov/niosh](http://www.cdc.gov/niosh) adresinde mevcuttur.
 3. Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
 4. Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
 5. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
- © 2015 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Tüm hakları saklıdır
Üretici: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Adres: Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, Shenzhen, 518057 Çin
E-posta adresi: service@mindray.com
Web sitesi: www.mindray.com
Tel: +86-755-26582888
Faks: +86-755-26582680
EC-Temsilcisi: Shanghai International Holding Corp.

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rs LM, Spencer CA. Laboratory medicine practice guidelines – Laboratory support for the diagnosis and monitoring of thyroid disease. The National Academy of Clinical Biochemistry. Available <http://www.cdc.gov/niosh>.

3. Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
4. Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037–1038.
5. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.

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© 2015 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Tous droits réservés
Fabricant : Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Représentant EC : Shanghai International Holding Corp. GmbH (Europe)