HSV-1+2 IgM

Herpes Simplex Virus 1+2 IgM(CLIA)

INSTRUCTIONS FOR USE

mindray

HSV-1+2 IgM

Herpes Simplex Virus 1+2 IgM(CLIA)

Order Information

Catalog No.	Package Size
105-012488-00	2×50 tests
105-012489-00	2×100 tests
105-012491-00	2×50 tests(calibrators included)
105-012492-00	2×100 tests(calibrators included)

Intended Use

The CL-series Herpes Simplex Virus 1+2 IgM assay is a Chemiluminescent Immunoassay (CLIA) for the qualitative determination of Herpes Simplex Virus type 1 and/or 2 specific IgM in human serum and plasma. The results of this assay can be used as an aid in the diagnosis of Herpes Simplex Virus infection and in the assessment of immune status of individuals including pregnant women.

Summary

Herpes simplex virus (HSV), a member of the Herpesviridae family, containing a double-stranded DNA molecule, naturally only infects human beings ^[1]. HSV is mainly transmitted through direct contact with mucosal membrane or abraded skin. Its infection is quite common worldwide, and there is no significant seasonal difference. Most of the primary infections have no obvious symptoms, and a few manifests as localized herpes of mucosal membrane or skin. After the symptoms of HSV primary infection are alleviated, the virus can reach the neurons through nerve fibers and remain in the neurons for a long time. Once immunocompromised, the virus can proliferate again, and virus excretion can be detected at the lesion site, along with or without near point infection symptoms. The incubation and reactivation of HSV is an important reason for its high prevalence ^[2-3].

According to the serotype, herpes simplex virus can be divided into herpes simplex virus type 1 (HSV-1) and herpes simplex virus type 2 (HSV-2). HSV-1 is normally associated with orofacial herpes and keratitis in immunocompetent individuals and viral encephalitis in infants and immunodeficiency adults ^[4-5]. In less developed countries, about one-third of children have seroconversion of HSV-1 at the age of 5 years, and about 70-80% of children have seroconversion of HSV-1 by adolescence. Genital herpes is the most common sexually transmitted disease, and HSV-2 is the most common cause of genital herpes. It is estimated that the global prevalence of HSV-2 is about 11.3%, and about 20 million primary HSV-2 infection cases are increased every year. HSV-2 has

synergistic effect with HIV infection and transmission. Patients infected with HSV-2 are about three times more likely to be infected with HIV than patients not infected with HSV-2, and patients infected with both HSV-2 and HIV are more likely to transmit HIV to others. It is reported that HSV-1 induced genital herpes has been increasingly found in recent years, especially in high-income countries. HSV-1 induced genital herpes has less symptoms and lower recurrence rate compared with HSV-2 induced genital herpes in general [6-9].

Neonates may be exposed to HSV in the birth canal during delivery, resulting in neonatal herpes, which is manifested as herpes in the skin, mouth or eyes. In severe cases, even disseminated diseases or central nervous system lesion may occur, with a mortality rate of up to 30%. In addition, about 20% of the surviving neonatal herpes patients may have long-term complications of the nervous system [10].

The detection of HSV specific antibody is very important for the evaluation of patients' immune status. IgM antibody generally starts to rise after one week of HSV infection, and its positive results suggest that patients may have acute infection. IgG antibody generally rises after 2 weeks of infection, which is a sign of previous infection. The first exposure of susceptible individuals to HSV-1 or HSV-2 will lead to "primary infection". "Non-primary Initial infection" occurs when an individual who has antibodies to either HSV-1 or HSV-2 is infected with the other virus type for the first time. Because of the cross reactivity between the IgG antibodies of HSV-1 and HSV-2, "non-primary initial infection causes milder symptoms and fewer complications than "primary infection" [11-12].

Assay Principles

The CL-series HSV-1+2 IgM assay is a μ -capture method to qualitatively detect IgM antibodies to Herpes simplex virus.

In the first step, sample, sample treatment solution and paramagnetic microparticles coated with anti-human IgM monoclonal antibody are added into a reaction cuvette. After incubation, HSV IgM antibodies in the sample will be bound to anti-human IgM monoclonal antibody coated microparticles. Afterwards, microparticles are magnetically captured while other unbound substances are removed by washing.

In the second step, test diluent solution and Herpes simplex virus antigens -alkaline phosphatase conjugates are added to the reaction cuvette for incubation. After incubation, Herpes simplex virus antigens -alkaline phosphatase conjugates will form sandwich structure with microparticles captured HSV IgM antibodies. Afterwards, microparticles are magnetically captured while other unbound substances are removed by washing.

Then, the substrate solution is added to the reaction cuvette. It is catalyzed by Herpes simplex virus antigens-alkaline phosphatase conjugates in the immune-complex retained on the microparticles. The resulting chemiluminescent reaction is measured as relative light units (RLUs) by a photomultiplier in the instrument.

A direct relationship exists between the amount of HSV IgM antibodies in the sample and the RLUs generated during the reaction. The presence or absence of analytes in the sample is determined by comparing the chemiluminescent signal in the reaction to the cutoff signal determined from calibration. A cutoff index (COI) is calculated by Sample RLU/Cutoff RLU.

Reagent Components

The kit is composed of Ra, Rb, Rc, and Rd, as well as calibrators C0 and C1. Calibrators C0 and C1 are optional (refer to the "Order Information" section). Components in different lots of reagent kit cannot be interchanged and the detailed information of each component is listed below:

Ra	Paramagnetic microparticles coated with anti-human IgM monoclonal antibody. Minimum concentration: 0.3 mg/mL solids. Tris buffer: 50 mmol/L. Preservative: 0.048% Proclin 300. 3.8 mL/50 tests; 6.6 mL/100 tests.
Rb	Herpes simplex virus antigens-alkaline phosphatase conjugates in MES buffer. Minimum concentration: 20 ng/mL. MES buffer: 50 mmol/L. Preservative: 0.048% Proclin 300. 3.5 mL/50 tests; 6.3 mL/100 tests.
Rc	Sample treatment solution. 50 mmol/L Tris buffer. Preservative: 0.048% Proclin 300. 3.5 mL/50 tests; 6.3 mL/100 tests.
Rd	Test diluent solution. 50 mmol/L Tris buffer. Preservative: 0.048% Proclin 300. 7.8 mL/50 tests; 14.8 mL/100 tests.
Calibrators C0, C1 (optional)	C0: IgM antibodies to Herpes simplex virus negative in buffer with 0.05% Proclin 300; C1: IgM antibodies to Herpes simplex virus positive in buffer with 0.05% Proclin 300; C0: 1×0.6 mL/vial; C1: 1×0.5 mL/vial.

Storage and Stability

The unopened reagents (Ra, Rb, Rc and Rd) and calibrators (C0 and C1) are stable up to the expiration date indicated on the label, when stored at 2-8°C.

Once opened in use, the reagents (Ra, Rb, Rc and Rd) are stable for 28 days at 2-8°C.

Once opened in use, the calibrators are stable for 28 days when refrigerated at 2-8°C and stable for 90 days when frozen at -20°C (Frozen only once).

Contamination must be avoided.

Materials Required but not Provided

Mindray CL-series Chemiluminescence Immunoassay Analyzer

Mindray ToRCH IgM Positive Control:

1×3.0 mL/vial (REF: 105-024840-00),

3×3.0 mL/vial (REF: 105-024841-00),

3×1.0 mL/vial (REF: 105-024842-00),

6×1.0 mL/vial (REF: 105-024843-00).

Mindray ToRCH IgG/IgM Negative Control:

1×3.0 mL/vial (REF: 105-024824-00),

3×3.0 mL/vial (REF: 105-024825-00),

3×1.0 mL/vial (REF: 105-024826-00),

6×1.0 mL/vial (REF: 105-024827-00).

Mindray Substrate Solution: 4x115 mL (REF: 105-004274-00)

Mindray Substrate Solution: 4x75 mL (REF: 105-009044-00)

Mindray Wash Buffer, 1×10 L (REF: 105-004552-00)

Mindray Reaction cuvette

Instrument System

Mindray CL-series Chemiluminescence Immunoassay Analyzers

Specimen Collection and Preparation

Specimen Types

- Human serum and plasma (K₂-EDTA, K₃-EDTA, Na₂-EDTA, sodium heparin, lithium heparin and sodium citrate) are recommended for this assay.
- Blood collection tubes from various manufacturers may contain different raw materials and additives which could affect the test results in some cases. Not all available tubes of all manufacturers were tested by Mindray. Each laboratory needs to judge the applicability of the tubes and plasma separation products. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Specimen Conditions

- Do not use:
 - heat-inactivated specimens
 - grossly hemolyzed specimens
 - specimens with obvious microbial contamination
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot

formation.

Preparation for Analysis

- Follow blood collection tube manufacturer's recommendations for centrifugation. Centrifuge the specimens after clot formation is complete. Ensure residual fibrin and cellular matter has been removed prior to analysis.
- For optimal results, inspect all samples for bubbles. Remove bubbles with a pipette tip prior to analysis. Specimens must be mixed thoroughly after thawing. Thawed samples should be centrifuged prior to use.
- If the sample is covered with lipid layer after centrifugation, the sample should be transferred to a clean tube and centrifuged before testing. Do not transfer the lipid layer. Handle carefully to prevent cross contamination.

Specimen Storage

- Specimens should be tested as soon as possible after sample collection. If testing is not completed within 48 hours, specimens should be tightly capped and refrigerated at 2-8°C. The sample is stable for 7 days at 2-8°C and 180 days at -20°C or below.
- The samples may be frozen no more than 5 times.

Frozen or refrigerated samples, controls and calibrators must be restored to room temperature before testing.

Assay Procedure

For optimal performance of this assay, operators should read the related Operator's Manual carefully to get enough information such as operation instructions, sample preservation and management, safety precautions, and maintenance. Prepare all required materials for the assay as well.

Before loading the HSV-1+2 IgM reagent kit on the instrument for the first time, invert unopened reagent bottle gently for at least 30 times to resuspend the microparticles, which have settled during shipment or storage. Visually inspect the bottle to ensure the microparticles have been well mixed. If the microparticles remain adhered to the bottle, continue inverting until suspended completely. Avoid foam formation. If the microparticles cannot be homogenized, it is recommended not to use this bottle of reagent. Contact Mindray Customer Service for help. Avoid inverting opened reagent bottle.

This assay requires 8 µL of sample volume for a single test not including dead volume of the sample container. Additional volume is required when performing additional tests from the same sample. Operators should refer to the Operator's Manual and follow the specific requirements of the assay to determine the minimum sample volume.

Calibration

The calibrators are traceable to Mindray internal reference.

The calibration information is stored in the barcodes attached inside the reagent and calibrator pack. When performing the calibration, please scan the barcodes into the system.

Take the vial out of refrigerator and let it stand for no less than 30 minutes at room temperature. Mix contents by gently inverting the vial no less than 10 times before use. Avoid bubble formation. Test C0 three times and C1 twice for calibration. The original calibrator vials may be used if the entire provided volume is used. Because of possible evaporation effects, the calibrator on the analyzer should be measured as soon as possible and not more than 3 calibration procedures per original vial should be performed. After use, close the bottles immediately and store upright at 2-8°C. The original calibrators can also be transferred into empty labeled bottles in aliquots. Use the aliquots immediately or store them immediately at -20°C. One aliquot can be used only once.

A valid calibration is required before any HSV-1+2 IgM testing.

Renewed calibration is required as follows:

1. The previous calibration was performed 28 days before.
2. Quality control results fall outside of the expected ranges.
3. A reagent kit with a new lot number is used.

For detailed instruction of calibration, refer to the Operator's Manual.

Quality Control

Users can prepare quality controls with clinical samples. Reference ranges can be established according to protocol approved by individual laboratories.

It is recommended that quality controls should be run once every 24 hours if the tests are in use, and following each calibration. The quality control frequency should be adapted to each laboratory's requirements.

Quality control results should be within the acceptable ranges. If a control is out of its specified range, the associated test results are invalid and the samples must be retested. Recalibration may be required. Refer to the Operator's Manual to check up the system. If quality control findings are still outside the identified ranges, please contact Mindray Customer Service for help.

Calculation

The analyzer automatically calculates the HSV-1+2 IgM result of each sample, on the basis of calibration master curve read from the barcode, and the relative light units (RLUs) generated from 2 levels of product calibrators of defined values. The results are shown in the unit of COI.

Cutoff RLU = [(Mean RLU of C1 – Mean RLU of C0) x Calibration Coefficient] + Mean RLU of C0

Calibration Coefficient is specific for each reagent and calibrator lot.

COI = Sample RLU/Cutoff RLU

Specimen Dilution

Specimens cannot be diluted for Mindray HSV-1+2 IgM assay.

Interpretation of Results

Specimens with a COI<0.90 are considered to be non-reactive for IgM antibodies to HSV.

Specimens with a COI≥0.90 to <1.10 are considered to be indeterminate for IgM antibodies to HSV. A second sample should be collected for testing within an appropriate time (e.g. 2-4 weeks) or other test methods should be used for confirmation.

Specimens with a COI≥1.10 are considered to be reactive for IgM antibodies to HSV.

The results of different detection systems may be inconsistent because of the detection methods and traceability.

Limitation of Measurement

The assay results should not be used solely for confirmation or exclusion of HSV infection. Clinical decisions should be made in conjunction with other evidences, such as symptoms, clinical history, results of nucleic acid test, etc. It is recommended to test the HSV IgM and HSV IgG antibodies simultaneously to evaluate and monitor the patient's immune status.

A negative HSV-1+2 IgM test result cannot completely exclude the possibility of acute HSV infection. The body may have not produced enough antibody in the early stage of acute infection. If acute infection is suspected, nucleic acid testin should be performed, and IgG antibody monitoring should be performed after a period of time (for example, 2-4 weeks). If the negative conversion of herpes simplex virus IgG in the second sample or positive PCR supports the diagnosis of herpes simplex virus acute infection. It is recommended that first sample should be stored at -20°C or below, so that two samples can be tested at the same time and same conditions^[12].

For pregnant women with negative IgM and IgG antibodies, if their sexual partners have or have been infected with herpes simplex virus or are likely to be exposed to other source of HSV infection, they should be monitored for antibodies. If seroconversion of IgG antibodies occurs, acute infection is indicated.

Because there is a high risk of false positive in the serological test of herpes simplex virus (HSV) IgM antibodies of pregnant women, prenatal screening is not recommended for asymptomatic

pregnant women, the laboratory test of pregnant women cannot reliably identify the risk of fetus, and the test results of this reagent should not be used as the basis for terminating pregnancy alone.

It is recommended to perform nucleic acid test in one week before delivery to determine whether there is a risk of neonatal infection through the birth canal, because positive and negative IgM do not directly indicate the active status of the virion.

For patients with immunodeficiency disease or receiving immunosuppressive therapy, the antibody test results should be interpreted carefully.

Exogenous interference should be considered in case a patient with HSV IgM positive result has recently received blood transfusion.

HSV IgM is not only increased in the primary infection, but also in the recurrent infection and reinfection, elevated IgM may persist even for years after primary infection in some individuals. Therefore, IgM level alone cannot be used to determine the time of initial infection.

Due to the existence of heterophilic antibodies, patients who often contact animals or use animal derived immunoglobulin for treatment may have abnormal detection results^[13].

The samples of patients who have been treated with mouse monoclonal antibody preparation may contain human anti-mouse antibody (HAMA), which may interfere with the test results of kits containing mouse monoclonal antibodies^[14-15]. Mindray chemiluminescence kit contains anti-interference components, which can effectively reduce the impact of HAMA in samples, but the interference of HAMA may still exist in a few samples. Only combined with the clinical information or other laboratory test results can the patient's condition be correctly determined.

Performance Characteristics

Precision

Precision was determined by following Clinical & Laboratory Standards Institute (CLSI) protocol EP05-A3^[16]. 3 serum samples and 2 controls were tested 2 runs each day in duplicate for totally 20 days (n=80). Results are shown in the table below (criterion: 4.5%CV or 0.045SD for repeatability and 8.0%CV or 0.08SD for intermediate precision)*.

Sample	Mean (COI)	Repeatability		Intermediate precision	
		SD(COI)	CV(%)	SD(COI)	CV(%)
HS1	0.44	0.01	1.55%	0.02	3.46%
HS2	0.99	0.02	1.62%	0.05	5.46%
HS3	2.35	0.04	1.85%	0.17	7.07%
NC	0.08	0.00	1.58%	0.00	2.90%
PC	2.50	0.03	1.34%	0.17	6.63%

*Representative data; results in individual laboratories may vary.

HS=Human serum. NC=negative control. PC=positive control.

Hook effect

When a sample containing a high titer of determinant is tested, the hook effect may lead to false negative results. This assay adopts two-step mode to keep the analysis signal at a high value, which could effectively avoid false negative results caused by hook effect.

2 samples with high titer of HSV-1+2 IgM were detected, all results showed the expected high concentration and there was no misjudgment of results.

Interference

The test results of Mindray HSV-1+2 IgM are not observed to be interfered with the potential endogenous substances up to the following concentrations (criterion: recovery within $\pm 15\%$ of initial value or absolute deviation within ± 0.15 COI)*.

Substance	Concentration
Hemoglobin	1000 mg/dL
Bilirubin	40 mg/dL
Triglyceride	37mmol/L
Total protein	12 g/dL
Albumin	6 g/dL
γ -Globulin	6 g/dL
Biotin	3500 ng/mL

**Representative data; results in individual laboratories may vary.*

The test results of Mindray HSV-1+2 IgM are not observed to be interfered with the potential exogenous substances up to the following concentrations (criterion: recovery within $\pm 15\%$ of initial value or absolute deviation within ± 0.15 COI)*.

Substance	Concentration
Acetaminophen	200 mg/L
Amoxicillin	700 mg/L
Cefalor	500 mg/L
Vitamin C	250 mg/L
Aspirin	1000 mg/L
Ibuprofen	500 mg/L
Acyclovir	6 mg/L
Famciclovir	12 mg/L
Valaciclovir	15 mg/L

**Representative data; results in individual laboratories may vary.*

Cross-reactivity

Studies were performed to evaluate the effect of potential cross reacting diseases on HSV-1+2 IgM assay. A total of 197 specimens were tested by Mindray HSV-1+2 IgM and a comparator HSV-1+2 IgM kit.

Substance	n	Substance	n
Treponema pallidum	10	parainfluenza	5
Influenza A	1	Influenza B	5
Chlamydia trachomatis	2	Mycoplasma pneumoniae	4
Candida albicans	5	Human herpes Virus 6	6
Gardella vaginalis	5	Human herpes Virus 8	6
Neisseria	5	Human papillomavirus	4
Escherichia coli	3	Hepatitis C virus	8
Rubella virus	6	HIV	5
Cytomegalo virus	9	Influenza vaccine	6
Hepers simplex virus 1+2 IgG	3	Human T cell leukemia virus	8
Toxoplasma gondii	7	SLE	7
Epstein-Barr virus	10	AMA	8
Parvo virus B19	9	RF	10
Hepatitis B virus	5	ANA	10
Hepatitis A virus	9	HAMA	8
Varicella-zoster virus	8	/	/

The results showed that the overall agreement of Mindray HSV-1+2 IgM and commercial HSV-1+2 IgM kits is 96.45%.

Method comparison

781 clinical specimens from pregnant women, women of childbearing age, neonates and patients with corresponding clinical symptoms were tested to evaluate the coincidence rate of Mindray HSV-1+2 IgM and a mainstream commercial chemiluminescence HSV-1+2 IgM kit. Specimens with inconsistent test results were confirmed by another commercial HSV-1+2 IgM kit. The positive coincidence rate was 96.80% (121/125, 95% CI 92.06%-98.75%); the negative coincidence rate was 99.37% (627/631, 95% CI 98.38%-99.75%); the total coincidence rate was 98.94% (748/767, 95% CI 97.93%-99.46%).

Mindray HSV-1+2 IgM	Comparator HSV-1+2 IgM			Total
	Reactive	Non-reactive	Indeterminate	
Reactive	121	4*	1	126
Non-reactive	4**	627	19	650
Indeterminate	1***	3***	1	5
Total	126	634	21	781

4*: 4 specimens with reactive on Mindray HSV-1+2 IgM and non-reactive on comparator HSV-1+2 IgM were found non-reactive on another commercial HSV-1+2 IgM kit.

4**: 4 specimens with non-reactive on Mindray HSV-1+2 IgM and reactive on comparator HSV-1+2 IgM were found reactive on another commercial HSV-1+2 IgM kit.

*1*** and 3***: 4 specimens with indeterminate on Mindray HSV-1+2 IgM and non-indeterminate on comparator HSV-1+2 IgM were found indeterminate on another commercial HSV-1+2 IgM kit.*

Warnings and Precautions

1. For *in vitro* diagnostic use only. For laboratory professional use.
2. Follow all the rules in handling laboratory reagents and take necessary safety precautions.
3. Due to the differences in methodology and antibody specificity, test results of the same sample may be different when using reagent kits from different manufacturers on Mindray system, or using Mindray reagent kits on other systems.
4. Do not use reagent kits beyond the expiration date.
5. Do not use reagents mixed from different reagent lots.
6. Please confirm the integrity of the package before use. Do not use the reagents with damaged packages. The results can't be assured when the reagents are stored at inappropriate condition.
7. If the reagents are opened unintentionally before use, they shall be used as soon as possible.
8. Wear clean gloves to avoid contaminating the reagents.
9. Always keep the reagent pack in the upright position to ensure no microparticle has been lost prior to use.
10. The calibrator was tested with CE-marked methods and shown to be negative for antibodies to HIV and HCV and nonreactive for HBsAg. However, as no testing method can rule out the potential risk of infection with absolute certainty, this material should be handled as a patient specimen.
11. Reagent pack opened for more than 14 days is not recommended for use.
12. Instability or deterioration should be suspected if there are visible signs of leakage, turbidity, precipitates or microbial growth.
13. Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
14. All the specimen and reaction wastes should be considered potentially biohazard. Specimens and reaction wastes should be handled in accordance with the local regulations and guidelines.
15. The Material Safety Data Sheet (MSDS) is available upon request.
16. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the local competent authority.
17. This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008.

**Warning**

H317 May cause an allergic skin reaction.

H412 Harmful to aquatic life with long lasting effects.

Prevention:

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P272 Contaminated work clothing should not be allowed out of the workplace.

P273 Avoid release to the environment.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

Response:

P302 + P352 IF ON SKIN: Wash with plenty of water.

P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.

P362 + P364 Take off contaminated clothing and wash it before reuse.

Disposal:

P501 Dispose of contents/ container in accordance with local regulation.

Graphical Symbols

In vitro
diagnostic
medical
device

Authorized
representative
in the European
Community

Consult
instructions
for use

European
Conformity

Catalogue
number

Temperature limit

Manufacturer

Use-by date

Biological
risks

Batch code

This way
up

Unique
device
identifier

Caution

Bibliography

1. Whitley R J, Roizman B. Herpes simplex virus infections[J]. 2001, 357(9267):1513-1518.
2. Centers for Disease Control and Prevention. Sexually Transmitted Infections Treatment Guidelines, 2021. MMWR Recomm Rep 2021;70(4).
3. World Health Organization. WHO guidelines for the treatment of genital herpes simplex virus. Geneva Switzerland WHO, 2016.
4. Sili U, Kaya A, Mert A. Herpes simplex virus encephalitis: Clinical manifestations, diagnosis and outcome in 106 adult patients[J]. Journal of Clinical Virology, 2014, 60(2):112-118.
5. Kaye S, Choudhary A. Herpes simplex keratitis[J]. Progress in Retinal and Eye Research, 2006, 25(4):355-380.
6. Fujie, Xu, Maya, et al. Trends in Herpes Simplex Virus Type 1 and Type 2 Seroprevalence in the United States[J]. Jama, 2006.
7. McQuillan G, Kruszon-Moran D, Flagg EW, Paulose-Ram R. Prevalence of herpes simplex virus type 1 and type 2 in persons aged 14-49: United States, 2015-2016. NCHS Data Brief 2018; 304:1-8.
8. Fleming D T, Mcquillan G M, Johnson R E, et al. Herpes simplex virus type 2 in the United States, 1976 to 1994.[J]. N Engl J Med, 1997, 337(16):1105-1111.
9. World Health Organization. Herpes simplex virus type 2: programmatic and research priorities in developing countries. Report of a WHO / UNAIDS / LSHTM workshop (London 14-16 February 2001). Geneva Switzerland WHO, 2001.
10. Obstetricians R. Management of genital herpes in pregnancy[J]. RCOG Green-top Guideline, 2007.
11. Bibbins-Domingo K, Grossman D C, Curry S J, et al. Serologic Screening for Genital Herpes Infection: US Preventive Services Task Force Recommendation Statement[J]. Jama, 2016, 316(23):2525.
12. Mendelson E, Aboudy Y, Smetana Z, et al. Laboratory assessment and diagnosis of congenital viral infections: Rubella, cytomegalovirus (CMV), varicella-zoster virus (VZV), herpes simplex virus (HSV), parvovirus B19 and

- human immunodeficiency virus (HIV).[J]. Reproductive Toxicology, 2006, 21(4):350-382.
13. Boscato Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
 14. Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
 15. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
 16. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline-Third Edition.

© 2022-2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

All rights Reserved.

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

E-mail Address: service@mindray.com

Website: www.mindray.com

Tel: +86-755-81888998

Fax: +86-755-26582680

EC-Representative: Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

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1. Whitley R J, Roizman B. Herpes simplex virus infections[J] (Уитли Р. Дж., Ройзман Б. Вирус простого герпеса [J]). 2001, 357(9267):1513-1518.
2. Centers for Disease Control and Prevention. Sexually Transmitted Infections Treatment Guidelines, 2021. MMWR Recomm Rep 2021;70(4).
3. World Health Organization. WHO guidelines for the treatment of genital herpes simplex virus. Geneva Switzerland WHO, 2016.
4. Sili U, Kaya A, Mert A. Herpes simplex virus encephalitis: Clinical manifestations, diagnosis and outcome in 106 adult patients[J]. Journal of Clinical Virology, 2014, 60(2):112-118.
5. Kaye S, Choudhary A. Herpes simplex keratitis[J]. Progress in Retinal and Eye Research, 2006, 25(4):355-380.
6. Fujie, Xu, Maya, et al. Trends in Herpes Simplex Virus Type 1 and Type 2 Seroprevalence in the United States[J]. Jama, 2006.
7. McQuillan G, Kruszon-Moran D, Flagg EW, Paulose-Ram R. Prevalence of herpes simplex virus type 1 and type 2 in persons aged 14-49: United States, 2015-2016. NCHS Data Brief 2018; 304:1-8.
8. Fleming D T, Mcquillan G M, Johnson R E, et al. Herpes simplex virus type 2 in the United States, 1976 to 1994.[J]. N Engl J Med, 1997, 337(16):1105-1111.
9. World Health Organization. Herpes simplex virus type 2: programmatic and research priorities in developing countries. Report of a WHO / UNAIDS / LSHTM workshop (London 14-16 February 2001). Geneva Switzerland WHO, 2001.
10. Obstetricians R. Management of genital herpes in pregnancy[J]. RCOG Green-top Guideline, 2007.
11. Bibbins-Domingo K, Grossman D C, Curry S J, et al. Serologic Screening for Genital Herpes Infection: US Preventive Services Task Force Recommendation Statement[J]. Jama, 2016, 316(23):2525.
12. Mendelson E, Aboudy Y, Smetana Z, et al. Laboratory assessment and diagnosis of congenital viral infections: Rubella, cytomegalovirus (CMV), varicella-zoster virus (VZV), herpes simplex virus (HSV), parvovirus B19 and human immunodeficiency virus (HIV).[J]. Reproductive Toxicology, 2006, 21(4):350-382.
13. Boscatto Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
14. Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.

15. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
16. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline-Third Edition.

© 2022-2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Все права защищены.

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Адрес: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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Atenção

Bibliografia

1. Whitley R J, Roizman B. Herpes simplex virus infections[J]. 2001, 357(9267):1513-1518.
2. Centers for Disease Control and Prevention. Sexually Transmitted Infections Treatment Guidelines, 2021. MMWR Recomm Rep 2021;70(4).
3. World Health Organization. WHO guidelines for the treatment of genital herpes simplex virus. Geneva Switzerland WHO, 2016.
4. Sili U, Kaya A, Mert A. Herpes simplex virus encephalitis: Clinical manifestations, diagnosis and outcome in 106 adult patients[J]. Journal of Clinical Virology, 2014, 60(2):112-118.
5. Kaye S, Choudhary A. Herpes simplex keratitis[J]. Progress in Retinal and Eye Research, 2006, 25(4):355-380.
6. Fujie, Xu, Maya, et al. Trends in Herpes Simplex Virus Type 1 and Type 2 Seroprevalence in the United States[J]. Jama, 2006.
7. McQuillan G, Kruszon-Moran D, Flagg EW, Paulose-Ram R. Prevalence of herpes simplex virus type 1 and type 2 in persons aged 14-49: United States, 2015-2016. NCHS Data Brief 2018; 304:1-8.
8. Fleming D T, Mcquillan G M, Johnson R E, et al. Herpes simplex virus type 2 in the United States, 1976 to 1994.[J]. N Engl J Med, 1997, 337(16):1105-1111.
9. World Health Organization. Herpes simplex virus type 2: programmatic and research priorities in developing countries. Report of a WHO / UNAIDS / LSHTM workshop (London 14-16 February 2001). Geneva Switzerland WHO, 2001.
10. Obstetricians R. Management of genital herpes in pregnancy[J]. RCOG Green-top Guideline, 2007.
11. Bibbins-Domingo K, Grossman D C, Curry S J, et al. Serologic Screening for Genital Herpes Infection: US

- Preventive Services Task Force Recommendation Statement[J]. Jama, 2016, 316(23):2525.
12. Mendelson E, Aboudy Y, Smetana Z, et al. Laboratory assessment and diagnosis of congenital viral infections: Rubella, cytomegalovirus (CMV), varicella-zoster virus (VZV), herpes simplex virus (HSV), parvovirus B19 and human immunodeficiency virus (HIV).[J]. Reproductive Toxicology, 2006, 21(4):350-382.
 13. Boscatto Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
 14. Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
 15. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
 16. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline-Third Edition.

© 2022-2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

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Fabricante: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Endereço: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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Precaución

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1. Whitley R J, Roizman B. Herpes simplex virus infections[J]. 2001, 357(9267):1513-1518.
2. Centers for Disease Control and Prevention. Sexually Transmitted Infections Treatment Guidelines, 2021. MMWR Recomm Rep 2021;70(4).
3. World Health Organization. WHO guidelines for the treatment of genital herpes simplex virus. Geneva Switzerland WHO, 2016.
4. Sili U, Kaya A, Mert A. Herpes simplex virus encephalitis: Clinical manifestations, diagnosis and outcome in 106 adult patients[J]. Journal of Clinical Virology, 2014, 60(2):112-118.
5. Kaye S, Choudhary A. Herpes simplex keratitis[J]. Progress in Retinal and Eye Research, 2006, 25(4):355-380.
6. Fujie, Xu, Maya, et al. Trends in Herpes Simplex Virus Type 1 and Type 2 Seroprevalence in the United States[J]. Jama, 2006.
7. McQuillan G, Kruszon-Moran D, Flagg EW, Paulose-Ram R. Prevalence of herpes simplex virus type 1 and type 2 in

- persons aged 14-49: United States, 2015-2016. NCHS Data Brief 2018; 304:1-8.
8. Fleming D T, Mcquillan G M, Johnson R E, et al. Herpes simplex virus type 2 in the United States, 1976 to 1994.[J]. N Engl J Med, 1997, 337(16):1105-1111.
 9. World Health Organization. Herpes simplex virus type 2: programmatic and research priorities in developing countries. Report of a WHO / UNAIDS / LSHTM workshop (London 14-16 February 2001). Geneva Switzerland WHO, 2001.
 10. Obstetricians R. Management of genital herpes in pregnancy[J]. RCOG Green-top Guideline, 2007.
 11. Bibbins-Domingo K, Grossman D C, Curry S J, et al. Serologic Screening for Genital Herpes Infection: US Preventive Services Task Force Recommendation Statement[J]. Jama, 2016, 316(23):2525.
 12. Mendelson E, Aboudy Y, Smetana Z, et al. Laboratory assessment and diagnosis of congenital viral infections: Rubella, cytomegalovirus (CMV), varicella-zoster virus (VZV), herpes simplex virus (HSV), parvovirus B19 and human immunodeficiency virus (HIV).[J]. Reproductive Toxicology, 2006, 21(4):350-382.
 13. Boscatto Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
 14. Kricka L. Interferences in immunoassays—still a threat. Clin Chem 2000; 46: 1037-1038.
 15. Bjerner J y otros. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
 16. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline-Third Edition.

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Attenzione

Bibliografia

1. Whitley R J, Roizman B. Herpes simplex virus infections[J]. 2001, 357(9267):1513-1518.
2. Centers for Disease Control and Prevention. Sexually Transmitted Infections Treatment Guidelines, 2021. MMWR Recomm Rep 2021;70(4).
3. World Health Organization. WHO guidelines for the treatment of genital herpes simplex virus. Geneva Switzerland WHO, 2016.
4. Sili U, Kaya A, Mert A. Herpes simplex virus encephalitis: Clinical manifestations, diagnosis and outcome in 106 adult patients[J]. Journal of Clinical Virology, 2014, 60(2):112-118.
5. Kaye S, Choudhary A. Herpes simplex keratitis[J]. Progress in Retinal and Eye Research, 2006, 25(4):355-380.
6. Fujie, Xu, Maya, et al. Trends in Herpes Simplex Virus Type 1 and Type 2 Seroprevalence in the United States[J]. Jama, 2006.
7. McQuillan G, Kruszon-Moran D, Flagg EW, Paulose-Ram R. Prevalence of herpes simplex virus type 1 and type 2 in persons aged 14-49: United States, 2015-2016. NCHS Data Brief 2018; 304:1-8.
8. Fleming D T, Mcquillan G M, Johnson R E, et al. Herpes simplex virus type 2 in the United States, 1976 to 1994.[J]. N Engl J Med, 1997, 337(16):1105-1111.
9. World Health Organization. Herpes simplex virus type 2: programmatic and research priorities in developing countries. Report of a WHO / UNAIDS / LSHTM workshop (London 14-16 February 2001). Geneva Switzerland WHO, 2001.

10. Obstetricians R. Management of genital herpes in pregnancy[J]. RCOG Green-top Guideline, 2007.
11. Bibbins-Domingo K, Grossman D C, Curry S J, et al. Serologic Screening for Genital Herpes Infection: US Preventive Services Task Force Recommendation Statement[J]. Jama, 2016, 316(23):2525.
12. Mendelson E, Aboudy Y, Smetana Z, et al. Laboratory assessment and diagnosis of congenital viral infections: Rubella, cytomegalovirus (CMV), varicella-zoster virus (VZV), herpes simplex virus (HSV), parvovirus B19 and human immunodeficiency virus (HIV).[J]. Reproductive Toxicology, 2006, 21(4):350-382.
13. Boscatto LM., Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
14. Kricka L. Interferences in immunoassays—still a threat. Clin Chem 2000; 46: 1037-1038.
15. Bjerner J. et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
16. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline –Third Edition.

© 2022-2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Tutti i diritti riservati.

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Indirizzo: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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10. Obstetricians R. Management of genital herpes in pregnancy[J]. RCOG Green-top Guideline, 2007.
11. Bibbins-Domingo K, Grossman D C, Curry S J, et al. Serologic Screening for Genital Herpes Infection: US Preventive Services Task Force Recommendation Statement[J]. Jama, 2016, 316(23):2525.
12. Mendelson E, Aboudy Y, Smetana Z, et al. Laboratory assessment and diagnosis of congenital viral infections: Rubella, cytomegalovirus (CMV), varicella-zoster virus (VZV), herpes simplex virus (HSV), parvovirus B19 and human immunodeficiency virus (HIV).[J]. Reproductive Toxicology, 2006, 21(4):350-382.
13. Boscatto Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
14. Kricka L. Interferences in immunoassays–still a threat. Clin Chem 2000; 46: 1037-1038.
15. Bjerner J et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
16. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline –Third Edition.

© 2022-2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

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Üretici: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Adres: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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2. Centers for Disease Control and Prevention. Sexually Transmitted Infections Treatment Guidelines, 2021. MMWR Recomm Rep 2021;70(4).
3. World Health Organization. WHO guidelines for the treatment of genital herpes simplex virus. Geneva Switzerland WHO, 2016.
4. Sili U, Kaya A, Mert A. Herpes simplex virus encephalitis: Clinical manifestations, diagnosis and outcome in 106 adult patients[J]. Journal of Clinical Virology, 2014, 60(2):112-118.
5. Kaye S, Choudhary A. Herpes simplex keratitis[J]. Progress in Retinal and Eye Research, 2006, 25(4):355-380.
6. Fujie, Xu, Maya, et al. Trends in Herpes Simplex Virus Type 1 and Type 2 Seroprevalence in the United States[J]. Jama, 2006.
7. McQuillan G, Kruszon-Moran D, Flagg EW, Paulose-Ram R. Prevalence of herpes simplex virus type 1 and type 2 in persons aged 14-49: United States, 2015-2016. NCHS Data Brief 2018; 304:1-8.
8. Fleming D T, Mcquillan G M, Johnson R E, et al. Herpes simplex virus type 2 in the United States, 1976 to 1994.[J]. N Engl J Med, 1997, 337(16):1105-1111.
9. World Health Organization. Herpes simplex virus type 2: programmatic and research priorities in developing countries. Report of a WHO / UNAIDS / LSHTM workshop (London 14-16 February 2001). Geneva Switzerland WHO, 2001.
10. Obstetricians R. Management of genital herpes in pregnancy[J]. RCOG Green-top Guideline, 2007.
11. Bibbins-Domingo K, Grossman D C, Curry S J, et al. Serologic Screening for Genital Herpes Infection: US Preventive Services Task Force Recommendation Statement[J]. Jama, 2016, 316(23):2525.
12. Mendelson E, Aboudy Y, Smetana Z, et al. Laboratory assessment and diagnosis of congenital viral infections: Rubella, cytomegalovirus (CMV), varicella-zoster virus (VZV), herpes simplex virus (HSV), parvovirus B19 and human immunodeficiency virus (HIV).[J]. Reproductive Toxicology, 2006, 21(4):350-382.
13. Boscatto Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
14. Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
15. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.

16. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline-Third Edition.

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Fabricant : Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Représentant en Europe : Shanghai International Holding Corp. GmbH (Europe)