RV IgG

Rubella Virus IgG (CLIA)

INSTRUCTIONS FOR USE

mindray

RV IgG

Rubella Virus IgG (CLIA)

Order Information

Catalog No.	Package Size
105-012518-00	2×50 tests (calibrators not included)
105-012519-00	2×100 tests (calibrators not included)
105-012521-00	2×50 tests (calibrators included)
105-012522-00	2×100 tests (calibrators included)

Intended Use

The CL-series Rubella virus IgG assay is a Chemiluminescent Immunoassay (CLIA) for the quantitative detection of Rubella virus IgG in human serum and plasma. The results of this assay can be used as an aid in the diagnosis of Rubella virus infection and in the assessment of immune status of individuals including pregnant women.

Summary

Rubella virus (RV) is the pathogen of rubella, which belongs to the genus rubella virus in the family of Toga virus. Rubella virus is a single-stranded RNA virus with human being as the only natural host. Rubella is an acute respiratory infectious disease, which is mainly transmitted through respiratory tract. The infection rate is more than 90%. Rubella virus infection generally has mild symptoms and short duration, characterized by fever, rash, lymphadenopathy, etc. Infection confers lifelong immunity^[1].

In the first trimester of pregnancy, Rubella virus may infect the fetus through the placental barrier and cause congenital rubella syndrome (CRS). The proliferation of Rubella virus in fetus leads to continuous systemic infection. By inhibiting cell division, the virus may delay or even arrest the embryo development. Intrauterine death may occur in some fetuses. The overall defect rate of neonates is up to 90%, of which the mild ones show growth retardation, and the severe ones present systemic organ damage, such as cataract, nerve deafness, congenital heart disease, hematological disease, etc. CRS results in 10-20% of neonatal deaths within one year after birth^[2-3].

Detection of Rubella virus specific antibody is indispensable for diagnosis of acute infection and assessment of immune status. Rubella virus IgM antibody usually begins to rise within 1 week after infection and persists for 1-3 months. IgM positive results suggest that patients may have acute infection. IgG antibody generally rises after 2 weeks after infection and persists for years or even decades,

which is a sign of previous infection. For individuals who were not vaccinated against Rubella virus within one month, Rubella virus antibody seroconversion or Rubella virus IgG of the second sample was significantly higher than the first one, which may support the diagnosis of acute Rubella virus infection^[4].

Assay Principles

The CL-series Rubella virus IgG assay is a indirect method to quantitatively detect IgG antibodies to Rubella virus.

In the first step, the sample, sample treatment solution and superparamagnetic microparticles coated with Rubella virus antigens are added into a reaction cuvette. After incubation, Rubella virus IgG antibodies in the sample will bind to Rubella virus antigens coated microparticles. Afterwards, microparticles are magnetically captured while other unbound substances are removed by washing.

In the second step, test diluent solution and ALP labeled anti-human IgG monoclonal antibody are added to the reaction cuvette. After incubation, ALP labeled anti-human IgG monoclonal antibody will form sandwich structure with microparticles captured Rubella virus IgG antibodies. Afterwards, microparticles are magnetically captured while other unbound substances are removed by washing.

Then, the substrate solution is added to the reaction cuvette. It is catalyzed by anti-human IgG antibody -ALP conjugate in the immune-complex retained on the microparticles. The resulting chemiluminescent reaction is measured as relative light units (RLUs) by a photomultiplier in the instrument.

The amount of Rubella virus IgG antibodies present in the sample is proportional to the relative light units (RLUs) generated during the reaction. The Rubella virus IgG antibodies results can be determined via a calibration master curve, which is established on an encoded Master Curve and three-level product calibrators.

Reagent Components

The kit is composed of Ra, Rb, Rc, and Rd, as well as calibrators C0, C1 and C2. Calibrators C0, C1 and C2 are optional (refer to the "Order Information" section). Components in different lots of reagent kit cannot be interchanged, and the detailed information of each component is listed below:

Ra	Paramagnetic microparticles coated with Rubella virus specific antigens in Tris buffer. Minimum concentration: 0.3 mg/mL solids. Tris buffer: 25 mmol/L. Preservative: 0.048% Proclin 300. 3.8 mL/50 tests; 6.6 mL/100 tests.
----	---

Rb	ALP labeled anti-human IgG monoclonal antibody in MES buffer. Minimum concentration: 10ng/mL. MES buffer: 50mmol/L. Preservative: 0.048% Proclin 300. 3.5 mL/50 tests; 6.3 mL/100 tests.
Rc	Sample treatment solution. 25mmol/L Tris buffer. Preservative: 0.048% Proclin 300. 3.5 mL/50 tests; 6.3 mL/100 tests.
Rd	Test diluent solution. 25mmol/L Tris buffer. Preservative: 0.048% Proclin 300. 7.8 mL/50 tests; 14.8 mL/100 tests.
Calibrators C0, C1 and C2 (optional)	IgG antibodies to Rubella virus of different levels (0, ~5 and ~100 IU/mL) in buffer with 0.048% Proclin 300. C0:1×0.6 mL/vial; C1:1×0.5 mL/vial; C2:1×0.5 mL/vial.

Storage and Stability

The unopened reagents (Ra, Rb, Rc and Rd) are stable up to the stated expiration date, when stored at 2-8°C. They can be stored onboard and used for a maximum of 28 days at 2-8°C after opening.

The calibrators C0, C1 and C2 are stable up to the expiration date as indicated on the label when stored unopened at 2-8°C. Once opened in use, the calibrators are stable for 28 days when refrigerated at 2-8°C and stable for 90 days when frozen at -20°C (Frozen only once).

Materials Required but not Provided

Mindray CL-series Chemiluminescence Immunoassay Analyzer

Mindray ToRCH IgG Positive Control:

1×3.0 mL/vial (REF: 105-024832-00),

3×3.0 mL/vial (REF: 105-024833-00),

3×1.0 mL/vial (REF: 105-024834-00),

6×1.0 mL/vial (REF: 105-024835-00).

Mindray ToRCH IgG/IgM Negative Control:

1×3.0 mL/vial (REF: 105-024824-00),

3×3.0 mL/vial (REF: 105-024825-00),

3×1.0 mL/vial (REF: 105-024826-00),

6×1.0 mL/vial (REF: 105-024827-00).

Mindray Substrate Solution, 4×115mL (REF: 105-004274-00).

Mindray Substrate Solution, 4×75 mL (REF: 105-009044-00).

Mindray Sample Diluent, 6×8 mL (REF: 105-004628-00).

Mindray Sample Diluent, 2×30 mL (REF: 105-004276-00).

Mindray Wash Buffer, 1×10 L (REF: 105-004552-00).

Mindray Reaction Cuvette

Instrument System

Mindray CL-series Chemiluminescence Immunoassay Analyzers

Specimen Collection and Preparation**Specimen Types**

- Human serum and plasma (K₂-EDTA, K₃-EDTA, Na₂-EDTA, sodium heparin, lithium heparin and sodium citrate) are recommended for this assay.
- Blood collection tubes from various manufacturers may contain different raw materials and additives which could affect the test results in some cases. Not all available tubes of all manufacturers were tested by Mindray. Each laboratory needs to judge the applicability of the tubes and plasma separation products. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Specimen Conditions

- Do not use:
 - heat-inactivated specimens
 - grossly hemolyzed specimens
 - specimens with obvious microbial contamination
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.

Preparation for Analysis

- Follow blood collection tube manufacturer's recommendations for centrifugation. Centrifuge the specimens after clot formation is complete. Ensure residual fibrin and cellular matter has been removed prior to analysis.
- For optimal results, inspect all samples for bubbles. Remove bubbles with a pipette tip prior to analysis. Specimens must be mixed thoroughly after thawing. Thawed samples should be centrifuged prior to use.
- If the sample is covered with lipid layer after centrifugation, the sample should be transferred to a clean tube and centrifuged before testing. Do not transfer the lipid layer. Handle carefully to prevent cross contamination.

Specimen Storage

- Specimens should be tested as soon as possible after sample collection. If testing is not completed within 48 hours, specimens should be tightly capped and refrigerated at 2-8°C. The sample is stable for 7 days at 2-8°C and 180 days at -20°C or below.
- The samples may be frozen no more than 5 times.
- Samples, controls and calibrators that are frozen or refrigerated must be restored to room temperature before analysis.

Assay Procedure

For optimal performance of this assay, operators should read the related Operator's Manual carefully to get enough information such as operation instructions, sample preservation and management, safety precautions, and maintenance. Prepare all required materials for the assay as well.

Before loading the Rubella IgG reagent kit on the instrument for the first time, invert unopened reagent bottle gently for at least 30 times to resuspend the microparticles, which have settled during shipment or storage. Visually inspect the bottle to ensure the microparticles have been well mixed. If the microparticles remain adhered to the bottle, continue inverting until suspended completely. Avoid foam formation. If the microparticles cannot be homogenized, it is recommended not to use this bottle of reagent. Contact Mindray Customer Service for help. Avoid inverting opened reagent bottle.

This assay requires 6 μ L of sample volume for a single test not including dead volume of the sample container. Additional volume is required when performing additional tests from the same sample. Operators should refer to the Operator's Manual and follow the specific requirements of the assay to determine the minimum sample volume.

Calibration

The calibrators are traceable to WHO 1st International Standard (RUBI-1-94).

The calibration information is stored in the barcodes attached inside the reagent and calibrator pack. When performing the calibration, please scan the barcodes into the system.

Take the vial out of refrigerator and let it stand for no less than 30 minutes at room temperature. Mix contents by gently inverting the vial no less than 10 times before use. Avoid bubble formation. Test C0 three times and C1, C2 twice for calibration. The original calibrator vials may be used if the entire provided volume is used. Because of possible evaporation effects, the calibrator on the analyzer should be measured as soon as possible and not more than 3 calibration procedures per original vial should be performed. After use, close the bottles immediately and store upright at 2-8°C. The original calibrators can also be transferred into empty labeled bottles in aliquots. Use the aliquots immediately or store them immediately at -20°C. One aliquot can be used only once.

A valid calibration is required before any Rubella virus IgG testing.

Renewed calibration is required as follows:

1. The previous calibration was performed 28 days before.
2. Quality control results lie outside of the expected ranges.
3. A reagent kit with a new lot number is used.

For detailed instruction of calibration, refer to the Operator's Manual.

Quality Control

Users can prepare quality controls with clinical samples. Reference ranges can be established according to protocol approved by individual laboratories.

It is recommended that quality controls should be run once every 24 hours if the tests are in use, and following each calibration. The quality control frequency should be adapted to each laboratory's requirements.

Quality control results should be within the acceptable ranges. If a control is out of its specified range, the associated test results are invalid and the samples must be retested. Recalibration may be required. Refer to the Operator's Manual to check up the system. If quality control findings are still outside the defined ranges, please contact Mindray Customer Service for help.

Calculation

The analyzer automatically calculates the RV IgG result of each sample, on the basis of calibration master curve read from the barcode, and the mathematical relationship between the measured responses and the known analyte concentration establishes the calibration curve. The results are shown in the unit of IU/mL.

Specimen Dilution

Samples with RV IgG concentrations above the upper limit can be diluted with Mindray Sample Diluent. The recommended dilution is 1:10 (diluted either automatically by the analyzer or manually). The concentration of the diluted sample must be > 10 IU/mL. After manual dilution, multiply the result by the dilution factor. After automated dilution by the analyzer, the system automatically multiplies the result by the dilution factor when calculating the sample concentration.

Note: Antibodies to Rubella virus are heterogeneous. This may lead to non-linear dilution behavior.

Interpretation of Results

Results	Interpretation
<5.00 IU/mL	Non-reactive
≥5.00 to <10.00 IU/mL	Indeterminate
≥10.00 IU/mL	Reactive

The NCCLS subcommittee on Rubella Serology recommended 10 IU/mL as the cutoff level ^[5].

If the result is indeterminate, a second sample should be collected for testing within an appropriate time (e.g. 2-4 weeks) or other test methods should be used for confirmation.

The results of different detection systems may be inconsistent because of the detection methods and traceability.

The assay results should not be used solely for confirmation or exclusion of Rubella virus infection. Clinical decisions should be made in conjunction with other evidences, such as symptoms, clinical history, results of nucleic acid test, etc.

Limitation of Measurement

A negative RV IgG test result cannot completely exclude the possibility of acute Rubella virus infection; the body may have not produced enough antibody in the early stage of acute infection. If acute infection is suspected, a second sample should be collected for testing within an appropriate period of time (e.g. 2-4 weeks). If the second sample shows seroconversion of primary infection or significant increase of RV IgG, the diagnosis of acute Rubella virus infection is supported. It is recommended that the first sample should be stored at -20°C or below, so that two samples can be tested at the same time and same conditions^[6].

For patients with immunodeficiency disease or receiving immunosuppressive therapy, the antibody test results should be interpreted with caution.

Exogenous interference should be considered in case a patient with Rubella virus IgG positive result has recently received blood transfusion.

Due to the existence of heterophilic antibodies, patients who are often exposed to animals or use animal derived immunoglobulin for treatment may have abnormal detection results^[7].

The samples of patients who have been treated with mouse monoclonal antibody preparation may contain human anti-mouse antibody (HAMA), which may interfere with the test results of kits containing mouse monoclonal antibodies^[8-9]. Mindray chemiluminescence kit contain anti-interference components, which can effectively reduce the impact of HAMA in samples, but the interference of HAMA may still exist in a few samples. Only combined with the clinical information or other laboratory test results can the patient's condition be correctly determined.

Performance Characteristics**Analytical Sensitivity**

Limit of Blank (LoB) = 0.17 IU/mL

Limit of Detection (LoD) = 0.40 IU/mL

Limit of Quantitation (LoQ) = 1.00 IU/mL

The Limit of Blank, Limit of Detection and Limit of Quantitation were determined in accordance with the CLSI EP17-A2 requirements^[10].

The Limit of Blank is the 95th percentile value from $n \geq 60$ measurements of analyte-free samples over several independent series. The Limit of Blank corresponds to the concentration below which analyte-free samples are found with a probability of 95%.

The Limit of Detection is determined based on the Limit of Blank and the standard deviation of low concentration samples. The Limit of Detection corresponds to the lowest analyte concentration which can be detected (value above the Limit of Blank with a probability of 95%).

The LoQ was determined from $n \geq 60$ replicates of low-analyte level samples and is defined as the lowest concentration at which a

TEa of 25% was met.

Linearity

A study was performed based on guidance from CLSI EP06-A2^[11]. A high concentration RV IgG sample (>350 IU/mL) was mixed with a low concentration sample at different ratios, generating a series of dilutions. The Linearity of Mindray RV IgG was demonstrated in the range of 1.00 IU/mL to 350.00 IU/mL, the correlation coefficient $r \geq 0.9900$. The linearity data are summarized in the table below*.

No.	Expected value (IU/mL)	Measured value (IU/mL)	Predicted value (IU/mL)	Predicted deviation
1	0.00	0.00	/	/
2	73.28	76.15	75.27	1.16%
3	146.56	139.87	142.72	-2.00%
4	219.84	210.03	210.17	-0.07%
5	293.12	266.20	277.61	-4.11%
6	366.40	366.40	345.06	6.18%

*Representative data; results in individual laboratories may vary.

Precision

Precision was determined by following CLSI Protocol EP05-A3^[12]. Four serum samples and two controls were tested 2 runs each day in duplicate for 20 days (n=80). Results are shown in the table below*.

Concentration range	The acceptable criteria of Within laboratory precision
≥ 0.17 to < 10.00 IU/mL	$SD \leq 0.8$ IU/mL
≥ 10.00 to < 150.00 IU/mL	$CV \leq 8.0\%$
≥ 150.00 to ≤ 350.00 IU/mL	$CV \leq 12.0\%$

Sample	Mean (IU/mL)	Repeatability		Within laboratory	
		SD (IU/mL)	CV(%)	SD (IU/mL)	CV(%)
S1	4.23	0.09	2.13%	0.20	4.84%
S2	10.46	0.19	1.83%	0.54	5.21%
S3	28.63	0.39	1.38%	1.21	4.23%
S4	257.20	10.89	4.24%	23.95	9.31%
NC	0.00	0.00	N/A**	0.03	N/A
PC	31.13	0.561	1.80%	1.49	4.79%

*Representative data; results in individual laboratories may vary.

** N/A = Not Applicable.

Hook effect

When a sample containing a high concentration of determinand is tested, the hook effect may lead to false negative results. This assay adopts two-step mode to keep the analysis signal at a high value, which could effectively avoid false negative results caused by hook effect.

3 samples with high concentration of RV IgG were detected, all results showed the expected high concentration and there was no misjudgement of results.

Accuracy

Two trueness controls with traceable and defined values were used to verify the accuracy of RV IgG assay. Results are shown in the table below (criterion: relative deviation within $\pm 10\%$ of defined value or absolute deviation within ± 1.00 IU/mL)*.

Sample	Defined Value (IU/mL)	Measured Value (IU/mL)	Absolute Deviation (IU/mL)	Relative Deviation
Defined control L	7.60	7.16	-0.44	-5.79%
Defined control H	37.35	37.66	0.11	0.29%

**Representative data; results in individual laboratories may vary.*

Interference

The test results of Mindray RV IgG are not observed to be interfered with the potential endogenous substances up to the following concentrations (criterion: recovery within $\pm 15\%$ of initial value or absolute deviation within ± 1.50 IU/mL)*.

Substance	Concentration
Hemoglobin	1000 mg/dL
Bilirubin	40 mg/dL
Triglyceride	37 mmol/L
Total protein	12 g/dL
Albumin	6 g/dL
γ globulin	6 g/dL
Total IgM	4 g/L
Biotin	3500 ng/mL
Rheumatoid Factor(RF)	1200 IU/mL

**Representative data; results in individual laboratories may vary.*

The test results of Mindray RV IgG are not observed to be interfered with the potential exogenous substances up to the following concentrations (criterion: recovery within $\pm 15\%$ of initial value or absolute deviation within ± 1.50 IU/mL)*.

Substance	Concentration
Acetaminophen	200 mg/L
Amoxicilin	700 mg/L
Cefaclor	500 mg/L
Vitamin C	250 mg/L
Aspirin	1000 mg/L
Ibuprofen	500 mg/L

**Representative data; results in individual laboratories may vary.*

Cross-reactivity

Studies were performed to evaluate the effect of potential cross reacting diseases on Mindray RV IgG assay. A total of 300 specimens were tested by Mindray RV IgG and a comparator RV IgG kit.

Substance	n
Toxoplasma gondii	8
Cytomegalo virus	14
Herpes simplex virus infection	5
Morbillivirus measles virus	4
Coxsackie virus-16	3
Enterovirus-71	3
Dengue fever virus	3
Chikungunya fever virus	4
Zika Virus	6
Epstein-Barr virus	17
Parvovirus B19	4
HBV	54
HCV	5
HIV	11
TP	15
VZV-virus	6
Influenza A virus	7
Influenza B virus	6
Parainfluenza virus	7
Mycoplasma pneumoniae	17
AMA	3
ANA	45
RF	37
HAMA	4
SLE	4
Rubella IgM	5
Treponema pallidum	3

The results showed that the overall agreement of Mindray RV IgG and commercial RV IgG kits is 97.67%.

Method comparison

1073 clinical specimens from pregnant women, women of childbearing age, hospitalized individuals, neonates and patients with corresponding clinical symptoms were tested to evaluate the coincidence rate of Mindray RV IgG and a mainstream commercial RV IgG kit. The positive coincidence rate was 99.54% (650/653, 95%CI 98.66%-99.84%); the negative coincidence rate was 99.45% (364/366, 95%CI 98.03%-99.85%); and the total coincidence rate was 99.51% (1014/1019, 95%CI 98.86%-99.79%). Specimens with inconsistent test results were confirmed by another commercial RV

IgG kit.

Mindray RV IgG	Comparator RV IgG		Total
	Reactive	Non-reactive	
Reactive	650	2*	652
Non-reactive	3**	364	367
Indeterminate	0	54	54
Total	653	420	1073

*: Of 2 specimens with reactive on Mindray RV IgG and non-reactive on comparator RV IgG, 2 specimens were found reactive on another commercial RV IgG kit.

** : Of 3 specimens with non-reactive on Mindray RV IgG and reactive on comparator RV IgG, 3 specimens were found non-reactive on another commercial RV IgG kit.

Seroconversion panels

To evaluate the sensitivity of seroconversion, 4 commercial seroconversion panels were tested by Mindray RV IgG and a mainstream commercial RV IgG kit, the results are shown in the table below.

Panel ID	Mindray RV IgG		Comparator RV IgG	
	Days*	PR**	Days	PR
SCP-RUB-0001	21	60% (9/15)	21	60% (9/15)
SCP-RUB-0002	19	53.85% (7/13)	19	53.85% (7/13)
SCP-RUB-0003	19	75% (15/20)	19	75% (15/20)
SCP-RUB-0004	21	80% (12/15)	21	80% (12/15)

*: Days: time from the first blood collection after positive conversion.

** : PR: positive rate.

Warnings and Precautions

1. For *in vitro* diagnostic use only. For laboratory professional use.
2. Follow all the rules in handling laboratory reagents and take necessary safety precautions.
3. Due to the differences in methodology and antibody specificity, test results of the same sample may be different when using reagent kits from different manufacturers on Mindray system, or using Mindray reagent kits on other systems.

4. Do not use reagent kits beyond the expiration date.
5. Do not use reagents mixed from different reagent lots.
6. Please confirm the integrity of the package before use. Do not use the reagents with damaged packages. The results can't be assured when the reagents are stored at inappropriate condition.
7. If the reagents are opened unintentionally before use, they shall be used as soon as possible.
8. Wear clean gloves to avoid contaminating the reagents.
9. Always keep the reagent pack in the upright position to ensure no microparticle has been lost prior to use.
10. The calibrator was tested with CE-marked methods and shown to be negative for antibodies to HIV and HCV and nonreactive for HBsAg. However, as no testing method can rule out the potential risk of infection with absolute certainty, so this material should be handled as a patient specimen.
11. Reagent pack opened for more than 28 days is not recommended for use.
12. Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
13. All the specimen and reaction wastes should be considered potentially biohazard. Specimens and reaction wastes should be handled in accordance with the local regulations and guidelines.
14. The Material Safety Data Sheet (MSDS) is available upon request.
15. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the local competent authority.
16. Instability or deterioration should be suspected if there are visible signs of leakage, turbidity, precipitates or microbial growth.
17. This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

**Warning**

H317 May cause an allergic skin reaction.

H412 Harmful to aquatic life with long lasting effects.

Prevention:

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P272 Contaminated work clothing should not be allowed out of the workplace.

P273 Avoid release to the environment.

P280 Wear protective gloves/protective clothing/eye

protection/face protection.

Response:

P302 + P352 IF ON SKIN: Wash with plenty of water.

P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.

P362 + P364 Take off contaminated clothing and wash it before reuse.

Disposal:

P501 Dispose of contents/ container in accordance with local regulation.

Graphical Symbols

In vitro
diagnostic
medical
device

Authorized
Representative
in the European
Community

Consult
instructions
for use

European
Conformity

Catalogue
number

Temperature
limit

Manufacturer

Use-by
date

Biological
risks

Batch code

This way up

Unique
device
identifier

Caution

Bibliography

1. Banatvala JE, Brown DWG. Rubella. Lancet 2004; 363:1127-1137.
2. Anonymous. Rubella vaccines: WHO position paper.[J]. Wkly Epidemiol Rec, 2011, 86(29):301-316.
3. Dimech W, Grangeot-Keros L, Vauloup-Fellous C. Standardization of Assays That Detect Anti-Rubella Virus

- IgG Antibodies[J]. Clinical Microbiology Reviews, 2016, 29(1):163.
4. Plotkin S A. Infectious Diseases of the Fetus and Newborn || Rubella[J]. 2011:861-898.
 5. Detection and Quantitation of Rubella IgG Antibody: Evaluation and Performance Criteria for Multiple Component Test Products, Specimen Handling, and Use of Test Products in the Clinical Laboratory; Approved Guideline. NCCLS document I/LA6-A (ISBN) 1-56238-335-3. NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 1997.
 6. Sensini A. Toxoplasma gondii infection in pregnancy: opportunities and pitfalls of serological diagnosis[J]. Clinical Microbiology and Infection, 2010, 12(6):504-512.
 7. Boscatto Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
 8. Kricka L. Interferences in immunoassays—still a threat. Clin Chem 2000; 46: 1037-1038.
 9. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
 10. CLSI. EP17-A2: Vol. 32, No. 8, Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition.
 11. CLSI. EP06-A2: Vol. 40, No. 16, Evaluation of Linearity of Quantitative Measurement Procedures; Approved Guideline-Second Edition.
 12. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline –Third Edition.

© 2022-2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
All rights Reserved.

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

E-mail Address: service@mindray.com

Website: www.mindray.com

Tel: +86-755-81888998

Fax: +86-755-26582680

EC-Representative: Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

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3. Dimech W, Grangeot-Keros L, Vauloup-Fellous C . Standardization of Assays That Detect Anti-Rubella Virus IgG Antibodies[J]. Clinical Microbiology Reviews, 2016, 29(1):163.
4. Plotkin S A. Infectious Diseases of the Fetus and Newborn || Rubella[J]. 2011:861-898.
5. Detection and Quantitation of Rubella IgG Antibody: Evaluation and Performance Criteria for Multiple Component Test Products, Specimen Handling, and Use of Test Products in the Clinical Laboratory; Approved Guideline. NCCLS document I/LA6-A (ISBN) 1-56238-335-3. NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 1997.
6. Sensini A . Toxoplasma gondii infection in pregnancy: opportunities and pitfalls of serological diagnosis[J]. Clinical Microbiology and Infection, 2010, 12(6):504-512.
7. Boscatto Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
8. Kricka L. Interferences in immunoassays—still a threat. Clin Chem 2000; 46: 1037-1038.
9. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
10. CLSI. EP17-A2: Vol. 32, No. 8, Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition.
11. CLSI. EP06-A2: Vol. 40, No. 16, Evaluation of Linearity of Quantitative Measurement Procedures; Approved Guideline-Second Edition.
12. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline –Third Edition.

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Адрес: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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5. Detection and Quantitation of Rubella IgG Antibody: Evaluation and Performance Criteria for Multiple Component Test Products, Specimen Handling, and Use of Test Products in the Clinical Laboratory; Approved Guideline. NCCLS document I/LA6-A (ISBN) 1-56238-335-3. NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 1997.
6. Sensini A . Toxoplasma gondii infection in pregnancy: opportunities and pitfalls of serological diagnosis[J]. Clinical Microbiology and Infection, 2010, 12(6):504-512.
7. Boscatto Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
8. Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
9. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
10. CLSI. EP17-A2: Vol. 32, No. 8, Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition.
11. CLSI. EP06-A2: Vol. 40, No. 16, Evaluation of Linearity of Quantitative Measurement Procedures; Approved Guideline-Second Edition.
12. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline -Third Edition.

© 2022-2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

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Fabricante: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Endereço: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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5. Detection and Quantitation of Rubella IgG Antibody: Evaluation and Performance Criteria for Multiple Component Test Products, Specimen Handling, and Use of Test Products in the Clinical Laboratory; Approved Guideline. NCCLS document I/LA6-A (ISBN) 1-56238-335-3. NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 1997.
6. Sensini A . Toxoplasma gondii infection in pregnancy: opportunities and pitfalls of serological diagnosis[J]. Clinical Microbiology and Infection, 2010, 12(6):504-512.
7. Boscatto Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
8. Kricka L. Interferences in immunoassays—still a threat. Clin Chem 2000; 46: 1037-1038.
9. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
10. CLSI. EP17-A2: Vol. 32, No. 8, Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition.
11. CLSI. EP06-A2: Vol. 40, No. 16, Evaluation of Linearity of Quantitative Measurement Procedures; Approved Guideline-Second Edition.
12. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline –Third Edition.

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5. Detection and Quantitation of Rubella IgG Antibody: Evaluation and Performance Criteria for Multiple Component Test Products, Specimen Handling, and Use of Test Products in the Clinical Laboratory; Approved Guideline. NCCLS document I/LA6-A (ISBN) 1-56238-335-3. NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 1997.
6. Sensini A . Toxoplasma gondii infection in pregnancy: opportunities and pitfalls of serological diagnosis[J]. Clinical Microbiology and Infection, 2010, 12(6):504-512.
7. Boscatto Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
8. Kricka L. Interferences in immunoassays—still a threat. Clin Chem 2000; 46: 1037–1038.
9. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
10. CLSI. EP17-A2: Vol. 32, No. 8, Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition.
11. CLSI. EP06-A2: Vol. 40, No. 16, Evaluation of Linearity of Quantitative Measurement Procedures; Approved Guideline-Second Edition.
12. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline –Third Edition.

© 2022-2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Tutti i diritti riservati.

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Indirizzo: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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6. Sensini A . Toxoplasma gondii infection in pregnancy: opportunities and pitfalls of serological diagnosis[J]. Clinical Microbiology and Infection, 2010, 12(6):504-512.
7. Boscatto Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
8. Kricka L. Interferences in immunoassays–still a threat. Clin Chem 2000; 46: 1037–1038.
9. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
10. CLSI. EP17-A2: Vol. 32, No. 8, Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition.
11. CLSI. EP06-A2: Vol. 40, No. 16, Evaluation of Linearity of Quantitative Measurement Procedures; Approved Guideline-Second Edition.
12. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline –Third Edition.

© 2022-2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

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Üretici: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Adres: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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- IgG Antibodies[J]. Clinical Microbiology Reviews, 2016, 29(1):163.
4. Plotkin S A. Infectious Diseases of the Fetus and Newborn || Rubella[J]. 2011:861-898.
 5. Detection and Quantitation of Rubella IgG Antibody: Evaluation and Performance Criteria for Multiple Component Test Products, Specimen Handling, and Use of Test Products in the Clinical Laboratory; Approved Guideline. NCCLS document I/LA6-A (ISBN) 1-56238-335-3. NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 1997.
 6. Sensini A . Toxoplasma gondii infection in pregnancy: opportunities and pitfalls of serological diagnosis[J]. Clinical Microbiology and Infection, 2010, 12(6):504-512.
 7. Boscatto Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
 8. Kricka L. Interferences in immunoassays—still a threat. Clin Chem 2000; 46: 1037-1038.
 9. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
 10. CLSI. EP17-A2: Vol. 32, No. 8, Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition.
 11. CLSI. EP06-A2: Vol. 40, No. 16, Evaluation of Linearity of Quantitative Measurement Procedures; Approved Guideline-Second Edition.
 12. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline –Third Edition.

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1. Banatvala JE, Brown DWG. Rubella. Lancet 2004; 363:1127-1137.
2. Anonymous. Rubella vaccines: WHO position paper.[J]. Wkly Epidemiol Rec, 2011, 86(29):301-316.
3. Dimech W, Grangeot-Keros L, Vauloup-Fellous C. Standardization of Assays That Detect Anti-Rubella Virus IgG Antibodies[J]. Clinical Microbiology Reviews, 2016, 29(1):163.
4. Plotkin S A. Infectious Diseases of the Fetus and Newborn || Rubella[J]. 2011:861-898.
5. Detection and Quantitation of Rubella IgG Antibody: Evaluation and Performance Criteria for Multiple Component Test Products, Specimen Handling, and Use of Test Products in the Clinical Laboratory; Approved Guideline. NCCLS document I/LA6-A (ISBN) 1-56238-335-3. NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 1997.
6. Sensini A. Toxoplasma gondii infection in pregnancy: opportunities and pitfalls of serological diagnosis[J]. Clinical Microbiology and Infection, 2010, 12(6):504-512.
7. Boscatto Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
8. Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
9. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
10. CLSI. EP17-A2: Vol. 32, No. 8, Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition.
11. CLSI. EP06-A2: Vol. 40, No. 16, Evaluation of Linearity of Quantitative Measurement Procedures; Approved Guideline-Second Edition.
12. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline -Third Edition.

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«Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

1203 Nanhuan Avenue, Guangming District, 518106 Shenzhen, P. R. China

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