Toxo IgG

Toxoplasma gondii IgG (CLIA)

INSTRUCTIONS FOR USE

mindray

Toxo IgG

Toxoplasma gondii IgG (CLIA)

Order Information

Catalog No.	Package Size
105-012524-00	2×50 tests (calibrators not included)
105-012525-00	2×100 tests (calibrators not included)
105-012527-00	2×50 tests (calibrators included)
105-012528-00	2×100 tests (calibrators included)

Intended Use

The CL-series Toxoplasma gondii IgG assay is a Chemiluminescent Immunoassay (CLIA) for the quantitative detection of Toxoplasma gondii specific IgG in human serum and plasma. The results of this assay can be used as an aid in the diagnosis of Toxoplasma gondii infection and in the assessment of immune status of individuals including pregnant women.

Summary

Toxoplasmosis is a world spread zoonotic infectious caused by Toxoplasma gondii. The infection can be transmitted through oral intake of food containing Toxoplasma gondii, direct contact with soil or water contaminated by oocysts, or vertical transmission of fetus through placenta^[1].

There are no obvious clinical symptoms and signs in healthy individuals infected with Toxoplasma gondii. However, in patients with immune deficiency or suppression, this parasite may invade various organs and lead to serious diseases, such as Toxoplasmosis encephalopathy, ophthalmopathy and liver damage.

If the infection occurs in pregnancy, Toxoplasma gondii can infect the fetus through placenta. In general, the later the pregnant period when pregnant women are infected, the higher the fetal infection rate, and the lighter clinical symptoms. In the first trimester, the fetal infection rate is 10-20%, which often causes spontaneous abortion. In the second trimester, the fetal infection rate increases and most of the fetuses suffer from stillbirth, premature delivery or severe damage of central nervous system. In the third trimester, the fetal infection rate increases to 60%-90%; although most fetuses have no obvious symptoms at birth, they may gradually develop chorioretinitis, severe visual impairment, nervous system growth retardation and other sequelae months or years after birth^[2-6].

Toxoplasmosis is commonly diagnosed by detecting anti

Toxoplasma specific IgG and IgM antibody. Toxo IgM begins to rise after one week of infection and reaches the peak within one month. The level of Toxo IgM reflects the recent or current Toxoplasma gondii infection. Toxo IgG begins to rise after 1 to 3 weeks of infection, reaches the peak after 2 to 3 months, then slowly decreases and persists for a lifetime. Toxo IgG detection can be used to evaluate the serological status of patients [7].

Assay Principles

The CL-series Toxo IgG assay is an indirect assay to quantitatively detect IgG antibodies to Toxoplasma gondii.

In the first step, sample, sample treatment solution and paramagnetic microparticles coated with Toxoplasma gondii antigens are added into a reaction cuvette. After incubation, Toxo IgG antibodies in the sample will be bound to Toxoplasma gondii antigens coated microparticles. Afterwards, microparticles are magnetically captured while other unbound substances are removed by washing.

In the second step, test diluent solution and ALP labeled anti-human IgG monoclonal antibody are added to the reaction cuvette for incubation. After incubation, ALP labeled anti-human IgG monoclonal antibody will form sandwich structure with microparticles captured Toxo IgG antibodies. Afterwards, microparticles are magnetically captured while other unbound substances are removed by washing.

Then, the substrate solution is added to the reaction cuvette. It is catalyzed by anti-human IgG antibody -ALP conjugate in the immune-complex retained on the microparticles. The resulting chemiluminescent reaction is measured as relative light units (RLUs) by a photomultiplier in the instrument.

The amount of Toxo IgG antibodies present in the sample is proportional to the relative light units (RLUs) generated during the reaction. The Toxo IgG antibodies results can be determined via a calibration curve, which is established on an encoded Calibration Master Curve and three-level product calibrators.

Reagent Components

The kit is composed of Ra, Rb, Rc, and Rd, as well as calibrators C0, C1 and C2. Calibrators C0, C1 and C2 are optional (refer to the "Order Information" section). Components in different lots of reagent kit cannot be interchanged. The detailed information of each component is listed below:

Ra	Paramagnetic microparticles coated with Toxoplasma gondii specific antigens. Minimum concentration: 0.3 mg/mL solids. Tris buffer: 50 mmol/L. Preservative: 0.048% Proclin 300. 3.8 mL/50 tests; 6.6 mL/100 tests.
----	--

Rb	ALP labeled anti-human IgG monoclonal antibody in MES buffer. Minimum concentration: 10 ng/mL. MES buffer: 50 mmol/L. Preservative: 0.048% Proclin 300. 5.3 mL/50 tests; 9.8 mL/100 tests.
Rc	Sample treatment solution 25 mmol/L Tris buffer. Preservative: 0.048% Proclin 300. 3.5 mL/50 tests; 6.3 mL/100 tests.
Rd	Test diluent solution 25 mmol/L Tris buffer. Preservative: 0.048% Proclin 300. 3.5 mL/50 tests; 6.3 mL/100 tests.
Calibrators C0, C1 and C2 (optional)	IgG antibodies to <i>Toxoplasma gondii</i> of different levels (0, ~6 and ~50 IU/mL) in buffer with 0.048% ProClin 300. C0:1×0.6 mL/vial; C1:1×0.5 mL/vial; C2:1×0.5 mL/vial.

Storage and Stability

The unopened Toxo IgG reagents (Ra, Rb, Rc and Rd) are stable up to the stated expiration date when stored at 2-8°C. They can be stored onboard and used for a maximum of 28 days at 2-8°C after opening.

The calibrators C0, C1 and C2 are stable up to the expiration date as indicated on the label when stored unopened at 2-8°C. Once opened in use, the calibrators are stable for 28 days when refrigerated at 2-8°C and stable for 90 days when frozen at -20°C (Frozen only once).

Materials Required but not Provided

Mindray CL-series Chemiluminescence Immunoassay Analyzer

Mindray ToRCH IgG Positive Control:

1×3.0 mL/vial (REF: 105-024832-00),

3×3.0 mL/vial (REF: 105-024833-00),

3×1.0 mL/vial (REF: 105-024834-00),

6×1.0 mL/vial (REF: 105-024835-00).

Mindray ToRCH IgG/IgM Negative Control:

1×3.0 mL/vial (REF: 105-024824-00),

3×3.0 mL/vial (REF: 105-024825-00),

3×1.0 mL/vial (REF: 105-024826-00),

6×1.0 mL/vial (REF: 105-024827-00).

Mindray Substrate Solution, 4×115mL (REF: 105-004274-00).

Mindray Substrate Solution, 4×75 mL (REF: 105-009044-00).

Mindray Sample Diluent, 6×8 mL (REF: 105-004628-00).

Mindray Sample Diluent, 2×30 mL (REF: 105-004276-00).

Mindray Wash Buffer, 1×10 L (REF: 105-004552-00).

Mindray Reaction Cuvette

Instrument System

Mindray CL-series Chemiluminescence Immunoassay Analyzers

Specimen Collection and Preparation**Specimen Types**

- Human serum and plasma (K₂-EDTA, K₃-EDTA, Na₂-EDTA, sodium heparin, lithium heparin and sodium citrate) are recommended for this assay. The average measured value of sodium citrate samples may be lower ($\leq 15\%$) than that of serum samples.
- Blood collection tubes from various manufacturers may contain different raw materials and additives which could affect the test results in some cases. Not all available tubes of all manufacturers were tested by Mindray. Each laboratory needs to judge the applicability of the tubes and plasma separation products. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Specimen Conditions

- Do not use:
 - heat-inactivated specimens
 - grossly hemolyzed specimens
 - specimens with obvious microbial contamination
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.

Preparation for Analysis

- Follow blood collection tube manufacturer's recommendations for centrifugation. Centrifuge the specimens after clot formation is complete. Ensure residual fibrin and cellular matter has been removed prior to analysis.
- For optimal results, inspect all samples for bubbles. Remove bubbles with a pipette tip prior to analysis. Specimens must be mixed thoroughly after thawing. Thawed samples should be centrifuged prior to use.
- If the sample was covered with lipid layer after centrifugation, the sample should be transferred to a clean tube and centrifuged before testing. Do not transfer the lipid layer. Handle carefully to prevent cross contamination.

Specimen Storage

- Specimens should be tested as soon as possible after sample collection. If testing is not completed within 48 hours, specimens should be tightly capped and refrigerated at 2-8°C. The sample is stable for 7 days at 2-8°C and 180 days at -20°C or below.
- The samples may be frozen no more than 5 times.
- Frozen or refrigerated samples, controls and calibrators must

be restored to room temperature before analysis.

Assay Procedure

For optimal performance of this assay, operators should read the related Operator's Manual carefully to get enough information such as operation instructions, sample preservation and management, safety precautions, and maintenance. Prepare all required materials for the assay as well.

Before loading the Toxo IgG reagent kit on the instrument for the first time, invert unopened reagent bottle gently for at least 30 times to resuspend the microparticles, which have settled during shipment or storage. Visually inspect the bottle to ensure the microparticles have been well mixed. If the microparticles remain adhered to the bottle, continue inverting until suspended completely. Avoid foam formation. If the microparticles cannot be homogenized, it is recommended not to use this bottle of reagent. Contact Mindray Customer Service for help. Avoid inverting opened reagent bottle.

This assay requires 6 μL of sample volume for a single test not including dead volume of the sample container. Additional volume is required when performing additional tests from the same sample. Operators should refer to the Operator's Manual and follow the specific requirements of the assay to determine the minimum sample volume.

Calibration

The calibrators are traceable to WHO 1st International Standard (Anti-Toxo IgG 01/600).

The calibration information is stored in the barcodes attached inside the reagent and calibrator pack. When performing the calibration, please scan the barcodes into the system.

Take the vial out of refrigerator and let it stand for no less than 30 minutes at room temperature, mix contents by gently inverting the vial no less than 10 times before use. Avoid bubble formation. Test C0 three times and C1, C2 twice for calibration. The original calibrator vials may be used if the entire provided volume is used. Because of possible evaporation effects, the calibrator on the analyzer should be measured as soon as possible and not more than 3 calibration procedures per original vial should be performed. After use, close the bottles immediately and store upright at 2-8°C. The original calibrators can also be transferred into empty labeled bottles in aliquots. Use the aliquots immediately or store them immediately at -20°C. One aliquot can be used only once.

A valid calibration is required before any Toxo IgG testing.

Renewed calibration is required as follows:

1. The previous calibration was performed 28 days before.
2. Quality controls lie outside of the expected ranges.
3. A reagent kit with a new lot number is used.

For detailed instructions of calibration, refer to the Operator's

Manual.

Quality Control

Users can prepare quality controls with clinical samples. Reference ranges can be established according to protocol approved by individual laboratories.

It is recommended that quality controls should be run once every 24 hours if the tests are in use, and following each calibration. The quality control frequency should be adapted to each laboratory's requirements.

Quality control results should be within the acceptable ranges. If a control is out of its specified range, the associated test results are invalid and the samples must be retested. Recalibration may be required. Refer to the Operator's Manual to check up the system. If quality control findings are still outside the defined ranges, please contact Mindray Customer Service for help.

Calculation

The analyzer automatically calculates the Toxo IgG result of each sample on the basis of calibration master curve read from the barcode, and the mathematical relationship between the measured responses and the known analyte concentration establishes the calibration curve. The results are shown in the unit of IU/mL.

Specimen Dilution

Samples with Toxo IgG concentrations above the upper limit can be diluted with Mindray Sample Diluent. The recommended dilution is 1:10 (diluted either automatically by the analyzer or manually). The concentration of the diluted sample must be > 10 IU/mL. After manual dilution, multiply the result by the dilution factor. After automated dilution by the analyzer, the system automatically multiplies the result by the dilution factor when calculating the sample concentration.

Note: Antibodies to *Toxoplasma gondii* are heterogeneous. This may lead to non-linear dilution behavior.

Interpretation of Results

Specimens with results < 1.60 IU/mL are considered to be non-reactive for IgG antibodies to *Toxoplasma gondii*.

Specimens with results from 1.60 IU/mL to < 2.00 IU/mL are considered to be indeterminate for IgG antibodies to *Toxoplasma gondii*. A second sample should be collected for testing within an appropriate time (e.g. 2-3 weeks) or other test methods should be used for confirmation.

Specimens with results ≥ 2.00 IU/mL are considered to be reactive for IgG antibodies to *Toxoplasma gondii*.

The assay results should not be used solely for confirmation or exclusion of Toxo infection. Clinical decisions should be made in conjunction with other evidences, such as symptoms, clinical history, results of nucleic acid test, etc.

Limitation of Measurement

A negative Toxo IgG test result cannot completely exclude the possibility of acute *Toxoplasma gondii* infection. The body may have not produced enough antibody in the early stage of acute infection. If acute infection is suspected, a second sample should be collected for testing within an appropriate time (e.g. 2-3 weeks). If the second sample shows seroconversion of primary infection or significant increase of Toxo IgG, the diagnosis of acute *Toxoplasma gondii* infection is supported. It is recommended that the first sample should be stored at -20°C or below, so that two samples can be tested at the same time and same conditions^[8-9].

For patients with immunodeficiency disease or receiving immunosuppressive therapy, the antibody test results should be interpreted with caution.

Exogenous interference should be considered in case a patient with Toxo IgG positive result has recently received blood transfusion.

Due to the existence of heterophilic antibodies, patients who are often exposed to animals or use animal derived immunoglobulin for treatment may have abnormal detection results^[10].

The samples of patients who have been treated with mouse monoclonal antibody preparation may contain human anti mouse antibody (HAMA), which may interfere with the test results of kits containing mouse monoclonal antibodies^[11-12]. Mindray chemiluminescence kit contains anti-interference components, which can effectively reduce the impact of HAMA in samples, but the interference of HAMA may still exist in a few samples. Only combined with the clinical information or other laboratory test results can the patient's condition be correctly determined.

Performance Characteristics**Analytical Sensitivity**

Limit of Blank (LoB) = 0.13 IU/mL

Limit of Detection (LoD) = 0.30 IU/mL

Limit of Quantitation (LoQ) = 0.50 IU/mL

The Limit of Blank, Limit of Detection and Limit of Quantitation were determined in accordance with the CLSI EP17-A2 requirements^[13].

The Limit of Blank is the 95th percentile value from $n \geq 60$ measurements of analyte-free samples over several independent series. The Limit of Blank corresponds to the concentration below which analyte-free samples are found with a probability of 95%.

The Limit of Detection is determined based on the Limit of Blank and the standard deviation of low concentration samples. The Limit of Detection corresponds to the lowest analyte concentration which can be detected (value above the Limit of Blank with a probability of 95%).

The LoQ was determined from $n \geq 60$ replicates of low-analyte level samples and is defined as the lowest concentration at which a TEa of 25% was met.

Linearity

Linearity study was performed based on guidance from CLSI EP06-A^[14]. A high concentration Toxo IgG sample (>200 IU/mL) was mixed with a low concentration sample at different ratios, generating a series of dilutions. The Linearity of Toxo IgG was demonstrated in the range of 0.50 to 200.00 IU/mL, the correlation coefficient $r \geq 0.9900$. The linearity data are summarized in the table below*.

No.	Expected value (IU/mL)	Measured value (IU/mL)	Predicted value (IU/mL)	Predicted deviation
1	0.00	0.00	/	/
2	41.50	37.70	37.41	0.78%
3	83.01	77.40	78.95	-1.96%
4	124.51	120.23	120.49	-0.22%
5	166.02	161.56	162.04	-0.30%
6	207.52	207.52	203.58	1.94%

*Representative data; results in individual laboratories may vary.

Precision

Precision was determined by following CLSI Protocol EP05-A3^[15]. Three serum samples and two controls were tested 2 runs each day in duplicate for totally 20 days (n=80). Results are shown in the table below*.

Sample	Mean (IU/mL)	Repeatability		Within-Laboratory	
		SD (IU/mL)	CV(%)	SD (IU/mL)	CV(%)
S1	0.22	0.03	N/A**	0.03	N/A**
S2	6.28	0.17	2.70%	0.29	4.70%
S3	143.89	2.74	1.90%	5.53	3.80%
NC	0.00	0.01	N/A	0.01	N/A
PC	6.11	0.12	1.90%	0.21	3.40%

*Representative data; results in individual laboratories may vary.

** N/A = Not Applicable.

Hook effect

When a sample containing a high concentration of determinand is tested, the hook effect may lead to false negative results. This assay adopts two-steps mode to keep the analysis signal at a high value, which can effectively avoid false negative caused by hook effect.

3 samples with high concentration of Toxo IgG were detected, all results showed the expected high concentration and there was no misjudgement of results.

Accuracy

Two trueness controls with traceable and defined values were tested to verify the accuracy of Toxo IgG assay. Results are shown in the table below (criterion: relative deviation within $\pm 10\%$ of

defined value or absolute deviation within ± 0.20 IU/mL)*.

Sample	Defined Value (IU/mL)	Measured Value (IU/mL)	Absolute Deviation (IU/mL)	Relative Deviation
Defined control L	1.52	1.43	0.09	6.23%
Defined control H	38.92	35.95	2.97	8.26%

**Representative data; results in individual laboratories may vary.*

Interference

The test results of Mindray Toxo IgG are not observed to be interfered with the potential endogenous substances up to the following concentrations (criterion: recovery within $\pm 15\%$ of initial value or absolute deviation within ± 0.30 IU/mL)*.

Substance	Concentration
Hemoglobin	1000 mg/dL
Bilirubin	40 mg/dL
Triglyceride	37 mmol/L
Total protein	12 g/dL
Albumin	6 g/dL
γ -globulin	6 g/dL
Total IgM	1 g/dL
Biotin	3500 ng/mL
Rheumatoid Factor(RF)	1200 IU/mL

**Representative data; results in individual laboratories may vary.*

The test results of Mindray Toxo IgG are not observed to be interfered with the potential exogenous substances up to the following concentrations (criterion: recovery within $\pm 15\%$ of initial value or absolute deviation within ± 0.30 IU/mL)*.

Substance	Concentration
Acetaminophen	200 mg/L
Amoxicillin	700 mg/L
Cefaclor	500 mg/L
Vitamin C	250 mg/L
Aspirin	1000 mg/L
Ibuprofen	500 mg/L
Acetylspiramycin	3000 mg/L
Sulfadiazine	2500 mg/L
Pyrimethamine	500 mg/L
Leucovorin	3 mg/L

**Representative data; results in individual laboratories may vary.*

Cross-reactivity

Studies were performed to evaluate the effect of potential cross

reacting diseases on Toxo IgG assay. A total of 250 specimens were tested by Mindray Toxo IgG and a comparator Toxo IgG kit.

Substance	n
High Concentration Toxo IgM	2
RV	4
CMV	8
HSV	3
E-B virus	23
Parvovirus B19	11
Hepatitis A virus	8
Hepatitis B virus	32
Varicella zoster virus	11
Influenza A virus	6
Influenza B virus	6
Parainfluenza virus	3
Mycoplasma pneumoniae	13
HIV	8
AMA	12
RF	45
ANA	52
HAMA	2
SLE	1

The results showed that the overall agreement of Mindray Toxo IgG and commercial Toxo IgG kit is 95.20%.

Method comparison

253 clinical specimens from pregnant women, women of childbearing age, hospitalized individuals, neonates and patients with corresponding clinical symptoms were tested to evaluate the coincidence rate of Mindray Toxo IgG and a mainstream commercial Toxo IgG kit. The positive coincidence rate was 96.72% (118/122, 95%CI 91.87%-98.72%); the negative coincidence rate was 98.39% (122/124, 95%CI 94.31%-99.56%); the total coincidence rate was 97.56% (240/246, 95%CI 94.78%-98.88%). Specimens with inconsistent test results were confirmed by another commercial Toxo IgG kit.

Mindray Toxo IgG	Comparator Toxo IgG			Total
	Reactive	Non-reactive	Indeterminate	
Reactive	118	2*	2	122
Non-reactive	4**	122	3	129
Indeterminate	0	0	2	2
Total	122	124	7	253

*: Of 2 specimens with reactive on Mindray Toxo IgG and non-reactive on comparator Toxo IgG, 2 specimens were found reactive on another commercial Toxo IgG kit.

****:** Of 4 specimens with non-reactive on Mindray Toxo IgG and reactive on comparator Toxo IgG, 2 specimens were found non-reactive on another commercial Toxo IgG kit.

Warnings and Precautions

1. For *in vitro* diagnostic use only. For laboratory professional use.
2. Follow all the rules in handling laboratory reagents and take necessary safety precautions.
3. Due to the differences in methodology and antibody specificity, test results of the same sample may be different when using reagent kits from different manufacturers on Mindray system, or using Mindray reagent kits on other systems.
4. Do not use reagent kits beyond the expiration date.
5. Do not use reagents mixed from different reagent lots.
6. Please confirm the integrity of the package before use. Do not use the reagents with damaged packages. The results can't be assured when the reagents are stored at inappropriate condition.
7. If the reagents are opened unintentionally before use, they shall be used as soon as possible.
8. Always keep the reagent pack in the upright position to ensure no microparticle has been lost prior to use.
9. The calibrator was tested with CE-marked methods and shown to be negative for antibodies to HIV and HCV and nonreactive for HBsAg. However, as no testing method can rule out the potential risk of infection with absolute certainty, so this material should be handled as a patient specimen.
10. Reagent pack opened for more than 28 days is not recommended for use.
11. Wear clean gloves to avoid contaminating the reagents.
12. Instability or deterioration should be suspected if there are visible signs of leakage, turbidity, precipitates or microbial growth.
13. Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
14. All the specimen and reaction wastes should be considered potentially biohazard. Specimens and reaction wastes should be handled in accordance with the local regulations and guidelines.
15. The Material Safety Data Sheet (MSDS) is available upon request.
16. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the local competent authority.
17. This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

Warning

H317 May cause an allergic skin reaction.

H412 Harmful to aquatic life with long lasting effects.

Prevention:

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P272 Contaminated work clothing should not be allowed out of the workplace.

P273 Avoid release to the environment.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

Response:

P302 + P352 IF ON SKIN: Wash with plenty of water.

P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.

P362 + P364 Take off contaminated clothing and wash it before reuse.

Disposal:

P501 Dispose of contents/ container in accordance with local regulation.

Graphical Symbols

In vitro
diagnostic
medical
device

Authorized
Representative
in the European
Community

Consult
instructions
for use

European
Conformity

Catalogue
number

Temperature
limit

Manufacturer

Use-by
date

Biological
risks

Batch code

This way up

Unique
device
identifier

Caution

Bibliography

1. National Health Commission of the People's Republic of China. WS/T486-2015 Diagnosis of *Toxoplasma gondii*.
2. Jones J, Lopez A, Wilson M. Congenital toxoplasmosis [J]. *Am Fam Physician*, 2003, 67 (10): 2131-2138.
3. Ortiz-Alegria LB, Caballero-Ortega H, Canedo-Solares I, et al. Congenital toxoplasmosis: candidate host immune genes relevant for vertical transmission and pathogenesis [J]. *Genes Immun*, 2010, 11 (5): 363-373.
4. Montoya JG, Liesenfeld O. Toxoplasmosis [J]. *Lancet*, 2004, 363 (9425): 1965-1976.
5. Montoya JG, Remington JS. Management of *Toxoplasmosis gondii* infection during pregnancy [J]. *Clin Infect Dis*, 2008, 47 (4): 554-566.
6. Montoya JG, Rosso F. Diagnosis and management of toxoplasmosis [J]. *ClinPerinatol*, 2005, 32 (3):705-726.
7. Robertgangneux F. Epidemiology of and diagnostic strategies for toxoplasmosis [J]. *Clinical Microbiology Reviews*, 2012, 25(2):264-296.
8. Maldonado Y A , Read J S . Diagnosis, Treatment, and Prevention of Congenital Toxoplasmosis in the United States[J]. *Pediatrics*, 2017, 139(2):e20163860.
9. National Committee for Clinical Laboratory Standards. Clinical use and interpretation of serologic tests for *toxoplasma gondii*; Approved guideline, M36-A, 2004.
10. Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem* 1988;34:261-264.
11. Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45:879-885
12. Boscatto, LM and Stuart, MC. Heterophilic antibodies; A problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
13. CLSI. EP17-A2: Vol. 32, No. 8, Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition.

14. CLSI. EP06-A2: Vol. 40, No. 16, Evaluation of Linearity of Quantitative Measurement Procedures; Approved Guideline-Second Edition.
15. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline-Third Edition.

© 2022-2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
All rights Reserved.

Manufacturer: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

E-mail Address: service@mindray.com

Website: www.mindray.com

Tel: +86-755-81888998

Fax: +86-755-26582680

EC-Representative: Shanghai International Holding Corp. GmbH
(Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

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1. National Health Commission of the People's Republic of China. WS/T486-2015 Diagnosis of *Toxoplasma gondii*.
2. Jones J, Lopez A, Wilson M. Congenital toxoplasmosis [J]. Am Fam Physician, 2003, 67 (10): 2131-2138.
3. Ortiz-Alegria LB, Caballero-Ortega H, Canedo-Solares I, et al. Congenital toxoplasmosis: candidate host immune genes relevant for vertical transmission and pathogenesis [J]. Genes Immun, 2010, 11 (5): 363-373.
4. Montoya JG, Liesenfeld O. Toxoplasmosis [J]. Lancet, 2004, 363 (9425): 1965-1976.
5. Montoya JG, Remington JS. Management of *Toxoplasmosis gondii* infection during pregnancy [J]. Clin Infect Dis, 2008, 47 (4): 554-566.
6. Montoya JG, Rosso F. Diagnosis and management of toxoplasmosis [J]. ClinPerinatol, 2005, 32 (3):705-726.
7. Robertgangneux F. Epidemiology of and diagnostic strategies for toxoplasmosis [J]. Clinical Microbiology Reviews, 2012, 25(2):264-296.
8. Maldonado Y A , Read J S . Diagnosis, Treatment, and Prevention of Congenital Toxoplasmosis in the United States[J]. Pediatrics, 2017, 139(2):e20163860.
9. National Committee for Clinical Laboratory Standards. Clinical use and interpretation of serologic tests for *toxoplasma gondii*; Approved guideline, M36-A, 2004.
10. Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. Clin Chem 1988;34:261-264.
11. Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobuline responses in patients receiving monoclonal antibody therapy. Cancer Res 1985;45:879-885
12. Boscato, LM and Stuart, MC. Heterophilic antibodies; A problem for all immunoassays. Clin Chem 1988;34(1):27-33.
13. CLSI. EP17-A2: Vol. 32, No. 8, Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition.

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Адрес: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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7. Robertgangneux F. Epidemiology of and diagnostic strategies for toxoplasmosis [J]. Clinical Microbiology Reviews, 2012, 25(2):264-296.
8. Maldonado Y A, Read J S. Diagnosis, Treatment, and Prevention of Congenital Toxoplasmosis in the United States[J]. Pediatrics, 2017, 139(2):e20163860.
9. National Committee for Clinical Laboratory Standards. Clinical use and interpretation of serologic tests for toxoplasma gondii; Approved guideline, M36-A, 2004.
10. Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. Clin Chem 1988;34:261-264.
11. Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobuline responses in patients receiving monoclonal antibody therapy. Cancer Res 1985;45:879-885
12. Boscato, LM and Stuart, MC. Heterophilic antibodies; A problem for all immunoassays. Clin Chem 1988;34(1):27-33.
13. CLSI. EP17-A2: Vol. 32, No. 8, Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition.
14. CLSI. EP06-A2: Vol. 40, No. 16, Evaluation of Linearity of Quantitative Measurement Procedures; Approved Guideline-Second Edition.
15. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline-Third Edition.

© 2022-2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

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Fabricante: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Endereço: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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Precaución

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1. National Health Commission of the People's Republic of China. WS/T486-2015 Diagnosis of Toxoplasma gondii.
2. Jones J, Lopez A, Wilson M. Congenital toxoplasmosis [J]. Am Fam Physician, 2003, 67 (10): 2131-2138.
3. Ortiz-Alegria LB, Caballero-Ortega H, Canedo-Solares I, et al. Congenital toxoplasmosis: candidate host immune genes relevant for vertical transmission and pathogenesis [J]. Genes Immun, 2010, 11 (5): 363-373.
4. Montoya JG, Liesenfeld O. Toxoplasmosis [J]. Lancet, 2004, 363 (9425): 1965-1976.
5. Montoya JG, Remington JS. Management of Toxoplasmosis gondii infection during pregnancy [J]. Clin Infect Dis, 2008, 47 (4): 554-566.
6. Montoya JG, Rosso F. Diagnosis and management of toxoplasmosis [J]. ClinPerinatol, 2005, 32 (3):705-726.
7. Robertgangneux F. Epidemiology of and diagnostic strategies for toxoplasmosis [J]. Clinical Microbiology Reviews, 2012, 25(2):264-296.
8. Maldonado Y A, Read J S. Diagnosis, Treatment, and Prevention of Congenital Toxoplasmosis in the United States[J]. Pediatrics, 2017, 139(2):e20163860.

9. National Committee for Clinical Laboratory Standards. Clinical use and interpretation of serologic tests for toxoplasma gondii; Approved guideline, M36-A, 2004.
10. Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. Clin Chem 1988;34:261-264.
11. Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobuline responses in patients receiving monoclonal antibody therapy. Cancer Res 1985;45:879-885
12. Boscato, LM and Stuart, MC. Heterophilic antibodies; A problem for all immunoassays. Clin Chem 1988;34(1):27-33.
13. CLSI. EP17-A2: Vol. 32, No. 8, Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition.
14. CLSI. EP06-A2: Vol. 40, No. 16, Evaluation of Linearity of Quantitative Measurement Procedures; Approved Guideline-Second Edition.
15. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline-Third Edition.

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6. Montoya JG, Rosso F. Diagnosis and management of toxoplasmosis [J]. ClinPerinatol, 2005, 32 (3):705-726.
7. Robertgangneux F. Epidemiology of and diagnostic strategies for toxoplasmosis [J]. Clinical Microbiology Reviews, 2012, 25(2):264-296.
8. Maldonado Y A , Read J S . Diagnosis, Treatment, and Prevention of Congenital Toxoplasmosis in the United States[J]. Pediatrics, 2017, 139(2):e20163860.
9. National Committee for Clinical Laboratory Standards. Clinical use and interpretation of serologic tests for toxoplasma gondii; Approved guideline,M36-A,2004.
10. Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. Clin Chem 1988;34:261-264.
11. Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobuline responses in patients receiving monoclonal antibody therapy. Cancer Res 1985;45:879-885.
12. Boscatto, LM and Stuart, MC. Heterophilic antibodies; A problem for all immunoassays. Clin Chem 1988;34(1):27-33.
13. CLSI. EP17-A2: Vol. 32, No. 8, Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition.
14. CLSI. EP06-A2: Vol. 40, No. 16, Evaluation of Linearity of Quantitative Measurement Procedures; Approved Guideline-Second Edition.
15. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline-Third Edition.

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Indirizzo: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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Seri kodu

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Dikkat

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1. National Health Commission of the People's Republic of China.WS/T486-2015 Diagnosis of Toxoplasma gondii.
2. Jones J, Lopez A, Wilson M. Congenital toxoplasmosis [J]. Am Fam Physician, 2003, 67 (10): 2131-2138.
3. Ortiz-Alegria LB, Caballero-Ortega H, Canedo-Solares I, et al. Congenital toxoplasmosis: candidate host immune genes relevant for vertical transmission and pathogenesis [J]. Genes Immun, 2010, 11 (5): 363-373.
4. Montoya JG, Liesenfeld O. Toxoplasmosis [J]. Lancet, 2004, 363 (9425): 1965-1976.
5. Montoya JG, Remington JS. Management of Toxoplasmosis gondii infection during pregnancy [J]. Clin Infect Dis, 2008, 47 (4): 554-566.
6. Montoya JG, Rosso F. Diagnosis and management of toxoplasmosis [J]. ClinPerinatol, 2005, 32 (3):705-726.
7. Robertgangneux F. Epidemiology of and diagnostic strategies for toxoplasmosis [J]. Clinical Microbiology Reviews, 2012, 25(2):264-296.
8. Maldonado Y A , Read J S . Diagnosis, Treatment, and Prevention of Congenital Toxoplasmosis in the United States[J]. Pediatrics, 2017, 139(2):e20163860.
9. National Committee for Clinical Laboratory Standards. Clinical use and interpretation of serologic tests for toxoplasma gondii; Approved guideline,M36-A,2004.
10. Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. Clin Chem 1988;34:261-264.
11. Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobuline responses in patients receiving monoclonal antibody therapy. Cancer Res 1985;45:879-885
12. Boscato, LM and Stuart, MC. Heterophilic antibodies; A problem for all immunoassays. Clin Chem 1988;34(1):27-33.
13. CLSI. EP17-A2: Vol. 32, No. 8, Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition.

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Üretici: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Adres: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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1. National Health Commission of the People's Republic of China. WS/T486-2015 Diagnosis of *Toxoplasma gondii*.
2. Jones J, Lopez A, Wilson M. Congenital toxoplasmosis [J]. Am Fam Physician, 2003, 67 (10): 2131-8.
3. Ortiz-Alegria LB, Caballero-Ortega H, Canedo-Solares I, et al. Congenital toxoplasmosis: candidate host immune genes relevant for vertical transmission and pathogenesis [J]. Genes Immun, 2010, 11 (5): 363-73.
4. Montoya JG, Liesenfeld O. Toxoplasmosis [J]. Lancet, 2004, 363 (9425): 1965-76.
5. Montoya JG, Remington JS. Management of *Toxoplasmosis gondii* infection during pregnancy [J]. Clin Infect Dis, 2008, 47 (4): 554-66.
6. Montoya JG, Rosso F. Diagnosis and management of toxoplasmosis [J]. ClinPerinatol, 2005, 32 (3):705-26.
7. Robertgangneux F. Epidemiology of and diagnostic strategies for toxoplasmosis [J]. Clinical Microbiology Reviews, 2012, 25(2):264-96.
8. Maldonado Y A , Read J S . Diagnosis, Treatment, and Prevention of Congenital Toxoplasmosis in the United States[J]. Pediatrics, 2017, 139(2):e20163860.
9. National Committee for Clinical Laboratory Standards. Clinical use and interpretation of serologic tests for *toxoplasma gondii*; Approved guideline, M36-A, 2004.
10. Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving

murine monoclonal antibodies for diagnosis and therapy. Clin Chem 1988;34:261-264.

11. Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobuline responses in patients receiving monoclonal antibody therapy. Cancer Res 1985;45:879-885

12. Boscatto, LM and Stuart, MC. Heterophilic antibodies; A problem for all immunoassays. Clin Chem 1988;34(1):27-33.

13. CLSI. EP17-A2: Vol. 32, No. 8, Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition.

14. CLSI. EP06-A2: Vol. 40, No. 16, Evaluation of Linearity of Quantitative Measurement Procedures; Approved Guideline-Second Edition.

15. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline – Third Edition.

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