IL-6

Interleukin-6 (CLIA)

Order Information

Catalog No.	Package Size
105-027046-00	2×50 tests
105-027047-00	2×100 tests

Intended Use

The CL-series Interleukin-6 (IL-6) assay is a Chemiluminescent Immunoassay (CLIA) for the quantitative detection of Interleukin-6 (IL-6) in human serum and plasma. This assay can be used to monitor the immune status and inflammatory reaction.

Summary

Interleukin-6 (IL-6) is a pleiotropic cytokine with a large range of biological activities. It was described as interferon-β2, hepatocyte stimulating factor, and later it was described as human B cell stimulating factor 2 (BSF2). In 1988, it was named IL-6^[1]. Further studies demonstrated that IL-6 shows activities on B cells, T cells, hematopoietic stem cells, hepatocytes and brain cells ^[2]. IL-6 contains 212 amino acids, and the N-terminal is easy to be cut into 184 amino acid polypeptide 2 with a molecular weight of 22-27 kDa.

IL-6 is produced by lymphocytes and non-lymphoid cells, and by normal and transformed cells. It regulates the immune response, acute phase response and hematopoiesis ^[3]. IL-6 production is rapidly induced in the process of acute inflammatory response caused by injury, trauma, stress, infection, brain death, tumor formation and other conditions ^[2]. IL-6 is an alarm marker in the early stage of sepsis ^[4-7] and its concentration also aids in the assessment of severity in inflammatory diseases ^[8-10]. IL-6 also plays a role in chronic inflammatory response of diseases, such as rheumatoid arthritis ^[11,12].

Assay Principles

The CL-series IL-6 assay is a sandwich immunoassay to quantitatively detect Interleukin-6.

In the first step, sample, sample diluent and paramagnetic microparticles coated with IL-6 antibodies and anti-IL-6 antibody-alkaline phosphatase (ALP) conjugates are added into a reaction cuvette. After incubation, IL-6 in the sample will bind to IL-6 antibodies coated microparticles and anti-IL-6 antibody-ALP conjugates. Afterwards, microparticles are magnetically captured while other unbound substances are removed by washing.

In the second step, chemiluminescent substrate (sodium 3-[(1R,3S,5,7r)- 4'-methoxy Spiro[adamantane-2,3'-[1,2] dioxetan]-4'-yl)]phenyl phosphate, AMPPD) is added to the reaction cuvette. AMPPD was removed a phosphate group by ALP to form an unstable intermediate structure, which produces methyl methoxybenzoate anion through intramolecular electron transfer. Chemiluminescence is generated when the methyl methoxybenzoate anion in the excited state returns to the ground state. The resulting chemiluminescent reaction is measured as relative light units (RLUs) by a photomultiplier in the instrument.

The amount of IL-6 present in the sample is proportional to the relative light units (RLUs) generated during the reaction. The IL-6 results can be determined via a calibration curve, which is established on an encoded Calibration Master Curve and three-level product calibrators.

Reagent Components

Components in different lots of reagent kit cannot be interchanged:

Ra	Paramagnetic microparticles coated with anti-IL-6 antibody (mouse) in 50mM TRIS [®] buffer. Minimum concentration: 0.5 g/L. Preservative: 0.02% BND and 0.048% Proclin 300.
Rb	Anti-IL-6 antibody (mouse)-alkaline phosphatase conjugate in 50mM MES [®] buffer. Minimum concentration: 0.5 mg/L. Preservative: 0.02% BND and 0.048% Proclin 300.
Rc	Sample diluents in 50mM TRIS buffer. Preservative: 0.09% sodium azide and 0.048% Proclin 300.

- a) TRIS=Tris (hydroxymethyl)-aminomethane
- b) MES=2-morpholino-ethane sulfonic acid

Storage and Stability

The unopened IL-6 (CLIA) reagent is stable up to the stated expiration date indicated on the label when stored at 2-8°C.

The IL-6 (CLIA) reagent can be stored onboard at 2-8°C and used for a maximum of 28 days after opening for use.

Reagent Preparation

The reagents in the kit have been assembled into an intact unit that cannot be separated.

Materials Required but not Provided

Mindray CL-series Chemiluminescence Immunoassay Analyzer
Mindray IL-6 Calibrators
Mindray IL-6 Control (L)
Mindray IL-6 Control (H)
Mindray Substrate Solution
Mindray Wash Buffer
Mindray Reaction Cuvette
Mindray Sample Diluent

Applicable System

Mindray CL-series Chemiluminescence Immunoassay Analyzer

Specimen Collection and Preparation

Specimen Types

- Human serum and plasma (K₂-EDTA, lithium heparin) is recommended for this assay.
- Blood collection tubes from different manufacturers may contain different materials and additives which could affect the test results in some cases. Not all available tubes of all manufacturers were tested by Mindray. Laboratories should determine the suitability of the blood collection tubes and serum separation products. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Specimen Conditions

Do not use specimen with the following conditions:

- heat inactivated
- grossly hemolyzed
- obvious microbial contamination
- visible fibrin or other cell debris

For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.

Preparation for Analysis

- Follow blood collection tube manufacturer's recommendations

for centrifugation. Centrifuge the specimens after clot formation is complete. Ensure residual fibrin and cellular matter has been removed prior to analysis.

- For optimal results, inspect all samples for bubbles. Remove bubbles with a pipette tip prior to analysis. Specimens must be mixed thoroughly after thawing. Thawed samples should be centrifuged prior to use.
- If the sample was covered with lipid layer after centrifugation, the sample should be transferred to a clean tube and centrifuged before testing. Do not transfer the lipid layer. Handle carefully to prevent cross contamination.

Specimen Storage

- Specimens should be tested as soon as possible after sample collection. If testing is not completed within 5 hours, specimens should be tightly capped and refrigerated at 2-8°C. The samples are stable for 1 day at 2-8°C and 90 days at -20°C or below.
- Avoid repeated freeze and thaw cycles. The samples may be frozen only once.
- Frozen or refrigerated samples, controls and calibrators must be restored to room temperature before analysis.

Assay Procedure

For optimal performance of this assay, operators should read the related operator's manual carefully to get enough information such as operation instructions, sample preservation and management, safety precautions, and maintenance. Prepare all required materials for the assay as well.

Before loading the IL-6 reagent kit on the instrument for the first time, invert unopened reagent bottle gently for at least 30 times to resuspend the microparticles, which have settled during shipment or storage. Avoid foam formation. Visually inspect the bottle to ensure the microparticles have been well mixed. If the microparticles remain adhered to the bottle, continue inverting until suspended completely. If the microparticles cannot be homogenized, it is recommended not to use this bottle of reagent. Contact Mindray Customer Service for help. Avoid inverting opened reagent bottle.

This assay requires 30 µL of sample volume for a single test not including dead volume of the sample container. Additional volume is required when performing additional tests from the same sample. Operators should refer to the operator's manual and follow the specific requirements of the assay to determine the minimum sample volume.

Calibration

CL-series IL-6 (CLIA) has been standardized against the NIBSC (National Institute for Biological Standards and Control) 1st IS 89/548 Standard.

The specific information of calibration master curve of IL-6 (CLIA) reagent kit is stored in the two-dimensional barcode attached in the reagent pack. It's used in combination with IL-6 calibrators for the calibration of the specific reagent lot. When performing the calibration, scan the information of calibration master curve from the barcode into the system first, and then use the IL-6 calibrators at three levels.

Recalibration is required as follows:

- The previous calibration was performed 28 days before.
- Quality controls lie outside the expected ranges.
- A reagent kit with a new lot number is used.

For detailed instruction of calibration, refer to the operator's manual.

Quality Control

It is recommended that quality controls should be run once every

24 hours if the tests are in use and following each calibration. The quality control frequency should be adapted to each laboratory's individual requirements.

Quality control results should be within the acceptable ranges. If control results fall out of the defined range, the associated test results are invalid and the samples must be retested. Recalibration may be required. Examine the assay system by referring to the operator's manual. If the quality control results are still out of the defined range, please contact Mindray Customer Service for help.

Calculation

The analyzer automatically calculates the analyte concentration of each sample on the calibration master curve read from the barcode, and a 4-Parameter Logistic Curve Fitting (4PLC) with the relative light units (RLUs) generated from three-level calibrators of defined concentration values. The results are shown in the unit of pg/mL.

Specimen Dilution

Samples with IL-6 concentrations above the upper limit (5,000pg/mL) can be diluted with Mindray Sample Diluent. The recommended dilution is 1:40 (manually) or 1:10 automatically by the analyzer. The concentration of the diluted sample must be >100 pg/mL. After manual dilution, multiply the result by the dilution factor. After automated dilution by the analyzers, the system automatically multiplies the result by the dilution factor when calculating the sample concentration.

Expected Values

A study on a cohort of 1022 apparently healthy individuals of different age and gender has determined the 95th Percentile is 7pg/mL.

Interpretation of Results

The assay results used in clinical decisions should be made in conjunction with other evidences, such as symptoms, clinical history, and other clinical information.

Limitation of Measurement

The test results are for clinical reference only and cannot be used alone as the basis for diagnosis or ruling out cases.

The measuring range of reagent kit is 2.5-5,000 pg/mL. A specimen with IL-6 concentration lower than the lower limit will be reported as <2.5 pg/mL. A specimen with IL-6 concentration higher than the upper limit will be reported as >5,000 pg/mL or repeat the test after dilution.

The IL-6 result in a given specimen, determined with assays from different manufacturers, can vary due to differences in assay methods, calibration, and reagent specificity.

Performance Characteristics

Lower limits of measurement

Limit of Blank, Limit of Detection and Limit of Quantitation
Limit of Blank = 1.0 pg/mL
Limit of Detection =1.5 pg/mL
Limit of Quantitation =2.5 pg/mL.

Precision

Precision was determined by following CLSI Protocol EP15-A3 ^[13]. Two levels of quality controls and two levels of samples were tested 2 runs each day in duplicate for totally 20 days (n=80). Results are shown in the table below*.

Sample	Mean (pg/mL)	Within-run Precision		Total Precision	
		SD (pg/mL)	CV	SD(pg/mL)	CV
QC-L	42.18	0.75	1.78%	0.76	1.80%
QC-H	259.00	5.14	1.98%	5.19	2.00%
Sample1	49.28	0.85	1.72%	0.85	1.73%
Sample 2	439.99	7.30	1.66%	7.35	1.67%

*Representative data; results in individual laboratories may vary.

Linearity

Linearity study was performed based on guidance from CLSI EP06-A^[4]. A high concentration IL-6 sample (>5,000 pg/mL) was mixed with a low concentration sample (<1.5 pg/mL) at different ratios, generating a series of dilutions. The Linearity of IL-6 was demonstrated in the range of 1.5 to 5000.00 pg/mL, the correlation coefficient r ≥0.9900. The linearity data are summarized in the table below*.

No.	Expected value (pg/mL)	Measured value (pg/mL)
A high concentration sample	5630.02	5630.02
1:5	1126.00	1091.99
1:10	563.00	559.72
1:20	281.50	286.71
1:40	140.75	145.53
1:80	70.38	75.58
1:160	35.19	39.16
1:320	17.59	20.41
A low concentration sample	0.00	0.05

*Representative data; results in individual laboratories may vary.

Accuracy

Two samples defined to the reference reagent were tested to verify the accuracy of IL-6 assay. Results are shown in the table below (criterion: relative deviation within ± 15% of defined value)*.

Sample	Defined Value (pg/mL)	Measured Value (pg/mL)	Relative Deviation
Accuracy control 1	44.07	42.68	3.26%
Accuracy control 2	265.92	264.67	0.47%

*Representative data; results in individual laboratories may vary.

High Dose Hook

For the Mindray CL-Series IL-6 assay, no high dose hook effect was observed when samples contain up to approximately 200,000 pg/mL of IL-6.

Analytical Specificity

Samples with hyperlipidemia, icterus, hemolysis, high total protein, high rheumatoid factor, or positive antinuclear antibody (ANA) may cause inaccurate test results. Heterophilic antibodies in samples can react with immunoglobulin in kit and interfere with the immunoassay results^[15] Patients who have been regularly exposed to animals or have been treated with immunoglobulins or immunoglobulin fragments may produce such antibodies, resulting in abnormal results. For example, specimens from individuals who have been exposed to mouse monoclonal antibodies may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values with assay kits employing mouse monoclonal antibodies. ^[16-17] There are anti-interference components in Mindray IL-6 (CLIA) reagent kit, which can effectively reduce the influence caused by HAMA in samples, but there may still be interference problems in a few samples.

The test results of Mindray IL-6 are not observed to be interfered with the potential endogenous substances up to the following concentrations (criterion: relative deviation ± 10% of initial value)*.

Substance	Concentration
Hemoglobin	1000 mg/dL
Bilirubin	40 mg/dL
Triglyceride	1500mg/dL
Biotin	1200 ng/mL
Alkaline Phosphatase (ALP)	1500 U/L

*Representative data; results in individual laboratories may vary.

The test results of Mindray IL-6 are not observed to be interfered with the potential interfering drugs up to the following concentrations (criterion: relative deviation ± 10% of initial value)*.

Substance	Non-interfering concentrations
Imipenem	1.18 mg/mL
Cefotaxime	0.9 mg/mL
Vancomycin	3.5 mg/mL
Dopamine	0.13 mg/mL
Noradrenaline	0.002 mg/mL
Furosemide	0.02 mg/mL
Paracetamol	20 mg/dL
Acetylsalicylic acid	50 mg/dL
Ibuprofen	40 mg/dL

*Representative data; results in individual laboratories may vary.

The test results of Mindray IL-6 does not show any significant cross-reactivity with the following substances (criterion: relative deviation within ±10% of initial value)*.

Substance	Concentration
Interleukin-1α	50 ng/mL
Interleukin-1β	50 ng/mL
Interleukin-2	50 ng/mL
Interleukin-3	50 ng/mL
Interleukin-4	50 ng/mL
Interleukin-8	50 ng/mL
Interferon-γ	50 ng/mL
TNF-α	50 ng/mL
G-CSF	50 ng/mL
Soluble IL-6 receptor (IL-6sR)	50 ng/mL
Soluble glycoprotein 130 (sgp130)	50 ng/mL
IL-6sR/sgp130 complex	50 ng/mL

*Representative data; results in individual laboratories may vary.

Method Comparison

The Mindray CL-Series IL-6 Assay was compared to a commercially available diagnostic kit in a correlation study with 399 specimens. The statistical data obtained by linear regression are shown in the table below.*

Concentration Range (pg/mL)	Slope	Intercept	Correlation Coefficient
1.5-5000	0.9989	-0.3699	0.9999

*Representative data; results in individual laboratories may vary.

Warnings and Precautions

- For in vitro diagnostic use only. For laboratory professional use.
- Follow all the rules in handling laboratory reagents and take necessary safety precautions.
- Due to the differences in methodology and antibody specificity, test results of the same sample may be different when using reagent kits from different manufacturers on Mindray system, or using Mindray reagent kits on other systems.
- Do not use reagent kits beyond the expiration date.
- Do not use reagents mixed from different reagent lots.
- Avoid bubble formation in all reagents. Wear clean gloves to avoid contaminating the reagents.
- Always keep the reagent pack in the upright position to ensure no microparticle has been lost prior to use.

- Please confirm the integrity of the package before use. Do not use the reagents with damaged packages. Inappropriate storage will lead to incorrect results.
- If the reagents are opened unintentionally before use, they shall be used as soon as possible.
- Reagent pack opened for more than 28 days is not recommended for use.
- Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
- All the specimen and reaction wastes should be considered potentially biohazard. The handling of specimens and reaction wastes should be in accordance with the local regulations and guidelines.
- Instability or deterioration should be suspected if there are visible signs of leakage, turbidity, precipitates or microbial growth.
- Do not freeze. The results can't be assured when the reagents are stored at inappropriate condition.
- Swallowing or getting the reagents on yourself must be avoided. Wash hands thoroughly after handling. The device must be used in a room with adequate ventilation. The contaminated clothing must be removed immediately.
- The Material Safety Data Sheet (MSDS) is available upon request.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the local competent authority.

Graphical Symbols

In vitro
Diagnostic
medical
device

Authorized
representative
in the European
Community

Consult
instructions
for use

European
Conformity

Catalogue
number

Temperature
limit

Manufacturer

Use-by date

Caution

Batch code

This way up

Unique
device
identifier

References

- Yamasaki, K., Taga, T., Hirata, Y. et al. Cloning and expression of the human interleukin-6 (BSF-2/IFN beta 2) receptor. Science 1988; 241:825.
- Song M, Kellum JA. Interleukin-6. Crit Care Med 2005;33(Suppl12):463-465.
- Hirano, T. Interleukin-6 and its relation to inflammation and disease. Clin Immunol And Immunopathol 1992; 62:S60
- Oda S, Hirasawa H, Shiga H, et al. Sequential measurement of IL-6 blood levels in patients with systemic inflammatory response syndrome (SIRS)/sepsis. Cytokine 2005;29:169-175.
- Jawa RS, Anillo S, Hutton K, et al. Interleukin-6 in Surgery, Trauma, and Critical Care Part II: Clinical Implications. J Intensive Care Med 2011;26(2):73-87.
- Bloos F, Reinhart K. Rapid diagnosis of sepsis. Virulence. 2014 Jan 1;5(1):154-60.
- Molano FD, Arevalo-Rodriguez I, Roqué I FM, et al. Plasma interleukin-6 concentration for the diagnosis of sepsis in critically ill adults. Cochrane Database Syst Rev. 2019 Apr

30;4(4):CD011811.

- Marti AL, Moreno CA. Prognostic value of interleukin-6 for death in patients with sepsis. Med Clin (Barc). 2017 Dec 7;149(11):507.
- Miguel-Bayarri V, Casanoves-Laparra EB, Pallás-Beneyto L, et al. Prognostic value of the biomarkers procalcitonin, interleukin-6 and C-reactive protein in severe sepsis. 2012 Nov;36(8):556-62.
- Jekarl DW, Lee SY, Lee J, et al. Procalcitonin as a diagnostic marker and IL-6 as a prognostic marker for sepsis. Diagn Microbiol Infect Dis 2013;75(4):342-347.
- Madhok, R. et al. Serum interleukin 6 levels in rheumatoid arthritis: correlations with clinical and laboratory indices of disease activity. Ann Rheum Dis 1993; 52: 232.
- Yao X, Huang J, Zhong H, et al. Targeting interleukin-6 in inflammatory autoimmune diseases and cancers. Pharmacol Ther. 2014 Feb;141(2):125-39.
- CLSI. EP15-A3: Vol. 34 No. 12, User Verification of Precision and Estimation of Bias; Approved Guideline—Third Edition.
- CLSI. EP06-A: Vol. 23, No. 16, Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline.
- Boscato Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
- Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037–1038.
- Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613–621.

© 2022-2026 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
All rights Reserved.
Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Address: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China
E-mail Address: service@mindray.com
Website: www.mindray.com
Tel: +86-755-81889998
Fax: +86-755-26582680
EC-Representative: Shanghai International Holding Corp. GmbH (Europe)
Address: Eiffestraße 80, 20537 Hamburg, Germany
Tel: 0049-40-2513175
Fax: 0049 -40 -255726

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- Yamasaki, K., Taga, T., Hirata, Y. et al. Cloning and expression of the human interleukin-6 (BSF-2/IFN beta 2) receptor. Science 1988; 241:825.
- Song M, Kellum JA. Interleukin-6. Crit Care Med 2005;33(Suppl12):463-465.
- Hirano, T. Interleukin-6 and its relation to inflammation and disease. Clin Immunol And Immunopathol 1992; 62:560
- Oda S, Hirasawa H, Shiga H, et al. Sequential measurement of IL-6 blood levels in patients with systemic inflammatory response syndrome (SIRS)/sepsis. Cytokine 2005;29:169-175.
- Jawa RS, Anillo S, Huntoon K, et al. Interleukin-6 in Surgery, Trauma, and Critical Care Part II: Clinical Implications. J Intensive Care Med 2011;26(2):73-87.
- Bloos F, Reinhart K. Rapid diagnosis of sepsis. Virulence. 2014 Jan 15;5(1):154-60.
- Molano FD, Arevalo-Rodriguez I, Roqué I FM, et al. Plasma interleukin-6 concentration for the diagnosis of sepsis in critically ill adults. Cochrane Database Syst Rev. 2019 Apr 30;4(4):CD0011811.
- Marti AL, Moreno CA. Prognostic value of interleukin-6 for death in patients with sepsis. Med Clin (Barc). 2017 Dec 7;149(11):507.
- Miguel-Bayarri V, Casanoves-Laparra EB, Pallás-Beneyto L, et al. Prognostic value of the biomarkers procalcitonin, interleukin-6 and C-reactive protein in severe sepsis. 2012 Nov;36(8):556-62.
- Jekarl DW, Lee SY, Lee J, et al. Procalcitonin as a diagnostic marker and IL-6 as a prognostic marker for sepsis. Diagn Microbiol Infect Dis 2013;75(4):342-347.
- Madhok R, et al. Serum interleukin 6 levels in rheumatoid arthritis: correlations with clinical and laboratory indices of disease activity. Ann Rheum Dis 1993; 52: 232.
- Yao X, Huang J, Zhong H, et al. Targeting interleukin-6 in inflammatory autoimmune diseases and cancers. Pharmacol Ther. 2014 Feb;141(2):125-39.
- CLSI. EP15-A3: Vol. 34 No. 12, User Verification of Precision and Estimation of Bias; Approved Guideline—Third Edition.
- CLSI. EP06-A: Vol. 23, No. 16, Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline.
- Boscato Ln, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
- Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037–1038.
- Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613–621.

© 2022-2026 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Все права защищены.
Производитель: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Адрес: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China
Адрес электронной почты: service@mindray.com
Веб-сайт: www.mindray.com
Тел.: +86-755-81888998
Факс: +86-755-26582680
Представитель в ЕС: Shanghai International Holding Corp. GmbH (Europe)
Адрес: Eiffestraße 80, 20537 Hamburg, Germany
Тел.: 0049-40-2513175
Факс: 0049-40-2557226

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Identificador exclusivo do dispositivo

Referências

- Yamasaki, K., Taga, T., Hirata, Y. et al. Cloning and expression of the human interleukin-6 (BSF-2/IFN beta 2) receptor. Science 1988; 241:825.
- Song M, Kellum JA. Interleukin-6. Crit Care Med 2005;33(Suppl12):463-465.
- Hirano, T. Interleukin-6 and its relation to inflammation and disease. Clin Immunol And Immunopathol 1992; 62:S60
- Oda S, Hirasawa H, Shiga H, et al. Sequential measurement of IL-6 blood levels in patients with systemic inflammatory response syndrome (SIRS)/sepsis. Cytokine 2005;29:169-175.
- Jawa RS, Anillo S, Huntoon K, et al. Interleukin-6 in Surgery, Trauma, and Critical Care Part II: Clinical Implications. J Intensive Care Med 2011;26(2):73-87.
- Bloos F, Reinhart K. Rapid diagnosis of sepsis. Virulence. 2014 Jan 1;5(1):154-60.
- Molano FD, Arevalo-Rodriguez I, Roqué I FM, et al. Plasma interleukin-6 concentration for the diagnosis of sepsis in critically ill adults. Cochrane Database Syst Rev. 2019 Apr 30;4(4):CD011811.
- Martí AL, Moreno CA. Prognostic value of interleukin-6 for death in patients with sepsis. Med Clin (Barc). 2017 Dec 7;149(11):507.
- Miguel-Bayarri V, Casanoves-Laparra EB, Pallás-Beneyto L, et al. Prognostic value of the biomarkers procalcitonin, interleukin-6 and C-reactive protein in severe sepsis. 2012 Nov;36(8):556-62.

- Jekarl DW, Lee SY, Lee J, et al. Procalcitonin as a diagnostic marker and IL-6 as a prognostic marker for sepsis. Diagn Microbiol Infect Dis 2013;75(4):342-347.
- Madhok, R. A. et al. Serum interleukin 6 levels in rheumatoid arthritis: correlations with clinical and laboratory indices of disease activity. Ann Rheum Dis 1993; 52: 232.
- Yao X, Huang J, Zhong H, et al. Targeting interleukin-6 in inflammatory autoimmune diseases and cancers. Pharmacol Ther. 2014 Feb;141(2):125-39.
- CLSI. EP15-A3: Vol. 30, 12, User Verification of Precision and Estimation of Bias; Approved Guideline—Third Edition.
- CLSI. EP06-A: Vol. 23, No. 16, Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline.
- Boscatto Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
- Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
- Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.

© 2022-2026 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Todos os direitos reservados.
Fabricante: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Endereço: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China
Endereço de e-mail: service@mindray.com
Site: www.mindray.com
Tel: +86-755-81888998
Fax: +86-755-26582680
Representante da EC: Shanghai International Holding Corp. GmbH (Europe)
Endereço: Eiffestraße 80, 20537 Hamburg, Germany
Tel: 0049-40-2513175
Fax: 0049-40-255726

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Referencias

- Yamasaki, K., Taga, T., Hirata, Y. et al. Cloning and expression of the human interleukin-6 (BSF-2/IFN beta 2) receptor. Science 1988; 241:825.
- Song M, Kellum JA. Interleukin-6. Crit Care Med 2005;33(Suppl12):463-465.
- Hirano, T. Interleukin-6 and its relation to inflammation and disease. Clin Immunol And Immunopathol 1992; 62:S60
- Oda S, Hirasawa H, Shiga H, et al. Sequential measurement of IL-6 blood levels in patients with systemic inflammatory response syndrome (SIRS)/sepsis. Cytokine 2005;29:169-175.
- Jawa RS, Anillo S, Huntoon K, et al. Interleukin-6 in Surgery, Trauma, and Critical Care Part II: Clinical Implications. J Intensive Care Med 2011;26(2):73-87.
- Bloos F, Reinhart K. Rapid diagnosis of sepsis. Virulence. 2014 Jan 1;5(1):154-60.
- Molano FD, Arevalo-Rodriguez I, Roqué I FM, et al. Plasma

- interleukin-6 concentration for the diagnosis of sepsis in critically ill adults. Cochrane Database Syst Rev. 2019 Apr 30;4(4):CD011811.
- Marti AL, Moreno CA. Prognostic value of interleukin-6 for death in patients with sepsis. Med Clin (Barc). 2017 Dec 7;149(11):507.
 - Miguel-Bayarri V, Casanoves-Laparra EB, Pallás-Beneyto L, et al. Prognostic value of the biomarkers procalcitonin, interleukin-6 and C-reactive protein in severe sepsis. 2012 Nov;36(8):556-62.
 - Jekarl DW, Lee SY, Lee J, et al. Procalcitonin as a diagnostic marker and IL-6 as a prognostic marker for sepsis. Diagn Microbiol Infect Dis 2013;75(4):342-347.
 - Madhok, R. et al. Serum interleukin 6 levels in rheumatoid arthritis: correlations with clinical and laboratory indices of disease activity. Ann Rheum Dis 1993; 52: 232.
 - Yao X, Huang J, Zhong H, et al. Targeting interleukin-6 in inflammatory autoimmune diseases and cancers. Pharmacol Ther. 2014 Feb;141(2):125-39.
 - CLSI. EP15-A3: Vol. 34 No. 12, User Verification of Precision and Estimation of Bias; Approved Guideline—Third Edition.
 - CLSI. EP06-A: Vol. 23, No. 16, Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline.
 - Boscato Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
 - Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
 - Bjerner J, et al. Immunoetric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.

© 2022-2026 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Todos los derechos reservados.
Fabricante: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Dirección: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China
Dirección de correo electrónico: service@mindray.com
Sitio web: www.mindray.com
Tel.: +86-755-81888998
Fax: +86-755-26582680
Representante en la CE: Shanghai International Holding Corp. GmbH (Europe)
Dirección: Eiffestraße 80, 20537 Hamburg, Germany
Tel.: 0049-40-2513175
Fax: 0049-40-255726

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Riferimenti

- Yamasaki, K., Taga, T., Hirata, Y. et al. Cloning and expression of the human interleukin-6 (BSF-2/IFN beta 2) receptor. Science 1988; 241:825.
- Song M, Kellum JA. Interleukin-6. Crit Care Med 2005;33(Suppl12):463-465.
- Hirano, T. Interleukin-6 and its relation to inflammation and disease. Clin Immunol And Immunopathol 1992; 62:S60
- Oda S, Hirasawa H, Shiga H, et al. Sequential measurement of IL-6 blood levels in patients with systemic inflammatory response syndrome (SIRS)/sepsis. Cytokine 2005;29:169-175.
- Jawa RS, Anillo S, Huntton K, et al. Interleukin-6 in Surgery, Trauma, and Critical Care Part II: Clinical Implications. J Intensive Care Med 2011;26(2):73-87.
- Bloos F, Reinhart K. Rapid diagnosis of sepsis. Virulence.

- 2014 Jan 1;5(1):154-60.
- Molano FD, Arevalo-Rodriguez I, Roqué I FM, et al. Plasma interleukin-6 concentration for the diagnosis of sepsis in critically ill adults. Cochrane Database Syst Rev. 2019 Apr 30;4(4):CD011811.
- Marti AL, Moreno CA. Prognostic value of interleukin-6 for death in patients with sepsis. Med Clin (Barc). 2017 Dec 7;149(11):507.
- Miguel-Bayarri V, Casanoves-Laparra EB, Pallás-Beneyto L, et al. Prognostic value of the biomarkers procalcitonin, interleukin-6 and C-reactive protein in severe sepsis. 2012 Nov;36(8):556-62.
- Jekarl DW, Lee SY, Lee J, et al. Procalcitonin as a diagnostic marker and IL-6 as a prognostic marker for sepsis. Diagn Microbiol Infect Dis 2013;75(4):342-347.
- Madhok, R. et al. Serum interleukin 6 levels in rheumatoid arthritis: correlations with clinical and laboratory indices of disease activity. Ann Rheum Dis 1993; 52: 232.
- Yao X, Huang J, Zhong H, et al. Targeting interleukin-6 in inflammatory autoimmune diseases and cancers. Pharmacol Ther. 2014 Feb;141(2):125-39.
- CLSI. EP15-A3: Vol. 34 No. 12, User Verification of Precision and Estimation of Bias; Approved Guideline—Third Edition.
- CLSI. EP06-A: Vol. 23, No. 16, Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline.
- Boscatto Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
- Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037–1038.
- Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613–621.

© 2022-2026 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Tutti i diritti riservati.
Produttore: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Indirizzo: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China
Indirizzo e-mail: service@mindray.com
Sito Web: www.mindray.com
Tel.: +86-755-81888998
Fax: +86-755-26582680
Rappresentante CE: Shanghai International Holding Corp. GmbH (Europe)
Indirizzo: Eiffestraße 80, 20537 Hamburg, Germany
Tel.: 0049-40-2513175
Fax: 0049 -40 -255726

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Referanslar

- Yamasaki, K., Taga, T., Hirata, Y. et al. Cloning and expression of the human interleukin-6 (BSF-2/IFN beta 2) receptor. Science 1988; 241:825.
- Song M, Kellum JA. Interleukin-6. Crit Care Med 2005;33(Suppl12):463-465.
- Hirano, T. Interleukin-6 and its relation to inflammation and disease. Clin Immunol And Immunopathol 1992; 62:S60
- Oda S, Hirasawa H, Shiga H, et al. Sequential measurement of IL-6 blood levels in patients with systemic inflammatory response syndrome (SIRS)/sepsis. Cytokine 2005;29:169-175.
- Jawa RS, Anillo S, Huntoon K, et al. Interleukin-6 in Surgery, Trauma, and Critical Care Part II: Clinical Implications. J Intensive Care Med 2011;26(2):73-87.
- Bloos F, Reinhart K. Rapid diagnosis of sepsis. Virulence. 2014 Jan 1;5(1):154-60.
- Molano FD, Arevalo-Rodriguez I, Roqué I FM, et al. Plasma interleukin-6 concentration for the diagnosis of sepsis in critically ill adults. Cochrane Database Syst Rev. 2019 Apr 30;4(4):CD011811.
- Marti AL, Moreno CA. Prognostic value of interleukin-6 for death in patients with sepsis. Med Clin (Barc). 2017 Dec 7;149(11):507.
- Miguel-Bayarri V, Casanoves-Laparra EB, Pallás-Beneyto L, et al. Prognostic value of the biomarkers procalcitonin, interleukin-6 and C-reactive protein in severe sepsis. 2012 Nov;36(8):556-62.

- Jekarl DW, Lee SY, Lee J, et al. Procalcitonin as a diagnostic marker and IL-6 as a prognostic marker for sepsis. Diagn Microbiol Infect Dis 2013;75(4):342-347.
- Madhok, R. et al. Serum interleukin 6 levels in rheumatoid arthritis: correlations with clinical and laboratory indices of disease activity. Ann Rheum Dis 1993; 52: 232.
- Yao X, Huang J, Zhong H, et al. Targeting interleukin-6 in inflammatory autoimmune diseases and cancers. Pharmacol Ther. 2014 Feb;141(2):125-39.
- CLSI. EP15-A3: Vol. 34 No. 12, User Verification of Precision and Estimation of Bias; Approved Guideline—Third Edition.
- CLSI. EP06-A: Vol. 23, No. 16, Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline.
- Boscato Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
- Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037–1038.
- Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613–621.

© 2022-2026 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

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Üretici: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Adres: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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- receptor. Science 1988; 241:825.
- Song M, Kellum JA. Interleukin-6. Crit Care Med 2005;33(Suppl12):463-465.
- Hirano, T. Interleukin-6 and its relation to inflammation and disease. Clin Immunol And Immunopathol 1992; 62:S60
- Oda S, Hirasawa H, Shiga H, et al. Sequential measurement of IL-6 blood levels in patients with systemic inflammatory response syndrome (SIRS)/sepsis. Cytokine 2005;29:169-175.
- Jawa RS, Anillo S, Huntoon K, et al. Interleukin-6 in Surgery, Trauma, and Critical Care Part II: Clinical Implications. J Intensive Care Med 2011;26(2):73-87.
- Bloos F, Reinhart K. Rapid diagnosis of sepsis. Virulence. 2014 Jan 1;5(1):154-60.
- Molano FD, Arevalo-Rodriguez I, Roqué I FM, et al. Plasma interleukin-6 concentration for the diagnosis of sepsis in critically ill adults. Cochrane Database Syst Rev. 2019 Apr 30;4(4):CD011811.
- Marti AL, Moreno CA. Prognostic value of interleukin-6 for death in patients with sepsis. Med Clin (Barc). 2017 Dec 7;149(11):507.
- Miguel-Bayarri V, Casanoves-Laparra EB, Pallás-Beneyto L, et al. Prognostic value of the biomarkers procalcitonin, interleukin-6 and C-reactive protein in severe sepsis. 2012 Nov;36(8):556-62.
- Jekarl DW, Lee SY, Lee J, et al. Procalcitonin as a diagnostic marker and IL-6 as a prognostic marker for sepsis. Diagn Microbiol Infect Dis 2013;75(4):342-347.
- Madhok, R. et al. Serum interleukin 6 levels in rheumatoid arthritis: correlations with clinical and laboratory indices of disease activity. Ann Rheum Dis 1993; 52: 232.
- Yao X, Huang J, Zhong H, et al. Targeting interleukin-6 in inflammatory autoimmune diseases and cancers. Pharmacol Ther. 2014 Feb;141(2):125-39.
- CLSI. EP15-A3: Vol. 34 No. 12, User Verification of Precision and Estimation of Bias; Approved Guideline—Third Edition.
- CLSI. EP06-A: Vol. 23, No. 16, Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline.
- Boscato Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-32.
- Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
- Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.

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Fabricant : Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Représentant en Europe : Shanghai International Holding Corp. GmbH (Europe)