sCD14-ST Control

Order Information

Product	Catalog No.	Package Size
sCD14-ST Control (L)	105-023924-00	3×0.5 mL/vial
	105-023932-00	6×0.5 mL/vial
sCD14-ST Control (H)	105-023925-00	3×0.5 mL/vial
	105-023933-00	6×0.5 mL/vial

Intended Use

Mindray sCD14-ST Control is used for quality control by monitoring the accuracy and precision of Mindray CL-series Chemiluminescence Immunoassay Analyzer and test ability of clinical laboratory in the quantitative measurement of soluble CD14 subtype analyte.

Summary

Mindray CL measurement system is composed of Mindray CL-series Chemiluminescence Immunoassay Analyzer, Mindray reagent kits, calibrators, and controls. Mindray sCD14-ST control can be used to check whether QC result is within specified range or not. It is used with Mindray sCD14-ST (CLIA).

Components

sCD14-ST Control contains soluble CD14 subtype analyte in buffered bovine serum albumin (BSA) with 0.09% sodium azide and 0.05% ProClin 300 as the preservative.

Warnings and Precautions

- For in vitro diagnostic use only. For laboratory professional use.
- Laboratory professionals must be trained to take all necessary precautions while handling laboratory reagents. Do not use reagents beyond the expiration date.
- The control values are lot-specific with the matched set listed in the value sheet.
- The controls should be run daily in parallel with the patient samples, after each calibration, or after the reagent lot has been changed.
- Please perform maintenance and standard operation as scheduled, including calibration and system analysis to assure the performance of measurement.
- The controls were tested with CE-approved kits for HIV-1/2 Antibody, HCV Antibody, Treponema pallidum Antibody and HBsAg. The results are all negative. However, as no testing method can rule

out the potential risk of infection with absolute certainty, this material should be handled as potential biohazard.¹

- If there is any evidence of microbial contamination or excessive turbidity in the product, contact Mindray Customer Service for help.
- Material safety data sheet is available upon request.
- Disposal of all waste material should be in accordance with local guidelines.
- The target values and ranges of Mindray sCD14-ST Control are for reference only. Each laboratory can establish its own internal quality control procedure and set up its own target values and ranges.
- Please confirm the integrity of the package before use. Do not use the controls with damaged packages. Inappropriate storage will lead to incorrect results.
- If the controls are opened unintentionally before use, they shall be tightly capped and stored at 2-8°C, and then used within in-use stability time.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the local competent authority.

Test Method

sCD14-ST controls have similar components as test samples. Therefore, the controls can be tested by following the same procedure as tested samples, the target values and ranges are then calculated and compared to see if the controls meet specifications.

Preparation

- Refer to Mindray Operator's Manual and /or Customer Service to acquire information about control principle, control use, and control validity evaluation.
- Take the vial out of fridge, and let it stand for no less than 30 minutes at room temperature.
- Mix contents by gently inverting the vial at least 10 times before use. Avoid bubble formation.
- Because of possible evaporation effects, the control on the analyzers should be measured as soon as possible.
- Run QC test you have requested or have it analyzed in the same way as patient samples.

Storage and Stability

sCD14-ST Control is stable up to the expiration date as indicated on the label when stored in unopened vial at 2-8°C. Once opened, it is stable for 30 days at 2-8°C or 90 days at -20°C.

Materials Required but not Provided

- Mindray CL-series Chemiluminescence Immunoassay Analyzer, reagent kits, and matched calibrators.
- General laboratory equipment.

Assay Values

The control values (the target value and range) determined by Mindray routine measurement procedure are listed in the target value sheet. The target value was obtained from Mindray measurement system and the range was calculated as the target values ± 3 standard deviations. The control values are lot-specific. Please check the lot number of each control before use.

Quality Control

Use Mindray controls as the samples to monitor the performance of system. The values of controls are lot-specific on Mindray CL-series Chemiluminescence Immunoassay Analyzer. The control result should be within the defined range as shown in the control target value sheet. If the controls are out of specified range, the measurement system and reagents should be checked, for instance, expiration date or storage condition of the calibrators, reagent and controls; possible contamination of reagent; position of reagent, sample placement on the analyzer; the parameter setting in software; or operation status of the analyzer. Each laboratory should establish its own internal quality control scheme and procedure for corrective and preventive actions.

Graphical Symbols

In vitro diagnostic medical device

Authorized representative in the European Community

Consult instructions for use

Catalogue number

Temperature limit

Manufacturer

Biological risks

Batch code

This way up

Use-by date

Caution

European conformity

CONTROL

UDI

Unique device identifier

References

- HHS Publication, 5th ed., December 2009. Biosafety in Microbiological and Biomedical Laboratories.
© 2022-2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
All rights Reserved.

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Address: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P.R. China
E-mail Address: service@mindray.com
Website: www.mindray.com
Tel: +86-755-81888998
Fax: +86-755-26582680
EC-Representative: Shanghai International Holding Corp. GmbH (Europe)
Address: Eiffestraße 80, Hamburg 20537, Germany
Tel: 0049-40-2513175
Fax: 0049-40-255726

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Fabricant : Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Représentant en Europe : Shanghai International Holding Corp. GmbH (Europe)