CMV IgM

Cytomegalovirus IgM (CLIA)

INSTRUCTIONS FOR USE

mindray

CMV IgM

Cytomegalovirus IgM (CLIA)

Order Information

Catalog No.	Package Size
105-012542-00	2×50 tests (calibrators not included)
105-012543-00	2×100 tests (calibrators not included)
105-012545-00	2×50 tests (calibrators included)
105-012546-00	2×100 tests (calibrators included)

Intended Use

The CL-series Cytomegalovirus IgM assay is a Chemiluminescent Immunoassay (CLIA) for the qualitative determination of Cytomegalovirus IgM in human serum and plasma. The test is intended as an aid in the diagnosis of recent or current Cytomegalovirus infection and in the assessment of immune status of individuals including pregnant women.

Summary

Cytomegalovirus (CMV), a member of the herpes virus family, is ubiquitous in all human populations, causing infections which are followed by life-long latency in the host with occasional reactivations. Serological survey results show that 90% of the population in Asia and Africa have been infected with cytomegalovirus. The serological positive rate increases with age, and there is no significant difference between men and women. Horizontal transmission is the most common route of infection of cytomegalovirus, which is mainly caused by contact with viral secretions, including saliva, urine, feces, cervical or vaginal secretions, semen, etc. In addition, cytomegalovirus can also spread vertically through placenta or cause iatrogenic infection^[1-2].

Most normal adults and children do not have obvious clinical symptoms after infection with cytomegalovirus. After infection, antibodies can be produced, but the virus can not be completely eliminated. Cytomegalovirus can enter tissue cells through fusion with cell membrane, showing a long-term carrying state. When the immune function of the body is damaged, such as receiving immunosuppressants, organ transplantation or combined with human immunodeficiency virus infection, the virus can be reactivated and replicated, and there are clinical manifestations such as pneumonia, hepatitis, peptic ulcer and central nervous system damage^[3-4].

During pregnancy, cytomegalovirus could vertically spreads to the

fetus through the placenta, which can cause congenital cytomegalovirus infection, and this is the most common congenital virus infection. Infected fetus may not have clinical manifestations, but can also suffer from systemic visceral damage due to viremia, such as jaundice, hepatosplenomegaly, fetal growth restriction, cerebellar malformation, respiratory distress syndrome, etc. Most patients die within 2 months after birth or are left with serious neurological disorders, such as deafness, intellectual deficiency or dyskinesia. The severity of congenital cytomegalovirus infection is related to the type of maternal infection. In general, fetal congenital infection caused by maternal primary infection is more serious than reactivation or reinfection^[5-6].

Detection of cytomegalovirus specific IgG and IgM antibodies is a common method for the diagnosis of acute cytomegalovirus infection. Cytomegalovirus IgM antibody generally appeared after 1 week of infection and gradually disappeared after 4-8 weeks, which is an indicator of recent infection; Cytomegalovirus IgG antibody generally appears after 2 weeks of infection and lasts for many years. It is an indicator of previous infection. The significant increase of cytomegalovirus IgG after an interval of 3-4 weeks supports the diagnosis of acute cytomegalovirus infection^[7].

Cytomegalovirus IgM antibody test should not be used in the antenatal screening of the asymptomatic population. The test results of this reagent should not be used as the basis for termination of pregnancy.

Assay Principles

The CL-series CMV IgM assay is a μ -capture assay to qualitatively detect IgM antibodies to Cytomegalovirus.

In the first step, sample, sample treatment solution and paramagnetic microparticles coated with anti-human IgM monoclonal antibody are added into a reaction cuvette. After incubation, CMV IgM antibodies in the sample will bind to anti-human IgM monoclonal antibody coated microparticles. Afterwards, microparticles are magnetically captured while other unbound substances are removed by washing.

In the second step, test diluent solution and ALP labeled CMV antigens are added to the reaction cuvette. After incubation, ALP labeled CMV antigens will form sandwich structure with microparticles captured CMV IgM antibodies. Afterwards, microparticles are magnetically captured while other unbound substances are removed by washing.

Then, the substrate solution is added to the reaction cuvette. It is catalyzed by CMV antigens -ALP conjugates in the immune-complex retained on the microparticles. The resulting chemiluminescent reaction is measured as relative light units (RLUs) by a photomultiplier in the instrument.

A direct relationship exists between the amount of CMV IgM antibodies in the sample and the RLU's generated during the reaction. The presence or absence of analytes in the sample is determined by comparing the chemiluminescent signal in the reaction to the cutoff signal determined from calibration. A cutoff index (COI) is calculated by Sample RLU/Cutoff RLU.

Reagent Components

The kit is composed of Ra, Rb, Rc, and Rd, as well as calibrators C0 and C1. Calibrators C0 and C1 are optional (refer to the "Order Information" section). Components in different lots of reagent kit cannot be interchanged and the detailed information of each component is listed below:

Ra	Paramagnetic microparticles coated with anti-human IgM monoclonal antibody in Tris buffer. Minimum concentration: 0.3 mg/ml solids. Tris buffer: 50 mmol/L. Preservative: 0.048% Proclin 300. 3.8 mL/50 tests; 6.6 mL/100 tests.
Rb	ALP labeled Cytomegalovirus specific antigens in MES buffer. Minimum concentration: 2 ng/ml. MES buffer: 50 mmol/L. Preservative: 0.048% Proclin 300. 3.5 mL/50 tests; 6.3 mL/100 tests.
Rc	Sample treatment solution. 25 mmol/L Tris buffer. Preservative: 0.048% Proclin 300. 3.5 mL/50 tests; 6.3 mL/100 tests.
Rd	Test diluent solution. 25 mmol/L Tris buffer. Preservative: 0.048% Proclin 300. 5.3 mL/50 tests; 9.8 mL/100 tests.
Calibrators C0 and C1 (optional)	C0: IgM antibodies to Cytomegalovirus negative in buffer with 0.048% Proclin 300; C1: IgM antibodies to Cytomegalovirus positive in human serum with 0.1% Proclin 300. C0: 1×0.6 mL/vial; C1: 1×0.5 mL/vial.

Storage and Stability

The unopened reagents (Ra, Rb, Rc and Rd) are stable up to expiration date indicated on the label, when stored at 2-8°C. They can be stored onboard and used for a maximum of 28 days at 2-8°C after opening.

The calibrators C0 and C1 are stable up to the expiration date as indicated on the label when stored unopened at 2-8°C. Once opened in use, the calibrators are stable for 28 days when

refrigerated at 2-8°C or 90 days when frozen at -20°C (Frozen only once).

Materials Required but not Provided

Mindray CL-series Chemiluminescence Immunoassay Analyzer

Mindray ToRCH IgM Positive Control:

1×3.0 mL/vial (REF: 105-024840-00),

3×3.0 mL/vial (REF: 105-024841-00),

3×1.0 mL/vial (REF: 105-024842-00),

6×1.0 mL/vial (REF: 105-024843-00).

Mindray ToRCH IgG/IgM Negative Control:

1×3.0 mL/vial (REF: 105-024824-00),

3×3.0 mL/vial (REF: 105-024825-00),

3×1.0 mL/vial (REF: 105-024826-00),

6×1.0 mL/vial (REF: 105-024827-00).

Mindray Substrate Solution: 4x115 mL (REF: 105-004274-00).

Mindray Substrate Solution: 4x75 mL (REF: 105-009044-00).

Mindray Wash Buffer, 1×10 L (REF: 105-004552-00).

Mindray Reaction cuvette.

Instrument System

Mindray CL-series Chemiluminescence Immunoassay Analyzers

Specimen Collection and Preparation**Specimen Types**

- Human serum and plasma (K₂-EDTA, K₃-EDTA, Na₂-EDTA, sodium heparin, lithium heparin and sodium citrate) are recommended for this assay.
- Blood collection tubes from various manufacturers may contain different raw materials and additives which could affect the test results in some cases. Not all available tubes of all manufacturers were tested by Mindray. Each laboratory needs to judge the applicability of the tubes and plasma separation products. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Specimen Conditions

- Do not use:
 - heat-inactivated specimens
 - grossly hemolyzed specimens
 - specimens with obvious microbial contamination
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.

Preparation for Analysis

- Follow blood collection tube manufacturer's recommendations for centrifugation. Centrifuge the specimens after clot formation is complete. Ensure residual fibrin and cellular matter has been removed prior to analysis.

- For optimal results, inspect all samples for bubbles. Remove bubbles with a pipette tip prior to analysis. Specimens must be mixed thoroughly after thawing. Thawed samples should be centrifuged prior to use.
- If the sample was covered with lipid layer after centrifugation, the sample should be transferred to a clean tube and centrifuged before testing. Do not transfer the lipid layer. Handle carefully to prevent cross contamination.

Specimen Storage

- Specimens should be tested as soon as possible after sample collection. If testing is not completed within 48 hours, specimens should be tightly capped and refrigerated at 2-8°C. The sample is stable for 7 days at 2-8°C and 180 days at -20°C or below.
- The samples may be frozen no more than 5 times.
- Frozen or refrigerated samples, controls and calibrators must be restored to room temperature before analysis.

Assay Procedure

For optimal performance of this assay, operators should read the related Operator's Manual carefully to get enough information such as operation instructions, sample preservation and management, safety precautions, and maintenance. Prepare all required materials for the assay as well.

Before loading the CMV IgM reagent kit on the instrument for the first time, invert unopened reagent bottle gently for at least 30 times to resuspend the microparticles, which have settled during shipment or storage. Visually inspect the bottle to ensure the microparticles have been well mixed. If the microparticles remain adhered to the bottle, continue inverting until suspended completely. Avoid foam formation. If the microparticles cannot be homogenized, it is recommended not to use this bottle of reagent. Contact Mindray Customer Service for help. Avoid inverting opened reagent bottle, and avoid air bubbles.

This assay requires 10 µL of sample volume for a single test not including dead volume of the sample container. Additional volume is required when performing additional tests from the same sample. Operators should refer to the Operator's Manual and follow the specific requirements of the assay to determine the minimum sample volume.

Calibration

The calibrators are traceable to Mindray internal reference.

The calibration information is stored in the barcodes attached into the reagent and calibrator pack. When performing the calibration, please scan the barcodes into the system.

Take the vial out of refrigerator and let it stand for no less than 30 minutes at room temperature. Mix contents by gently inverting the vial no less than 10 times before use. Avoid bubble formation. Test

C0 three times and C1 twice for calibration. The original calibrator vials may be used if the entire provided volume is used. Because of possible evaporation effects, the calibrator on the analyzer should be measured as soon as possible and not more than 3 calibration procedures per original vial should be performed. After use, close the bottles immediately and store upright at 2-8°C. The original calibrators can also be transferred into empty labeled bottles in aliquots. Use the aliquots immediately or store them immediately at -20°C. One aliquot can be used only once.

A valid calibration is required before any CMV IgM testing.

Renewed calibration is required as follows:

1. The previous calibration was performed 28 days before.
2. Quality control results lie outside of the expected ranges.
3. A reagent kit with a new lot number is used.

For detailed instruction of calibration, refer to the system Operator's Manual.

Quality Control

Users can prepare quality controls with clinical samples. Reference ranges can be established according to protocol approved by individual laboratories.

It is recommended that quality controls should be run once every 24 hours if the tests are in use, and following each calibration. The quality control frequency should be adapted to each laboratory's requirements.

Quality control results should be within the acceptable ranges. If a control is out of its specified range, the associated test results are invalid and the samples must be retested. Recalibration may be required. Refer to the system Operator's Manual to check up the system. If quality control findings are still outside of the identified ranges, please contact Mindray Customer Service for help.

Calculation

The CL-series analyzer calculates the mean RLU signals of both calibrator C1 and C0.

The CL-series analyzer calculates CMV IgM result based on the ratio of the sample RLU to the cutoff RLU (Cutoff Index, COI) for each specimen and control.

$$\text{Cutoff RLU} = [(\text{Mean RLU of C1} - \text{Mean RLU of C0}) \times \text{Calibration Coefficient}] + \text{Mean RLU of C0}$$

Calibration Coefficient is specific for each reagent and calibrator lot.

$$\text{COI} = \text{Sample RLU} / \text{Cutoff RLU}$$

Specimen Dilution

Specimens cannot be diluted for Mindray CMV IgM assay.

Interpretation of Results

Specimens with a COI < 0.80 are considered to be non-reactive for IgM antibodies to CMV.

Specimens with a $COI \geq 0.80$ to < 1.00 are considered to be indeterminate for IgM antibodies to CMV. A second sample should be collected for testing within an appropriate time (e.g. 2-3 weeks) or other test methods should be used for confirmation.

Specimens with a $COI \geq 1.00$ are considered to be reactive for IgM antibodies to CMV, suggesting recent infection.

The results of different detection systems may be inconsistent because of the detection methods and traceability.

The assay results should not be used solely for confirmation or exclusion of CMV infection. Clinical decisions should be made in conjunction with other evidences, such as symptoms, clinical history, results of nucleic acid test, etc.

Limitation of Measurement

A negative CMV IgM test result cannot completely exclude the possibility of acute CMV infection. If acute infection is suspected, a second sample should be collected for testing within an appropriate time (e.g. 2-3 weeks). If the second sample shows seroconversion of primary infection or significant increase of CMV IgG, the diagnosis of acute CMV infection is supported. It is recommended that the first sample should be stored at -20°C or below, so that two samples can be tested at the same time and same conditions^[8].

For patients with immunodeficiency disease or receiving immunosuppressive therapy, the antibody test results should be interpreted with caution.

Exogenous interference should be considered in case a patient with CMV IgM positive result has recently received blood transfusion.

CMV IgM is not only increased in the primary infection, but also in the recurrent infection and reinfection, elevated IgM may persist even for years after primary infection in some individuals. Therefore, IgM level alone cannot be used to determine the time of initial infection.

Due to the existence of heterophilic antibodies, patients who are often exposed to animals or use animal derived immunoglobulin for treatment may have abnormal detection results ^[9].

The samples of patients who have been treated with mouse monoclonal antibody preparation may contain human anti mouse antibody (HAMA), which may interfere with the test results of kits containing mouse monoclonal antibodies ^[10-11]. Mindray chemiluminescence kit contains anti-interference components, which can effectively reduce the impact of HAMA in samples, but the interference of HAMA may still exist in a few samples. Only combined with the clinical information or other laboratory test results can the patient's condition be correctly determined.

Specimens from patients with high levels of non-specific IgM, e.g.

from patients with multiple myeloma, may show depressed values when tested with μ -capture format assays.

Performance Characteristics

Precision

Precision was determined by following Clinical & Laboratory Standards Institute (CLSI) protocol EP05-A3^[12]. 3 serum samples and 2 controls were tested 2 runs each day in duplicate for totally 20 days (n=80). Results are shown in the table below*.

Sample	Mean (COI)	Repeatability		Intermediate precision	
		SD (COI)	CV (%)	SD (COI)	CV (%)
S1	0.18	0.00	N/A**	0.02	N/A
S2	0.94	0.02	1.92%	0.04	4.10%
S3	2.64	0.06	2.10%	0.11	4.24%
NC	0.08	0.00	N/A	0.01	N/A
PC	2.61	0.04	1.35%	0.09	3.59%

*Representative data; results in individual laboratories may vary.

**N/A = Not Applicable.

Hook effect

When a sample containing a high concentration of determinand is tested, the hook effect may lead to false negative results. This assay adopts two-steps mode to keep the analysis signal at a high value, which could effectively avoid false negative results caused by hook effect.

3 samples with high concentration of CMV IgM were detected, all results showed the expected high concentration and there was no misjudgement of results.

Interference

The test results of Mindray CMV IgM are not observed to be interfered with the potential endogenous substances up to the following concentrations (criterion: recovery within $\pm 15\%$ of initial value or absolute deviation within ± 0.15 COI)*.

Substance	Concentration
Hemoglobin	1000 mg/dL
Bilirubin	40 mg/dL
Triglyceride	37 mmol/L
Total protein	12 g/dL
Albumin	6 g/dL
γ -globulin	6 g/dL
Biotin	3500 ng/mL

*Representative data; results in individual laboratories may vary.

The test results of Mindray CMV IgM are not observed to be interfered with the potential exogenous substances up to the following concentrations (criterion: recovery within $\pm 15\%$ of initial value or absolute deviation within ± 0.15 COI)*.

Substance	Concentration
Acetaminophen	200 mg/L

Substance	Concentration
Amoxicillin	700 mg/L
Cefaclor	500 mg/L
Vitamin C	250 mg/L
Aspirin	1000 mg/L
Ibuprofen	500 mg/L
Ganciclovir	40 mg/L
Valganciclovir	40 mg/L
Foscarnet sodium	400 mg/L

**Representative data; results in individual laboratories may vary.*

Cross-reactivity

Studies were performed to evaluate the effect of potential cross reacting diseases on CMV IgM assay. A total of 265 specimens were tested by Mindray CMV IgM and a comparator CMV IgM kit.

Substance	n
RV	4
Toxo	3
HSV-1	2
HSV-2	2
E-B virus	19
Parvovirus B19	3
HBV	61
HAV	4
HCV	2
HIV	20
HTLV	2
VZV	4
Measles virus	2
Mumps virus	2
Borrelia burgdorferi	2
Influenza A virus	5
Influenza B virus	4
Parainfluenza virus	4
TP	17
E.coli	2
AMA	4
RF	40
ANA	51
HAMA	2
SLE	2
CMV IgG	2

The results showed that the overall agreement of Mindray CMV IgM and commercial CMV IgM kits is 95.09%.

Method comparison

1032 clinical specimens from pregnant women, women of childbearing age, neonates and patients with corresponding clinical symptoms were tested to evaluate the coincidence rate of Mindray CMV IgM and a mainstream commercial chemiluminescence CMV IgM kit. The positive coincidence rate was 96.67% (116/120, 95%CI 91.74%-98.70%); the negative coincidence rate was 99.78% (901/903, 95%CI 99.20%-99.94%); the total coincidence rate was 99.41% (1017/1023 95%CI 98.73%-99.73%). Specimens with inconsistent test results were confirmed by another commercial CMV IgM kit.

Mindray CMV IgM	Comparator CMV IgM			Total
	Reactive	Non-reactive	Indeterminate	
Reactive	116	2*	4	122
Non-reactive	4**	901	4	909
Indeterminate	0	0	1	1
Total	120	903	9	1032

*: Of 2 specimens with reactive on Mindray CMV IgM and non-reactive on comparator CMV IgM, 2 specimens were found non-reactive on another commercial CMV IgM kit.

** : Of 4 specimens with non-reactive on Mindray CMV IgM and reactive on comparator CMV IgM, 4 specimens were found non-reactive on another commercial CMV IgM kit.

Seroconversion panels

To evaluate the sensitivity of seroconversion, 6 commercial seroconversion panels were tested by Mindray CMV IgM and a mainstream commercial CMV IgM kit, the results are shown in the table below.

Panel ID	Timepoint of seroconversion (days)*	
	Mindray CMV IgM	Comparator CMV IgM
SCP-CMV-001	8	8
SCP-CMV-002	36	43
SCP-CMV-003	25	25
SCP-CMV-005	6	6
SCP-CMV-006	6	6
SCP-CMV-007	27	27

*: Timepoint of seroconversion: time from the first blood collection after positive conversion.

Warnings and Precautions

- For *in vitro* diagnostic use only. For laboratory professional use.
- Follow all the rules in handling laboratory reagents and take necessary safety precautions.
- Due to the differences in methodology and antibody specificity, test results of the same sample may be different when using reagent kits from different

- manufacturers on Mindray system, or using Mindray reagent kits on other systems.
4. Do not use reagent kits beyond the expiration date.
 5. Do not use reagents mixed from different reagent lots.
 6. Please confirm the integrity of the package before use. Do not use the reagents with damaged packages. The results can't be assured when the reagents are stored at inappropriate condition.
 7. If the reagents are opened unintentionally before use, they shall be used as soon as possible.
 8. Wear clean gloves to avoid contaminating the reagents.
 9. Always keep the reagent pack in the upright position to ensure no microparticle has been lost prior to use.
 10. The calibrator was tested with CE-marked methods and shown to be negative for antibodies to HIV and HCV and nonreactive for HBsAg. However, as no testing method can rule out the potential risk of infection with absolute certainty, so this material should be handled as a patient specimen.
 11. Reagent pack opened for more than 14 days is not recommended for use.
 12. Instability or deterioration should be suspected if there are visible signs of leakage, turbidity, precipitates or microbial growth.
 13. Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
 14. All the specimens and reaction wastes should be considered potentially biohazard. Specimens and reaction wastes should be handled in accordance with the local regulations and guidelines.
 15. The Material Safety Data Sheet (MSDS) is available upon request.
 16. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the local competent authority.
 17. This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008.

**Warning**

H317 May cause an allergic skin reaction.

H412 Harmful to aquatic life with long lasting effects.

Prevention:

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P272 Contaminated work clothing should not be allowed out of the workplace.

P273 Avoid release to the environment.

P280 Wear protective gloves/protective clothing/eye

protection/face protection.

Response:

P302 + P352 IF ON SKIN: Wash with plenty of water.

P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.

P362 + P364 Take off contaminated clothing and wash it before reuse.

Disposal:

P501 Dispose of contents/ container in accordance with local regulation.

Graphical Symbols

In vitro
diagnostic
medical
device

Authorized
Representative
in the European
Community

Consult
instructions
for use

European
Conformity

Catalogue
number

Temperature
limit

Manufacturer

Use-by
date

Biological
risks

Batch code

This way up

Unique
device
identifier

Caution

Bibliography

1. Cannon MJ, Schmid DS, Hyde TB. Review of cytomegalovirus seroprevalence and demographic characteristics associated with infection. Rev Med Virol 2010;20:202-213.

2. Johnson J, Anderson B, Pass RF. Prevention of maternal and congenital cytomegalovirus infection. Clin Obstet Gynecol 2012;55:521-530.
3. Yinon Y, Farine D, Yudin MH. Screening, diagnosis, and management of cytomegalovirus infection in pregnancy. Obstet Gynecol Surv 2010;65:736-743.
4. Buonsenso D, Serranti D, Gargiullo L, et al. Congenital cytomegalovirus infection: current strategies and future perspectives. Eur Rev Med Pharmacol Sci 2012;16:919-935.
5. Coll O, Benoist G, Ville Y, et al. WAPM Perinatal Infections Working Group. Guidelines on CMV congenital infection. J Perinat Med 2009;37:433-445.
6. Joseph S, Beliveau C, Muecke C, et al. Cytomegalovirus as an occupational risk in daycare educators. Paediatr Child Health 2006;11:401-407.
7. Lazzarotto T, Guerra B, Lanari M, et al. New advances in the diagnosis of congenital cytomegalovirus infection. J Clin Virol 2008;41:192-197.
8. Revello M G, Gerna G. Diagnosis and management of human cytomegalovirus infection in the mother, fetus, and newborn infant. [J]. Clinical microbiology reviews, 2002, 15(4):680-715.
9. Boscatto Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
10. Kricka L. Interferences in immunoassays—still a threat. Clin Chem 2000; 46: 1037-1038.
11. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
12. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline –Third Edition.

© 2022-2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
All rights Reserved.

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

E-mail Address: service@mindray.com

Website: www.mindray.com

Tel: +86-755-81888998

Fax: +86-755-26582680

EC-Representative: Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffeustraße 80, Hamburg 20537, Germany

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1. Cannon MJ, Schmid DS, Hyde TB. Review of cytomegalovirus seroprevalence and demographic characteristics associated with infection. Rev Med Virol 2010;20:202-213.
2. Johnson J, Anderson B, Pass RF. Prevention of maternal and congenital cytomegalovirus infection. Clin Obstet Gynecol 2012;55:521-530.
3. Yinon Y, Farine D, Yudin MH. Screening, diagnosis, and management of cytomegalovirus infection in pregnancy. Obstet Gynecol Surv 2010;65:736-743.
4. Buonsenso D, Serranti D, Gargiullo L, et al. Congenital cytomegalovirus infection: current strategies and future perspectives. Eur Rev Med Pharmacol Sci 2012;16:919-935.
5. Coll O, Benoist G, Ville Y, et al. WAPM Perinatal Infections Working Group. Guidelines on CMV congenital infection. J Perinat Med 2009;37:433-445.
6. Joseph S, Beliveau C, Muecke C, et al. Cytomegalovirus as an occupational risk in daycare educators. Paediatr Child Health 2006;11:401-407.
7. Lazzarotto T, Guerra B, Lanari M, et al. New advances in the diagnosis of congenital cytomegalovirus infection. J Clin Virol 2008;41:192-197.

8. Revello M G , Gema G . Diagnosis and management of human cytomegalovirus infection in the mother, fetus, and newborn infant.[J]. Clinical microbiology reviews, 2002, 15(4):680-715.
9. Boscato Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
10. Kricka L. Interferences in immunoassays—still a threat. Clin Chem 2000; 46: 1037-1038.
11. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
12. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline –Third Edition.

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Адрес: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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Atenção

Bibliografia

1. Cannon MJ, Schmid DS, Hyde TB. Review of cytomegalovirus seroprevalence and demographic characteristics associated with infection. Rev Med Virol 2010;20:202-213.
2. Johnson J, Anderson B, Pass RF. Prevention of maternal and congenital cytomegalovirus infection. Clin Obstet Gynecol 2012;55:521-530.
3. Yinon Y, Farine D, Yudin MH. Screening, diagnosis, and management of cytomegalovirus infection in pregnancy. Obstet Gynecol Surv 2010;65:736-743.
4. Buonsenso D, Serranti D, Gargiullo L, et al. Congenital cytomegalovirus infection: current strategies and future perspectives. Eur Rev Med Pharmacol Sci 2012;16:919-935.
5. Coll O, Benoist G, Ville Y, et al. WAPM Perinatal Infections Working Group. Guidelines on CMV congenital infection. J Perinat Med 2009;37:433-445.
6. Joseph S, Beliveau C, Muecke C, et al. Cytomegalovirus as an occupational risk in daycare educators. Paediatr Child Health 2006;11:401-407.
7. Lazzarotto T, Guerra B, Lanari M, et al. New advances in the diagnosis of congenital cytomegalovirus infection. J Clin Virol 2008;41:192-197.
8. Revello M G, Gema G. Diagnosis and management of human cytomegalovirus infection in the mother, fetus, and newborn infant.[J]. Clinical microbiology reviews, 2002, 15(4):680-715.
9. Boscato Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
10. Kricka L. Interferences in immunoassays—still a threat. Clin Chem 2000; 46: 1037-1038.

11. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
12. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline –Third Edition.

© 2022-2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

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Fabricante: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Endereço: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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Precaución

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1. Cannon MJ, Schmid DS, Hyde TB. Review of cytomegalovirus seroprevalence and demographic characteristics associated with infection. Rev Med Virol 2010;20:202-213.
2. Johnson J, Anderson B, Pass RF. Prevention of maternal and congenital cytomegalovirus infection. Clin Obstet Gynecol 2012;55:521-530.
3. Yinon Y, Farine D, Yudin MH. Screening, diagnosis, and management of cytomegalovirus infection in pregnancy. Obstet Gynecol Surv 2010;65:736-743.
4. Buonsenso D, Serranti D, Gargiullo L, et al. Congenital cytomegalovirus infection: current strategies and future

- perspectives. Eur Rev Med Pharmacol Sci 2012;16:919-935.
5. Coll O, Benoist G, Ville Y, et al. WAPM Perinatal Infections Working Group. Guidelines on CMV congenital infection. J Perinat Med 2009;37:433-445.
 6. Joseph S, Beliveau C, Muecke C, et al. Cytomegalovirus as an occupational risk in daycare educators. Paediatr Child Health 2006;11:401-407.
 7. Lazzarotto T, Guerra B, Lanari M, et al. New advances in the diagnosis of congenital cytomegalovirus infection. J Clin Virol 2008;41:192-197.
 8. Revello M G, Gerna G. Diagnosis and management of human cytomegalovirus infection in the mother, fetus, and newborn infant. [J]. Clinical microbiology reviews, 2002, 15(4):680-715.
 9. Boscatto Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
 10. Kricka L. Interferences in immunoassays—still a threat. Clin Chem 2000; 46: 1037-1038.
 11. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
 12. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline –Third Edition.

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Attenzione

Bibliografia

1. Cannon MJ, Schmid DS, Hyde TB. Review of cytomegalovirus seroprevalence and demographic characteristics associated with infection. Rev Med Virol 2010;20:202-213.
2. Johnson J, Anderson B, Pass RF. Prevention of maternal and congenital cytomegalovirus infection. Clin Obstet Gynecol 2012;55:521-530.
3. Yinon Y, Farine D, Yudin MH. Screening, diagnosis, and management of cytomegalovirus infection in pregnancy. Obstet Gynecol Surv 2010;65:736-743.
4. Buonsenso D, Serranti D, Gargiullo L, et al. Congenital cytomegalovirus infection: current strategies and future perspectives. Eur Rev Med Pharmacol Sci 2012;16:919-935.
5. Coll O, Benoist G, Ville Y, et al. WAPM Perinatal Infections Working Group. Guidelines on CMV congenital infection. J Perinat Med 2009;37:433-445.
6. Joseph S, Beliveau C, Muecke C, et al. Cytomegalovirus as an occupational risk in daycare educators. Paediatr Child Health 2006;11:401-407.
7. Lazzarotto T, Guerra B, Lanari M, et al. New advances in the diagnosis of congenital cytomegalovirus infection. J Clin Virol 2008;41:192-197.
8. Revello M G, Gerna G. Diagnosis and management of human cytomegalovirus infection in the mother, fetus, and newborn infant.[J]. Clinical microbiology reviews, 2002, 15(4):680-715.

9. Boscatto Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
10. Kricka L. Interferences in immunoassays—still a threat. Clin Chem 2000; 46: 1037-1038.
11. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
12. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline—Third Edition.

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Indirizzo: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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Dikkat

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1. Cannon MJ, Schmid DS, Hyde TB. Review of cytomegalovirus seroprevalence and demographic characteristics associated with infection. Rev Med Virol 2010;20:202-213.
2. Johnson J, Anderson B, Pass RF. Prevention of maternal and congenital cytomegalovirus infection. Clin Obstet Gynecol 2012;55:521-530.
3. Yinon Y, Farine D, Yudin MH. Screening, diagnosis, and management of cytomegalovirus infection in pregnancy. Obstet Gynecol Surv 2010;65:736-743.
4. Buonsenso D, Serranti D, Gargiullo L, et al. Congenital cytomegalovirus infection: current strategies and future perspectives. Eur Rev Med Pharmacol Sci 2012;16:919-935.
5. Coll O, Benoist G, Ville Y, et al. WAPM Perinatal Infections Working Group. Guidelines on CMV congenital infection. J Perinat Med 2009;37:433-445.

6. Joseph S, Beliveau C, Muecke C, et al. Cytomegalovirus as an occupational risk in daycare educators. Paediatr Child Health 2006;11:401-407.
7. Lazzarotto T, Guerra B, Lanari M, et al. New advances in the diagnosis of congenital cytomegalovirus infection. J Clin Virol 2008;41:192-197.
8. Revello M G , Gema G . Diagnosis and management of human cytomegalovirus infection in the mother, fetus, and newborn infant.[J]. Clinical microbiology reviews, 2002, 15(4):680-715.
9. Boscatto Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
10. Kricka L. Interferences in immunoassays–still a threat. Clin Chem 2000; 46: 1037-1038.
11. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
12. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline–Third Edition.

© 2022-2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

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Üretici: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Adres: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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2. Johnson J, Anderson B, Pass RF. Prevention of maternal and congenital cytomegalovirus infection. Clin Obstet Gynecol 2012;55:521-530.
3. Yinon Y, Farine D, Yudin MH. Screening, diagnosis, and management of cytomegalovirus infection in pregnancy. Obstet Gynecol Surv 2010;65:736-743.
4. Buonsenso D, Serranti D, Gargiullo L, et al. Congenital cytomegalovirus infection: current strategies and future perspectives. Eur Rev Med Pharmacol Sci 2012;16:919-935.
5. Coll O, Benoist G, Ville Y, et al. WAPM Perinatal Infections Working Group. Guidelines on CMV congenital infection. J Perinat Med 2009;37:433-445.
6. Joseph S, Beliveau C, Muecke C, et al. Cytomegalovirus as an occupational risk in daycare educators. Paediatr Child Health 2006;11:401-407.
7. Lazzarotto T, Guerra B, Lanari M, et al. New advances in the diagnosis of congenital cytomegalovirus infection. J Clin Virol 2008;41:192-197.
8. Revello M G, Gerna G. Diagnosis and management of human cytomegalovirus infection in the mother, fetus, and newborn infant.[J]. Clinical microbiology reviews, 2002, 15(4):680-715.
9. Boscatto Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
10. Kricka L. Interferences in immunoassays—still a threat. Clin Chem 2000; 46: 1037-1038.
11. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
12. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline –Third Edition.

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6. Joseph S, Beliveau C, Muecke C, et al. Cytomegalovirus as an occupational risk in daycare educators. Paediatr Child Health 2006;11:401-407.
7. Lazzarotto T, Guerra B, Lanari M, et al. New advances in the diagnosis of congenital cytomegalovirus infection. J Clin Virol 2008;41:192-197.
8. Revello M G, Gerna G. Diagnosis and management of human cytomegalovirus infection in the mother, fetus, and newborn infant.[J]. Clinical microbiology reviews, 2002, 15(4):680-715.
9. Boscatto Lm, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
10. Kricka L. Interferences in immunoassays – still a threat. Clin Chem 2000; 46: 1037-1038.
11. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
12. CLSI. EP05-A3: Vol. 34, No. 13, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline – Third Edition.

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