ZINCO

Anvisa 80115310247

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ZINC

Anvisa 80115310247

BEFORE USING THE PRODUCT, CHECK THE INSTRUCTIONS FOR USE NUMBER AND THE CORRESPONDING VERSION ON THE PRODUCT PACKAGING.

TO OBTAIN A PRINTED COPY OF THE INSTRUCTIONS FOR USE AT NO ADDITIONAL COST, PLEASE CONTACT CUSTOMER SERVICE: SAC (21) 3907 2534 / 0800 015 1414 / sac@kovalent.com.br

PRESENTATION

Catalog No.	Package Size
1170050KWS	R1: 2 x 25 mL
1170100KWS	R1: 4 x 25 mL

INTENDED USE

In vitro diagnostic reagent for the quantitative determination of Zinc in serum, plasma, or urine using photometric systems.

SUMMARY[1]

Zinc is involved in many enzymatic reactions at the molecular level. It plays an important role in the synthesis of DNA and RNA and exerts a clearly enhancing effect on the immune system. Another important function of zinc is its involvement in the cellular protective function against free radicals and reactive oxygen compounds.
Causes of zinc deficiency may be – among others – malnutrition, malabsorption, diseases of the small intestine, alcoholism, diabetes mellitus, rheumatic disorders, acute and chronic infections, or chronic liver diseases.

METHOD

Colorimetric assay using 5-Br-PAPS.

PRINCIPLE

Zinc forms with 2-(5-Brom-2-pyridylazo)-5-(N-propyl-N-sulfo-propylamino)-phenol a red chelate complex.The increase in absorbance of this complex is proportional to the zinc concentration in the sample.

REAGENTS

Components and Concentrations

Bicarbonate buffer	pH 9,4	200 mmol/L
5-Br-PAPS		0,02 mmol/L
Sodium citrate		170 mmol/L
Dimethylglyoxime		4 mmol/L
Detergent		1%

REAGENT STORAGE AND STABILITY

The reagents are stable until the stated expiration date when contamination is avoided, protected from light, and stored at 15–25 °C.

After opening, the reagent is stable for 30 days when stored at 15–25 °C.

WARNINGS AND PRECAUTIONS

1. For *in vitro* diagnostic use only.
2. Consult the Safety Data Sheet and observe the precautions required for handling laboratory reagents.
3. For a final diagnosis, results must always be interpreted in conjunction with the patient's medical history, clinical examination, and other findings.
4. In the event of a malfunction or a change in product appearance that may affect performance, contact the manufacturer.
5. Any serious incident related to the product must be reported to the manufacturer and to the competent authority in the jurisdiction where the user and/or patient is located.
6. Although verification and validation studies evaluated the effects of potential interfering substances in this assay, interference from unknown substances may occur. Therefore, as with any diagnostic test, technical errors, procedural errors, and unlisted substances or factors may adversely affect the proper performance of the test kit.
7. For professional use only.

WASTE DISPOSAL

Comply with the applicable regulations governing the technical requirements for healthcare waste management, as well as with equivalent biosafety practices.

REAGENT PREPARATION

The reagent is ready to use.

MATERIALS REQUIRED BUT NOT PROVIDED

1. 9 g/L NaCl solution.
2. General laboratory equipment.

SPECIMEN[3]

Serum, heparinized plasma (do not use EDTA plasma), or urine.

Serum or plasma stability:	3 days	at	20 – 25 ºC
	7 days	at	4 – 8 ºC
	6 months	at	-20 ºC

Discard contaminated specimens. Freeze only once.

Urine stability:	4 days	at	20 – 25 ºC

ASSAY PROCEDURE

Applications for automated systems are available upon request or on our website: www.kovalent.com.br

Wavelength	560 nm
Temperature	25 °C / 37 ºC
Measurement	Against the reagent blank

Allow reagents and specimens to reach room temperature.

	Blank	Specimen or standard
Specimen or standard	-	50 µL
Distilled water	50 µL	-
Reagent	1000 µL	1000 µL
Mix and incubate for 8 min at 25 °C or for 5 min at 37 °C. Read the absorbance against the reagent blank.		

CALCULATIONS

Using a standard or calibrator

$$\text{Zinc [µg/dL]} = \frac{\Delta A_{\text{Specimen}}}{\Delta A_{\text{Standard/calibrator}}} \times \text{Standard/calibrator conc. [µg/dL]}$$

Conversion factor
Zinc [µg/dL] × 0.153 = Zinc [µmol/L]

CALIBRATORS AND CONTROLS

Use the Topkal Zinco calibrator to calibrate the assay. For internal quality control, Kovalent Topkon N and P controls should be measured. Each laboratory should establish corrective actions when control recovery deviates from the acceptable range.

WARRANTY

These Instructions for Use must be read carefully before using the product, and the information contained herein must be strictly followed. The reliability of assay results cannot be guaranteed if these instructions are not followed.

PERFORMANCE CHARACTERISTICS

Accuracy and Precision
CV ≤ 2% for within-run and between-run precision.

Sensitivity / Limit of Detection
The limit of detection is 2.9 µg/dL (0.444 µmol/L).

Linearity and Measuring Range

The assay was developed to determine zinc concentrations from 2.9 to 500 µg/dL (0.444 – 76.5 µmol/L). Samples with values above this range should be diluted 1 + 1 with 9 g/L NaCl solution, and the result should be multiplied by 2.

Specificity / Interferences

No interference was observed from bilirubin up to 15 mg/dL, lipemia up to 1000 mg/dL triglycerides, or hemoglobin up to 160 mg/dL.

Method Comparison

A method comparison between Zinco 5-Br-PAPS (y) and a commercially available assay (X), using 66 samples, yielded the following results: y = 0.9663x + 2.6613; r = 0.9911.

REFERENCE VALUES

Serum or plasma

	µg/dL (µmol/L)
Adults	46 – 150 (7 – 23)
Children	
< 4 months	65 – 137 (10 – 21)
4–12 months	65 – 130 (10 – 20)
1–5 years	65 – 118 (10 – 18)
6–9 years	78 – 105 (12 – 16)

	Male	**Female**
10–13 years	78 – 98 (12 – 15)	78 – 118 (12 – 18)
14–19 years	65 – 118 (10 – 18)	59 – 98 (9 – 15)

Urine

150 – 800	µg/24h	24-hour urine
15 – 120	µg/dL	Spot urine

*Each laboratory should verify whether these reference values are applicable to its patient population and establish its own reference values, if necessary.

REFERENCES

1. Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 347-9.
2. Johnsen and R.Eliasson. Evaluation of a commercially available kit for the colorimetric determination of zinc. International Journal of Andrology, 1987, April 10 (2): 435-440.
3. Guder WG, Zatwa B et al. The quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag, 2001: 38-39, 50-51.

INFORMATION FOR USERS

Symbols used

	Manufacturer
	Temperature limit
	*In vitro* diagnostic medical device
	Caution
	Consult instructions for use
	Recyclable material
	Do not dispose of directly into the environment
	Batch code
	Date of manufacture
	Use-by date
	Biological risks
	Highly toxic
	Corrosive
	Harmful

MANUFACTURER

Kovalent do Brasil Ltda.
110 Cristóvão Sardinha Street – Jardim Bom Retiro
São Gonçalo – RJ – Postal Code 24722-414 – Brazil
www.kovalent.com.br
CNPJ: 04.842.199/0001-56

SAC: sac@kovalent.com.br - (21) 3907-2534 / 0800 015 1414

Expiration date and batch number: SEE LABEL