Alfa-1 Antitripsina

Anvisa 80115310303

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Alpha-1 Antitrypsin

Anvisa 80115310303

BEFORE USING THE PRODUCT, CHECK THE INSTRUCTIONS FOR USE NUMBER AND THE CORRESPONDING VERSION ON THE PACKAGING.

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT AT NO ADDITIONAL COST, CONTACT CUSTOMER SERVICE: SAC (21) 3907 2534 / 0800 015 1414 / sac@kovalent.com.brORDER INFORMATION

PRESENTATION

Cat. No.	Kit size
4300055K	R1: 2 x 25 mL + R2: 1 x 5 mL
4300050KS	R1: 1 x 40 mL + R2: 1 x 10 mL

INTENDED USE

Reactive for the *in vitro* quantitative determination of α-1 antitrypsin (AAT) in human serum by immunoturbidimetric assay.

SUMMARY

α1-Antitrypsin is an acute phase protein. It inhibits proteinase and serine proteases. AAT has a strong binding constant for leukocyte elastase.
Increased serum levels are found in: acute infection and inflammation, acute malaria, pregnancy (in 100 %), anabolic steroid therapy, advanced malignant tumours. Decreased serum levels are found in: congenital deficiencies, juvenile cirrhosis, lung emphysema, testosterone administration.

METHOD

Measurement of antigen-antibody reaction by the end-point method.

REAGENTS

Components and Concentrations
R1 (Buffer):
Phosphate buffered saline
Enhancer
Sodium azide (0.95 g/L)

R2 (Antiserum):
Phosphate buffered saline
Polyclonal goat anti-human α1-Antitrypsin
Sodium azide (0.95 g/L)

STORAGE AND STABILITY

The reagents are stable until the expiration date indicated on the kit, as long as they are stored at a temperature of 2 to 8 °C and contamination is avoided. The stability of the equipment is at least 4 weeks, as long as contamination is avoided. Do not freeze.

WARNINGS AND PRECAUTIONS

1. For in vitro diagnostic use only.
2. Sodium azide has been reported to form lead or copper azide in laboratory plumbing, potentially causing explosions. After dispensing solutions containing sodium azide, flush with copious amounts of water to ensure complete dilution.
3. The enhancing polymer (polyethylene glycol) poses no biological risk.
4. Exercise appropriate caution when handling laboratory reagents.
5. Each donor used in the preparation of standards and controls tested negative for HIV-1 and HIV-2 antibodies, hepatitis B surface antigen, and anti-hepatitis C antibodies using FDA-approved methods. Only donors with negative results were used in manufacturing. However, all products derived from human body fluids must be handled with appropriate care, in accordance with recommended procedures for biohazardous materials, as the absence of infectious agents cannot be guaranteed.
6. In the event of product malfunction or a change in appearance that could affect performance, contact the manufacturer.
7. Any serious incident related to the product must be reported to the manufacturer and the competent authorities in the state where the user or patient is located.
8. Consult the Safety Data Sheet (SDS) and take necessary precautions for the use of laboratory reagents. For diagnostic purposes, results should always be evaluated in conjunction with the patient's medical history, clinical examinations, and other findings.
9. For professional use only.

WASTE MANAGEMENT

Follow the requirements of the current guidelines about technical regulation for the management of healthcare service waste, as well as other equivalent biosafety practices.

REAGENT PREPARATION

The reagents are liquid and ready for use.

MATERIALS REQUIRED, BUT NOT PROVIDED

General laboratory equipment
Sodium Chloride (0.9%)
Calibrators and controls

SPECIMEN

Human serum.

If the test cannot be performed on the same day, the serum must be stored at 2–8°C for up to 48 hours. If a longer storage period is required, the sample must be frozen.

Freeze only once. Discard contaminated samples.

ASSAY PROCEDURE

Applications for automated systems are available upon request.
Wavelength: 340 nm

Sample, control or calibrator	2 µL
Reagent 1	350 µL
Mix, incubate for 2 minutes, read the absorbance, and then add:	
Reagent 2	70 µL
Mix, incubate for 5 minutes, and read the absorbance.	

WARRANTY

These instructions for use should be read carefully before using the product and the information contained therein should be strictly adhered to. The reliability of the test results cannot be guaranteed if the instructions are not followed.

PERFORMANCE CHARACTERISTICS

To evaluate the performance characteristics of the Alpha-1 Antitrypsin kit, studies were conducted on a clinical chemistry analyzer.

Measurement Range: 0 – 400 mg/dL
Limit of Detection: 8 mg/dL
Hook Effect: > 800 mg/dL
Sensitivity: 0.0013 ABS units/concentration units

Precision [%CV]:

	Low	Medium	High
Within-run	1.51	1.66	4.79
Between-run	-	3.44	-

Accuracy [mg/dL]:

Control	Target value	Measured value
Bio-Rad 1	79 (63–94)	88
Bio-Rad 2	176 (141–211)	174

Specificity: Monospecific
Interferents: No interference from: Hemoglobin (1000 mg/dL), Sodium Citrate (1000 mg/dL), Heparin (50 mg/dL), Bilirubin (20 mg/dL), and Triglycerides (2500 mg/dL).
Limitations: None
Comparison with Nephelometry: y = 1.0041X + 2.9236 / r = 0.9934

REFERENCE VALUES

89 – 205 mg/dL (IFCC)
Each laboratory should verify whether the reference interval applies to its own patient population and, if necessary, determine its own reference intervals.

LITERATURE

1. Naito, H.K., J. Clin. Immunoassay, 9, 155 (1986)
2. Kottke, B.A., et. al., Mayo Clin. Proc.,61, 313 (1986)
3. Dati, F. et al., Lab. Med. 13, 87 (1989)

CONSUMER INFORMATION

Symbols used:

	Manufacturer
	Temperature limit
	In vitro diagnostic device
	Caution
	Operating instructions
	Recycling material
	Do not discard directly into the environment
	Batch code
	Date of manufacture
	Use by date
	Biological hazards
	Highly toxic
	Corrosive
	Harmful

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Customer service: sac@kovalent.com.br - (21) 3907-2534 / 0800 015 1414

Expiration date and Lot no.: See label